



# BioSim<sup>™</sup> Etanercept (Human) ELISA Kit

rev 12/20

(Catalog # E4374-100, 100 assays, Store at 4°C)

## I. Introduction:

Etanercept is a therapeutic fusion protein specific for Tumor Necrosis Factor-Alpha (TNF-α) and is used to treat rheumatic arthritis, intestinal disorders, dermatological diseases and cancer. Etanercept specifically binds to TNF alpha and blocks its interaction with cell surface TNF receptors and reduces inflammation and subsequently improves the patient's health. Drug level quantification can be important to adapt to patient prescription or to switch to an alternative TNF inhibitor drug. BioVision's BioSim™ Etanercept ELISA kit has been developed for specific quantification of Etanercept concentration in human serum or plasma with high sensitivity and reproducibility.

#### II. Application:

This ELISA kit is used for in vitro quantitative determination of Etanercept.

Detection Range: 0.1 - 3 µg/ml

Sensitivity: 0.1 µg/ml

Assay Precision: Intra-Assay: CV < 30%; Inter-Assay: CV < 30% (CV (%) = SD/mean X 100)

Recovery rate: < 100 ± 30% with normal human serum samples with known concentrations

Cross Reactivity: No significant cross-reactivity or interference with other proteins present in native human serum or other therapeutic immunoglobulins.

### III. Sample Type:

Human serum and plasma

#### IV. Kit Contents:

Components	E4374-100	Part No.
Micro ELISA Plate	1 plate	E4374-100-1
Etanercept Standards (S1 – S7)	0.3 ml X 7	E4374-100-2.x
Assay Buffer	50 ml	E4374-100-3
HRP-conjugate Probe	12 ml	E4374-100-4
TMB substrate (Avoid light)	12 ml	E4374-100-5
Stop Solution	12 ml	E4374-100-6
Wash buffer (20X)	50 ml	E4374-100-7
Plate sealers	2	E4374-100-8

# V. User Supplied Reagents and Equipment:

- Microplate reader capable of measuring absorbance at 450 nm
- Calibrated measures
- · Precision pipettes with disposable tips
- · Clean eppendorf tubes for preparing standards or sample dilutions
- Absorbent paper

#### VI. Storage and Handling:

The entire kit may be stored at 4°C for up to 12 months from the date of shipment.

### VII. Reagent and Sample Preparation:

Note: Prepare reagents within 30 minutes before the experiment.

Before using the kit, spin tubes and bring down all components to the bottom of tubes.

1. Wash Buffer: Dilute the 20X Wash Buffer to 1X solution in ddH<sub>2</sub>O (10 ml of Wash Buffer stock to 190 ml of ddH<sub>2</sub>O). Mix the 1X solution thoroughly by vortex manually. The working stock can be stable for 2 weeks after preparation at 4°C.

#### 2. Standard Preparation:

Ready to use

Name	S1	S2	S3	S4	<b>S</b> 5	S6	<b>S</b> 7
Conc. (µg/ml)	3	1	0.3	0.1	0	High Control	Low Control

#### 3. Sample Dilution:

- Serum/Plasma: Initially dilute samples at 1:20 (10 µl Serum/Plasma+ 190 µl Assay Buffer). (Dilution Buffer 1:20)
- Diluted samples should further be diluted if the concentration of Etanercept is higher than the measuring range.
- The usual precautions for venipuncture should be observed. Samples are stable at 4°C for 2 days and -20°C for 6 months. Avoid freeze-and-thaw cycle.

# VIII. Assay Protocol:

Note: Bring all reagents, microplate and samples to room temperature 15 minutes prior to the assay. It is recommended that all standards and samples be run at least in duplicate.



#### Gentaur Europe BVBA Voortstraat 49, 1910 Kampenhout BELGIUM Tel 0032 16 58 90 45 info@gentaur.com



A standard curve must be run with each assay.

- 1. Prepare all reagents, samples and standards as instructed in section VII.
- 2. Pipette 100µl of **Assay Buffer** non-exceptionally into each of the wells to be used.
- 3. Add 20 µl of standards, controls and diluted-samples into appropriate wells. Cover wells and incubate for 30 minutes at room temperature (RT).
- 4. Discard incubation solution. Wash plate 3 times each with 300 µl of diluted Wash Buffer. Remove excess solution by tapping the inverted plate on a paper towel.
- 5. Add 100 µl of HRP-conjugate into each well. Cover wells with adhesive plate sealer and incubate at RT for 30 minutes.
- 6. Discard the solution and wash the wells as step 4.
- 7. Add 100 µl of 1X TMB substrate solution and incubate the plate in dark at RT for 10 minutes
- 8. Add 100 µl of **Stop solution** to stop the reaction
- 9. Read the absorbance in micro plate reader set to 450 nm within 20 minutes. (reference wavelength to 650 nm)

#### **CALCULATION:** IX.

Using the standards disregarding zero standard, construct a standard curve by plotting the OD450/650 nm for each of 4 standards on the Yaxis versus the corresponding Etanercept concentration on the X-axis. Construct a standard curve of difference data using software capable of generating four parameter logistic (4PL) or point-to-point calculation curve fit. To obtain the exact values of the samples, the concentration determined from the standard-curve should be multiplied by the dilution factor.




Figure: Typical Standard Curve: These standard curves are for demonstration only. A standard curve must be run with each assay.

# X. RELATED PRODUCTS:

- BioSim<sup>™</sup> Rituximab (Human) ELISA Kit (Cat. No. E4371-100)
- BioSim™ Adalimumab (Human) ELISA Kit (Cat. No. E4372-100)
  BioSim™ Bevacizumab (Human) ELISA Kit (Cat. No. E4373-100)

- BioSim<sup>TM</sup> Etanercept (Human) ELISA Kit (Cat. No. E4374-100)
  BioSim<sup>TM</sup> anti-HER2 (Human) ELISA Kit (Cat. No. E4376-100)
  BioSim<sup>TM</sup> Golimumab (Human) ELISA Kit (Cat. No. E4377-100)
- BioSim<sup>™</sup> Cetuximab (Human) ELISA Kit (Cat. No. E4379-100)
- BioSim<sup>™</sup> Denosumab (Human) ELISA Kit (Cat. No. E4380-100)
- BioSim<sup>™</sup> Omalizumab (Human) ELISA Kit (Cat. No. E4381-100)
- BioSim<sup>™</sup> Nivolumab (Human) ÉLISA Kit (Cat. No. E4382-100)
- BioSim<sup>™</sup> Pembrolizumab (Human) ELISA Kit (Cat. No. E4383-100)
- BioSim<sup>™</sup> Ipilimumab (Human) ELISA Kit (Cat. No. E4384-100)