



# BioSim<sup>™</sup> anti- HER2 (Human) ELISA Kit

rev 06/20

# (Catalog # E4376-100, 100 assays, Store at 4°C)

#### I. Introduction:

Anti-HER2 is a recombinant DNA-derived humanized monoclonal antibody that selectively targets the extracellular domain of the human epidermal growth factor receptor 2 protein (HER2). It has antitumor activity against HER2-positive human breast tumor cells in laboratory models and is active for the treatment of women with HER2-overexpressing breast cancers. This antibody was approved in 1998 for clinical use for HER2 overexpressing metastatic breast cancer. In HER2 overexpressing cells, anti- HER2 antibody markedly down-regulates HER2 expression by accelerating receptor endocytosis and degradation and inhibits cell cycle progression by inducing the formation of p27Kip1/Cdk2 complexes. BioSim<sup>™</sup> anti-HER2 ELISA kit has been developed for specific quantification of anti-HER2 antibody concentration in human serum or plasma with high sensitivity and reproducibility.

#### II. Application:

This ELISA kit is used for *in vitro* quantitative determination of anti- HER2 antibody Detection Range: 10 - 300 ng/ml

Sensitivity: Quantitative limit - 10 ng/ml, Detection limit - 2 ng/ml

Assay Precision: Intra-Assay: CV < 30%; Inter-Assay: CV < 30% (CV (%) = SD/mean X 100)

Recovery rate: <100±30% with normal human serum samples with known concentrations

Cross Reactivity: No significant cross-reactivity or interference with other proteins present in native human serum or other therapeutic immunoglobulins.

#### III. Sample Type:

Human serum and plasma

#### IV. Kit Contents:

Components	E4376-100	Part No.	
Micro ELISA Plate	1 plate	E4376-100-1	
anti- HER2 Standards (S1 – S7)	0.3 ml X 7	E4376-100-2.x	
Assay Buffer	2 X 50 ml	E4376-100-3	
HRP-conjugate Probe	12 ml	E4376-100-4	
TMB substrate (Avoid light)	12 ml	E4376-100-5	
Stop Solution	12 ml	E4376-100-6	
Wash buffer (20X)	50 ml	E4376-100-7	
Plate sealers	2	E4376-100-8	

# V. User Supplied Reagents and Equipment:

- Microplate reader capable of measuring absorbance at 450 nm
- · Calibrated measures
- Precision pipettes with disposable tips
- Clean eppendorf tubes for preparing standards or sample dilutions
- · Absorbent paper

#### VI. Storage and Handling:

The entire kit may be stored at 4°C for up to 12 months from the date of shipment.

#### VII. Reagent and Sample Preparation:

Note: Prepare reagents within 30 minutes before the experiment.

Before using the kit, spin tubes and bring down all components to the bottom of tubes.

- 1. Assay Buffer: Dilute 5X assay buffer to 1X in ddH<sub>2</sub>O (25 ml of Assay Buffer stock to 100 ml of ddH<sub>2</sub>O)
- 2. Wash Buffer: Dilute the 20X Wash Buffer to 1X solution in ddH<sub>2</sub>O (10 ml of Wash Buffer stock to 190 ml of ddH<sub>2</sub>O). Mix the 1X solution thoroughly by vortex manually. The working stock can be stable for 2 weeks after preparation at 4°C.

#### 3. Standard Preparation:

Dilute 10X stock with Assay Buffer. (20 µl Standards + 180 µl Assay Buffer)

Name	S1	S2	<b>S</b> 3	S4	S5	S6	S7
Conc. (ng/ml)	3000	1000	300	100	0	High Control	Low Control
Working Conc. (ng/ml)	300	100	30	10	0	-	-





#### 4. Sample Dilution:

- Serum/Plasma: Initially dilute samples 1:10 (10 μl Serum/Plasma+ 90 μl Assay Buffer). Then Dilute another 1:100 (5 μl Standard + 495 μl Assay Buffer) to a total of 1:1000 dilution.
- Diluted samples should further be diluted if the concentration of anti- HER2 is higher than the measuring range.
- The usual precautions for venipuncture should be observed. Samples are stable at 4°C for 2 days and -20°C for 6 months. Avoid freeze-and-thaw cycle.

#### VIII. Assay Protocol:

Note: Bring all reagents, microplate and samples to room temperature 15 minutes prior to the assay.

- It is recommended that all standards and samples be run at least in duplicate.
- A standard curve must be run with each assay.
- 1. Prepare all reagents, samples and standards as instructed in section VII.
- 2. Add 100 µl of **standards** and **diluted-samples** into appropriate wells. Cover wells and incubate for 30 minutes at room temperature (RT).
- 3. Discard incubation solution. Wash plate 3 times each with 300 µl of diluted **Wash Buffer**. Remove excess solution by tapping the inverted plate on a paper towel.
- 4. Add 100 µl of HRP-conjugate into each well. Cover wells with adhesive plate sealer and incubate at RT for 30 minutes.
- 5. Discard the solution and wash the wells as step 3.
- 6. Add 100 µl of 1X TMB substrate solution and incubate the plate in dark at RT for 10 minutes
- 7. Add 100 µl of Stop solution to stop the reaction
- 8. Read the absorbance in micro plate reader set to 450 nm within 20 minutes. (reference wavelength to 650 nm)

### IX. CALCULATION:

Using the standards disregarding zero standard, construct a standard curve by plotting the OD450/650 nm for each of 4 standards on the Yaxis versus the corresponding anti- HER2 concentration on the X-axis. Construct a standard curve of difference data using software capable of generating four parameter logistic (4PL) or point-to-point calculation curve fit. To obtain the exact values of the samples, the concentration determined from the standard-curve should be multiplied by the dilution factor.




Figure: Typical Standard Curve: These standard curves are for demonstration only. A standard curve must be run with each assay.

# X. RELATED PRODUCTS:

- BioSim<sup>™</sup> Rituximab (Human) ELISA Kit (Cat. No. E4371-100)
- BioSim<sup>™</sup> Adalimumab (Human) ELISA Kit (Cat. No. E4372-100)
- BioSim<sup>™</sup> Bevacizumab (Human) ELISA Kit (Cat. No. E4373-100)
- BioSim<sup>™</sup> Etanercept (Human) ELISA Kit (Cat. No. E4374-100)
- BioSim™ Ipilimumab (Human) ELISA Kit (Cat. No. E4375-100)
- BioSim<sup>™</sup> Golimumab (Human) ELISA Kit (Cat. No. E4377-100)
- BioSim™ Infliximab (Human) ÉLISA Kit (Cat. No. E4378-100)
- BioSim<sup>™</sup> Cetuximab (Human) ELISA Kit (Cat. No. E4379-100)
- BioSim<sup>™</sup> Denosumab (Human) ELISA Kit (Cat. No. E4380-100)
- BioSim<sup>™</sup> Omalizumab (Human) ELISA Kit (Cat. No. E4381-100)
- BioSim<sup>™</sup> Nivolumab (Human) ELISA Kit (Cat. No. E4382-100)
- BioSim<sup>™</sup> Pembrolizumab (Human) ELISA Kit (Cat. No. E4383-100)
- BioSim<sup>™</sup> Ipilimumab (Human) ELISA Kit (Cat. No. E4384-100)