



3. After the incubation at 37 °C, incubate the sample at 85-90 °C for 10 min.
4. Dilute the sample 40 fold using the Sample Diluent. (For example, mix 5 µl of serum with 195 µl of Sample Diluent.)
5. Use 50 µl per well for the assay.

Note: Dilution factor: 40

• **Urine and Saliva**

1. Centrifuge 0.5 ml of urine or 0.2 ml of saliva at 10,000 g for 5 min and recover the supernatant.
2. Dilute the supernatant 40 fold using the Sample Diluent. (For example, mix 5 µl of urine with 195 µl of Sample Diluent.)
3. Use 50 µl per well for the assay.

Note: Dilution factor: 40

IX. Caffeine ELISA Assay Protocol:


Notes: It is recommended that all standards and samples should be run at least in duplicate. Standard curves must be run each time as reference for sample quantification.

1. Prepare all reagents, standards and samples as sections VII and VIII respectively.
2. Add 50 µl of Standards or Samples per well. Add 50 µl of conjugate working solution and 50 µl of Antibody to all wells containing standard or sample.
3. Cover the microtiter plate with plate sealer and mix well. Incubate the plate at room temperature (25 °C) for 45 min.
4. Aspirate all reagents and wash each well 4 times: add 250 µl of 1X Wash Buffer and incubate for 30 sec. Remove 1X Wash buffer completely before the next wash. (This is essential for accurate results.) Repeat this step 3 more times.
5. Add 100 µl of TMB Substrate to each well. Tap or shake the plate to ensure complete mixing.
6. Check the OD at 650 nm for the well containing no caffeine (S0). When its reading is approximately between 0.8 and 1.0 (usually between 5-30 min after addition of TMB Substrate), add 50 µl of Stop Solution and gently tap the plate to ensure thorough mixing.
7. Measure OD at 450 nm for the standards and samples.


X. Calculation:

The Standard Curve is prepared by plotting OD at 450 nm vs. caffeine concentrations. The concentration of caffeine in each sample (ng/ml), which can be read from the calibration curve, is multiplied by the corresponding dilution factor.

A.



B.



Figures. A. Caffeine Standard Curve (*This standard curve is for demonstration only. A standard curve must be run with each assay*). **B.** Spike recovery experiment: Human serum, urine and saliva samples were assayed alone or with a spike to a final expected concentration of 500 ng/ml of caffeine (80-95% recovery of spike).

XI. RELATED PRODUCTS:

Gentamicin (serum/urine) ELISA Kit (Cat. No. K4315-100)
Ampicillin ELISA Kit (Cat. No. E4350-100)
Enrofloxacin (ENR) ELISA Kit (Cat. No. E4277-100)
Quinolone ELISA Kit (Cat. No. E4530-100)

Folic Acid ELISA Kit (Cat. No. 4523-100)
Kanamycin ELISA Kit (Cat. No. K4210-100)
Cell Lysis Buffer (Cat. No. 1067-100)

FOR RESEARCH USE ONLY! Not to be used on humans.