



Sudan ELISA Kit

10/19

(Catalog # E4769-100; 96 assays, Storage at 4°C)

I. Introduction:

- II. Sudan dyes are synthetic chemical dyes of similar chemical structure. They are aromatic compounds containing azo group (- N=N -). Sudan dyes are widely used in commercial manufacturing to impart a rich red, red-orange or yellow-orange color to plastics and textiles. They can be generally applied for coloring substances such as hydrocarbon solvents, oils, fats, waxes and plastics. Although Sudan dyes have been reported as contact allergens and sensitizers, the greatest concern has been on their possible carcinogenicity. BioVision's Sudan ELISA Kit is based on Competitive ELISA method. The microtiter plate provided in the kit has been pre-coated with coupled antigen. During the reaction, Sudan dye in the samples or standard competes with coupled antigen on the solid phase supporter for sites of anti-Sudan antibody. Then Horseradish Peroxidase (HRP) conjugate is added to each microtiter plate well, and TMB substrate is added for color development. There is a negative correlation between the OD value of samples and the concentration of Sudan. The concentration of Sudan in the samples can be calculated by comparing the OD of the samples to the standard curve.

III. Applications:

In vitro, quantitative determination of Sudan

Detection Range: Tomato juice, Ketchup, Chilli sauce - 12ppb, Chilli powder, Feed -120ppb,

Eggs (Chicken egg, Duck egg, Goose egg) - 30ppb

Sensitivity: 0.3 ppb (ng/mL)

Sample recovery rate: Tomato juice, Ketchup, Chilli sauce - 80%±15%, Chilli powder, Feed - 95%±15%,

Eggs(Chicken egg, Duck egg, Goose egg) - 80%±15%

IV. Sample Type:

Tomato juice, Ketchup, Chilli sauce, Paprika, Feed, Eggs

V. Kit Contents:

Components	E4769-100	Part Number
Micro ELISA Plate	96 wells	E4769-100-1
High concentrated Standard (1.0 ppm)	1 ml	E4769-100-2
HRP Conjugate	11 ml	E4769-100-3
Antibody Working Solution	5.5 ml	E4769-100-4
Substrate Reagent A	6 ml	E4769-100-5
Substrate Reagent B	6 ml	E4769-100-6
Stop Solution	6 ml	E4769-100-7
Wash Buffer (20X)	40 ml	E4769-100-8
Empty vials	6	E4769-100-9
Plate Sealer	3	E4769-100-10

VI. User Supplied Reagents and Equipment:

- Microplate reader capable of measuring absorbance at 450 nm
- Methanol
- Clean eppendorf tubes for preparing standards or sample dilutions

VII. Storage and Handling:

Store at 4°C.

VIII. Reagent and Sample Preparation:

Bring all reagents to room temperature before use. Before using the kit, spin tubes and bring down all components to the bottom of tubes.

- **Wash Buffer (20X):** Dilute 20X Concentrated Wash Buffer to 1X with deionized water.
- **10% Methanol:** Add 10 ml Methanol to 90 ml deionized water, mix well.
- **Standard:** Prepare fresh standards each time. Label the empty bottles provided as 0 ppb, 0.3 ppb, 0.9 ppb, 2.7 ppb, 8.1 ppb and 24.3 ppb. Add 3 ml of 10% Methanol into vial label as 0 ppb. Add 2 ml of 10% Methanol into standard vial label as 0.3 ppb, 0.9 ppb, 2.7 ppb, and 8.1 ppb respectively. Add 2.93 ml of 10% Methanol into 24.3 ppb vial. To prepare standards, mix high concentrated standard (1.0 ppb) with methanol as described below:
 - Standard 6 (24.3 ppb):** Add 73 µl of 1.0 ppm high concentrated standard into 24.3 ppb vial, mix well.
 - Standard 5 (8.1 ppb):** Add 1 ml of Standard Solution 6 into 8.1 ppb vial, mix well.
 - Standard 4 (2.7ppb):** Add 1 ml of Standard Solution 5 into 2.7ppb vial, mix well.
 - Standard 3 (0.9 ppb):** Add 1 mL of Standard Solution 4 into 0.9 ppb vial, mix well.
 - Standard 2 (0.3 ppb):** Take 1 mL of Standard Solution 3 into 0.3 ppb vial, mix well.
 - Standard 1 (0 ppb):** Use 10% Methanol only.

