

BasicExoRNA™ Basic RNA Extraction Kit

(Cat#: K1222-50; Store at +4 °C)

I. Introduction:

Exosomes are small endosome derived lipid nanoparticles (50-120 nm) actively secreted by exocytosis by most living cells. Exosome release occurs either constitutively or upon induction, under both normal and pathological conditions, in a dynamic, regulated and functionally relevant manner. Both the amount and molecular composition of released exosomes depend on the state of a parent cell. Exosomes have been isolated from diverse cell lines (hematopoietic cells, tumor lines, primary cultures, and virus infected cells) as well as from biological fluids in particular blood (e.g. serum and plasma from cancer patients) and other body fluids (broncho alveolar lavage fluid, pleural effusions, synovial fluid, urine, amniotic fluid, semen, saliva etc). Exosomes have pleiotropic physiological and pathological functions and an emerging role in diverse pathological conditions such as cancer, infectious and neurodegenerative diseases. Exosomes shuttle functional RNA molecules in the target cell. Increasing evidence suggests a role for exosome-derived miRNAs in the development and/or progression of specific human diseases. Pathogenic miRNAs might be exploited as novel therapeutic targets or disease biomarkers in complex diseases, including cancer. In fact, miRNAs seem to play critical roles as transcriptional and post-transcriptional regulators of epigenetic mechanisms and cell processes and have been linked to the etiology, progression and prognosis of cancer. Similar miRNA expression patterns between tumor tissue samples and circulating exosomes have been observed.

BasicExoRNA™ Kit (basic RNA extraction kit) allows total RNA (miRNA + mRNAs) extraction from pre-isolated exosomes pre-isolated via different methods including ultracentrifugation, chemical precipitation, immunocapture, size-chromatography etc). *RNA basic kit does not contain immunobeads for exosome isolation. Pre-isolated exosomes are lysed with an optimized lysis buffer and total RNA is purified using spin columns with a fast and user-friendly process. Eluted RNA can be used for downstream analyses or stored at -80 °C. All our kits guarantee high specificity for exosomal RNA and high yield of total RNA (including small RNAs) than similar products.




Figure 1. Exosomes RNA Extraction Kit.

II. Applications:

- Exosome RNA extraction from pre-isolated exosomes pre-isolated via different methods including ultracentrifugation, chemical precipitation, immunocapture, size-chromatography etc
- Simultaneous miRNA and mRNA profiling (qRT-PCR, RT-PCR, microarray).

III. Sample Type:

- Pre-isolated exosomes from different samples (biofluids or cell culture media, either pellet after precipitation, frozen, or lyophilized).

IV. Kit Contents (for RNA extraction from Plasma, Serum, Urine, Cell culture media, pre-isolated exosomes):

Components	K1222-50	Part Number	Storage Temperature
	50 reactions		
Lysis buffer	1 bottle	K1222-50-1	+4 °C
RNA Washing buffer* (add Ethanol 96%)	2 bottles	K1222-50-2	+4 °C
Elution buffer	1 vial	K1222-50-3	+4 °C
Columns	52 columns	K1222-50-4	+4 °C or RT
Elution tubes	52 tubes	K1222-50-5	+4 °C or RT

*Add into the bottle containing RNA Washing Buffer, the volume of pure ethanol (96%) indicated on the bottle's label to get the final ethanol concentration of approximately 70%, which is around 48 ml of 96% Ethanol to the RNA Washing Buffer bottle.

V. User Supplied Reagents and Equipment:

- Single-use and/or pipettes with disposable tips 2-100 µl
- Pipettes 1 ml and 5 ml for reagent preparation
- PBS
- Disposable pipetting reservoirs
- Ethanol 96%
- Chloroform or BCP (1-bromo-3-chloropropane)
- Sample concentrator (urine and cell culture supernatant samples)

VI. Shipment and Storage:

All the reagents and buffers with the RNA extraction kit are shipped and stored at +4 °C for up to 24 months, if unopened. Spin columns and Elution tubes may be stored at room temperature (RT). DO NOT FREEZE!

VII. Reagent Preparation and Storage Conditions:

- RNA Washing Buffer. Add into the bottle containing RNA Washing Buffer the volume of pure ethanol (96%) mentioned above to get the final ethanol concentration of approximately 70% ethanol.
- Elution buffer and Lysis buffer are ready to use.
- The RNase free column and elution tubes should be stored at RT.
- All the opened buffers and diluted reagents including the bead washing buffer, RNA washing buffer, lysis buffer and the elution buffer should be stored at 4°C.

VIII. BasicExoRNA™ Assay Protocol:

1. RNA Extraction:

a) Lysis:

- Add 700 µl of Lysis buffer directly onto the exosome preparation.
- Resuspend by pipetting up and down until the lysate is clear.
- Incubate for 5 min at RT.

b) Extraction:

- Add 70 µl of 1-Bromo-3-chloropropane (BCP) or 140 µl of pure Chloroform.
- Shake for 30 sec.
- Incubate for 10 min at RT.
- Incubate for 1 min on ice and centrifuge at 12,000 g at 4 °C for 10 min. **Note:** Incubation on ice prior to centrifugation helps in reducing DNA contamination, which tend to remain in the interphase.
- Transfer the top phase (aqueous) to a fresh tube.
- Add 2X of ethanol 96%. Mix by gently inverting 4-5 times. If the top phase volume is 400 µl, add 800 µl of ethanol 96%.

c) Purification:

- Transfer the half volume of the mixture into the spin column.
- Spin at 14,000 g for 30 sec.
- Discard the flow-through.
- Add the remaining volume into the same spin column.
- Spin at 14,000 g for 30 sec. Discard the flow-through.
- Wash the column with RNA Washing buffer. Add 400 µl of RNA Washing buffer in the column. Spin at 14,000 g for 30 sec. Discard the flow-through. Perform the washing step twice more.
- Spin for 5 additional min at 14,000 g to eliminate any ethanol residues from the column. Discard the flow-through.
- Remove the tube and transfer the spin column into an elution tube.
- Elute the column with 15 µl of Elution buffer. Incubate for 5 min at RT. Spin for 2 min at 200g and 1 min at 14,000 g. Keep the flow-through. Eluted RNA is now ready for downstream analysis or for storage at -80 °C

IX. Sensitivity: Purified exosome RNA can be quantified and analyzed using NanoDrop spectrophotometer (Thermo Scientific), although the measured concentration values are likely to end toward the bottom limit of detection of the instrument. For better quantification, we recommend the concomitant use of electropherogram-based technologies (eg Bioanalyzer, Agilent Technologies) or fluorimetric technologies (Qubit nano; ThermoScientific). Since most of the RNA contained in extracellular vesicles are small-non-coding RNAs (eg. miRNA), the expected Nanodrop profile, purity and yield are as shown in the representative Figure 2. ExoRNA™ allows extraction of high quality of exosome-derived RNAs from low volumes of sample and better performance than competitors. Efficiency of ExoRNA™ kit was tested vs a competitor kit for RNA extraction from healthy donor plasma derived exosomes. RNA yield was quantified by Nanodrop (Figure 3) and extracted RNA was subsequently retrotranscribed using the miScript II RT kit (Qiagen). miRNA21 was amplified by qPCR (Figure 4).




Figure 2. Expected Nanodrop profile for RNA extracted from immunocaptured exosomes (100 µl of human plasma). Yield = 8.4 ng/µl; $A_{260/280} = 1.6$; $A_{260/230} = 1.85$




Figure 3. Nanodrop quantification of total RNA yield




Figure 4. miRNA 21 amplification by qPCR

