





## VI. Shipment and Storage:

All the reagents are shipped and stored at 4°C for up to 12 months, if unopened. Briefly centrifuge small vials prior to opening. **DO NOT FREEZE!**

## VII. Reagent Preparation and Storage Conditions:

- ExoFACS™ contains reagents and antibodies for 20 reactions.
- Exosome standards must be reconstituted in 100 µl of deionized water.
- Isolation Component is included in the kit for exosome isolation
- Beads are ready to use for exosome capture.
- Primary and secondary antibody must be appropriately diluted in sample buffer.
- 1 vial (100 µg) of exosome standards (lyophilized), from human serum (number of particles/ml  $1 \times 10^{10}$ ).
- Exosome standards: The remaining reconstituted standard stock solution should be aliquoted into polypropylene vials (preferably low binding) and stored at -20°C for up to one month or at -80°C for up to six months. Strictly avoid repeated freeze-and-thaw cycles.
- Store opened and diluted reagents at 4°C up to 12 months if unopened.

## VIII. ExoQuant™ Assay Protocol:

- 1. Human Serum sample preparation:** Prepare samples by 3 centrifugation steps to eliminate red blood cells and cellular debris. After each step, transfer the supernatant to a new tube and discard the pellet.
  - a) 10 min at 300g at 4°C (save supernatant; discard pellet).
  - b) 20 min at 1200g at 4°C (save supernatant; discard pellet).
  - c) 30 min at 10,000g at 4°C (save supernatant; discard pellet).

- 2. Exosome isolation from Human Serum:**

| Fluid | Minimum volume required | Volume suggested |
|-------|-------------------------|------------------|
| Serum | 200 µl                  | 250 µl -500 µl   |

- a) Add Isolation Component to your sample in ratio 1/4.
  - b) Mix well by pipetting and inverting tube.
  - c) Incubate on ice for 1 hr.
  - d) Centrifuge 20 min at 10,000g (centrifuge can be performed at 4°C or at RT).
  - e) Discard the supernatant.
  - f) Centrifuge for 2 min at 1500g to eliminate entirely the supernatant.
  - g) Resuspend the pellet in 100 µl\* of 1X PBS. \* Volume of resuspension can be defined by the user.
- 3. Lyophilized Exosome Standard reconstitution:**
    - a) Reconstitute lyophilized exosome standard by adding 100 µl of deionized water to get a final concentration of 1 µg/µl.
    - b) Resuspend exosomes pipetting the solution up and down 10-15 times, avoiding bubbles.
    - c) Vortex the reconstituted standard for 60 secs. Briefly centrifuge the tubes containing the standard to ensure that the solution is collected at the bottom of the tube. Pipette the solution up and down 10 times, avoiding the introduction of bubbles. After this step, the standard is ready to use.
    - d) Use 5 µl of reconstituted Exosome Standard for each reaction.
  - 4. Exosome binding onto latex FACS-Beads:**
    - a) It is recommended to prepare the complex Exosome-Beads (Exo-Beads) in one single tube, then to divide in single reactions before the antibody incubation.
    - b) Latex FACS-Beads are ready to use. Resuspend well FACS-Beads prior to use by vortexing or pipetting several times.
    - c) For each reaction mix together 5 µl of Exosome Standards and 5 µl of FACS-beads in an eppendorf tube (preferably low binding). Mix well by pipetting 5-6 times. Example: if you want to run 10 reactions, mix into the same eppendorf low binding tube 50 µl of Exosome Standards and 25 µl of FACS-Beads.
    - d) For exosome isolated using the Isolation Component, mix 25 µl of FACS-Beads with the volume of resuspended exosomes suggested (volumes are indicative only; the user should define the appropriate volumes on the base of exosome yield). Serum: 20 µl - 50 µl /reaction.
    - e) Incubate for 15 min at room temperature (RT).
    - f) Add 0.7 ml of 1X PBS and incubate in rotator or shaker for 2 hr at RT or overnight (ON) at 4°C.
    - g) Centrifuge the complex Exosomes-Beads (Exo-Beads) for 5 min at 4500g at 4°C and discard the supernatant.
    - h) Add 1 ml of Sample Buffer, resuspend Exo-Beads for 5-6 times and incubate for 15 min at RT.
    - i) Centrifuge for 5 min at 4500g at 4°C, discard the supernatant.
  - 5. Antibody Incubation:**
    - a) Prepare the Washing buffer (not provided in the kit) diluting 2% of FBS (or FCS) in 1X PBS (consider that you need 8 ml of washing buffer for each reaction). Alternatively, if you don't have FBS or FCS, prepare the Washing buffer diluting 0.5% of BSA in 1X PBS. Keep on ice.
    - b) Resuspend the Exo-Beads in 100 µl of Sample buffer for each reaction. Example: if you are running 10 reactions resuspend Exo-Beads in 1 ml of Sample buffer.

