



ExoPure™ Isolation Kit (Stem Cell Media)

(Catalog # K1239-2, -10; Store at RT)

I. Introduction:

Exosomes are small endosome derived lipid nanoparticles (50-120 nm) actively secreted by exocytosis by most living cells. Exosome release occurs either constitutively or upon induction, under both normal and pathological conditions, in a dynamic, regulated and functionally relevant manner. Both the amount and molecular composition of the released exosomes depend on the state of a parent cell. Exosomes have been isolated from diverse cell lines (hematopoietic cells, tumor lines, primary cultures, and virus infected cells) as well as from biological fluids in particular blood (e.g. serum and plasma from cancer patients) and other body fluids (broncho alveolar lavage fluid, pleural effusions, synovial fluid, urine, amniotic fluid, semen, saliva etc). Exosomes have pleiotropic physiological and pathological functions and an emerging role in diverse pathological conditions such as cancer, infectious and neurodegenerative diseases.

ExoPure™ Isolation Kit (Stem Cell Media) is a fast and convenient method of exosome isolation and purification at high yields from stem cell culture media. This kit yields highly pure exosomes using filtration method as compared to exosome precipitation methods.

II. Applications:

- Easy to use: No ultra-centrifugation (<2 hr)
- 10-fold higher yield as compared to other kits or ultracentrifuge method
- Cost effective as compared to antibody bead method
- **Isolates pure exosome (exosome purity >95%)**
- Intact exosomes (good morphology)
- Each reaction can process 20 mL stem cell culture medium. The yield of each reaction can yield pure exosome suspended in 50–200 µL PBS.
- Isolated exosomes are suitable for a wide range of downstream analyses, such as EM study, exosome label, exosome subpopulation, qRT-PCR profiling of exosomal miRNAs, and gel analysis of exosomal proteins western blotting
- Easy to store and ships at room temperature (RT).

III. Sample Type:

- Stem Cell Media

IV. Kit Contents: (exosome Isolation from stem cell media):

Components	K1239-2	K1239-10	Part Number
	2 reactions	10 reactions	
Solution A (Blue)	1.5 mL	7.5 mL	K1239-XX-1
Solution B	1.5 mL	7.5 mL	K1239-XX-2
Solution C	6 mL	10 mL x 3 bottles	K1239-XX-3
ExoPure™ Column	2	10	K1239-XX-4

V. User Supplied Reagents and Equipment:

- Glass tubes

VI. Shipment and Storage:

- ExoPure™ Isolation Kit is shipped at room temp. Keep all the bottles upright in a cool and dark place for up to 12 months. DO NOT FREEZE!

VII. Reagent Preparation and Storage Conditions:




- Cap all the bottles well immediately after each use, to prevent evaporation.
- The maximum medium volume of each reaction is 20 mL from at most 1×10^7 cells. Do not exceed the suggested sample volume or the cell number. Otherwise it may cause indistinct layer separation and column clogging.
- One ExoPure™ Column is only for one reaction.

VIII. ExoPure™ Assay Protocol:

Sample Preparation for Exosome Isolation:

Fetal bovine serum (FBS) contains high level exosomes which may contaminate the cell derived exosomes. Use serum-free conditioned media to **starve stem cells for 48 hr before media harvest**.


1. Collect 20 mL cell culture medium.
2. Centrifuge the cell media at 3,000g for 15 min at 4°C to remove cells & debris. *Imp: skipping this step may cause filter clog in step 14.*
3. Transfer 20 mL clear supernatant (cell-free culture media) to a new 50 mL centrifuge tube (tube 1) and keep it on ice. **The maximum medium volume of each reaction is 20 mL from at most 1×10^7 cells. Do not exceed the suggested sample volume or the cell number. Otherwise it may cause indistinct layer separation and column clogging. One Column can be used for only one reaction.*
4. In another 50 mL centrifuge tube (tube 2), add the solutions in the following order to prepare A/B/C mixture (always prepare A/B/C mixture right before use): 1st Solution A (0.75 mL); 2nd Solution B (0.75 mL); 3rd Solution C (3 mL). **Cap Solution A, B & C bottles immediately after each use.*
5. Vortex the tube 2 (4.5 mL A/B/C mixture) for 5–10 sec to obtain a homogenous solution.
6. Add the 4.5 mL A/B/C mixture from tube 2 to tube 1 (20 mL cell-free culture media).
7. Tightly cap tube 1, vigorously vortex for 30 sec, and then incubate at 4°C for 30 min.
- 8a. The mixture now appears as 3 layers: Top layer, medium color; Bottom layer, colorless; Middle fluffy layer (exosome enriched layer).



Transmission electron microscopy (TEM) micrography of exosomes recovered from NIH3T3 cell media with ExoPure™ Isolation kit. Isolated exosomes show spherical and membrane encapsulated particles with the diameters varying between 20-200 nm. Homogeneous Spherical Exosome Isolated by ExoPure™ Isolation Kit (EM Analysis).



Representative dynamic light scattering (DLS, 632.8 nm laser) number distribution measurement of isolated exosome population from NIH3T3 cells (1x10⁶ cells) demonstrates a single peak (~150 nm) diameter. Homogenous Size Distribution of Exosomes Isolated by ExoPure™ Kit.



Products/Catalog Number
ExoPure™ Isolation Kit (Cell Media) # K1237-2, -10
ExoPure™ Isolation Kit (Serum, Plasma) # K1238-2, -10
ExoPure™ Isolation Kit (Stem Cell Media) # K1239-2, -10
ExoPure™ Isolation Kit (Urine) # K1240-2, -10
ExoPure™ Isolation Kit (Bio Fluids) # K1241-2, -10

- Our kit cannot isolate vesicles bigger than 300 nm.
- K1239 kit can be used to isolate exosome from stem cell media but not for isolating microparticles from cells.

