



### Concentrator improves the quality of circulating and EV-associated DNA, suitable for dPCR and NGS:

Circulating and EVs associated DNA was extracted from different samples, including human plasma, serum and four cell culture media using both ExoDNA™ (EXO-DNA) and ExoDNA™ Conc (EXO-DNAc), optimized with the concentrator columns. DNA yields were finally analyzed by Agilent Bioanalyzer and by Nanodrop



ExoDNA™ Conc Kit (EXO-DNAc) is able to concentrate around 4X the DNA yield extracted with the standard protocol (EXO-DNA). DNA yield measured by Nanodrop.

ExoDNA™ Conc Kit (EXO-DNAc) improves the DNA purity, compared to the standard protocol (EXO-DNA). DNA purity is expressed by 260/280 ratio.



Detection of the mutation BRAFV600E by digital PCR in genomic DNA isolated from serum of a metastatic prostate cancer patient treated with abiraterone. Circulating and EVs associated DNA was extracted by ExoDNA™ Conc Kit (EXO-DNAc) and a competitor kit. Digital PCR analysis revealed mutations in DNA extracted with (EXO-DNAc), whereas no mutation was detectable in DNA extracted with competitor kit. (Validated by Exosomes Siena).

### XI. Related Products:

Product Name	Cat. No.
ExoDNAPS™ Circulating and Exosome associated DNA (plasma/serum)	*K1230
ExoDNAPS™ Conc Exosome DNA Extraction Kit (plasma/serum)	*K1243
ExoDNAUC™ Circulating and Exosome associated DNA (urine/cell media)	*K1231
ExoDNAUC™ Conc Exosome DNA Extraction Kit (urine/cell media)	*K1244
ExoDNA™ Extraction Kit	K1242

\*Not available for sale in USA

### XII. Trouble Shooting:

Problems/Cause	Solution
<b>Low DNA Yield</b>	<ul style="list-style-type: none"> <li>. Be sure to add Proteinase K in the mixture of lysis buffer</li> <li>. Increase the incubation at RT during the elution step</li> <li>. Do not use water to elute DNA but use only the Elution buffer provided in the kit</li> <li>. Avoid to mix the lysate too vigorously</li> <li>. Avoid to form bubbles during mixing steps</li> </ul>
<b>DNA is Sheared or Degraded</b>	<ul style="list-style-type: none"> <li>. Do not touch the membrane of the column with the tip</li> <li>. Treatment with DNase must be done before to lyse the vesicles. Be careful to deactivate the DNase before to proceed to the lysis</li> <li>. Avoid repeated freeze and thaw cycles</li> </ul>
<b>Incomplete Elution</b>	<ul style="list-style-type: none"> <li>. Prolong the incubation time with Elution Buffer to 5-10 min or repeat elution step once again (25 μl + 25 μl)</li> </ul>
<b>Ethanol Contamination</b>	<ul style="list-style-type: none"> <li>. After the second washing step, centrifuge once again for 2 min at 15,000g. Dry the membrane of the column by incubation at RT (no flow hood)</li> </ul>