

- Discard the flow through then place the VB Column back in the 2 ml Collection Tube.
- Add 600 µl of Wash Buffer (make sure ethanol was added) to the VB Column.
- Centrifuge at 14-16,000g for 30 sec.
- Discard the flow-through and place the VB Column back in the 2 ml Collection Tube.
- Centrifuge at 14-16,000g for 3 min to dry the column matrix.

Step 4. Nucleic Acid Elution:

- Place the dried VB Column in a clean 1.5 ml microcentrifuge tube.
- Add 50 µl of RNase-free Water to the CENTER of the VB Column matrix.
- Let stand for at least 3 min to ensure the RNase-free Water is absorbed by the matrix.
- Centrifuge at 14-16,000g for 1 min to elute the purified nucleic acid.

X. Sensitivity:




Figure 1. Virus RNA was purified from 10E1-10E4 copy number of Red Spotted Grouper Nervous Necrosis Virus (RGNNV) using the BV Virus DNA/RNA Kit II (3 replications of each copy number). The purified RNA was eluted with 30 µl RNase-free Water. cDNA synthesis was carried out with a 10 µl aliquot of purified RNA using a Transcriptor First Strand cDNA Synthesis Kit (Roche) in a final volume of 20 µl. A Real time PCR assay was then performed with 3 µl of synthesized cDNA as template, primers (designed to amplify the T4 region on the RNA2 segment), and Fast SYBR Green PCR Master Mix using the StepOnePlus™ Real-Time PCR system (Applied Biosystems). The results confirmed that virus RNA can be successfully extracted and detected from as low as 10E1 copy number of RGNNV. The average cycle threshold (C_t): 10E4 = 23.88, 10E3 = 27.72, 10E2 = 31.22, 10E1 = 34.62. The low C_t values indicate a high number of target nucleic acid in the sample.




Figure 2. Hepatitis A Virus (HAV) RNA was extracted using the Viral Nucleic Acid Extraction Kit II. The purified RNA was analyzed by electrophoresis on a 1% agarose gel.

- M: BV 1 Kb DNA Ladder
1: HAV from original serum
2: HAV from 10X diluted serum
3: HAV from 100X diluted serum




Figure 3. HBV (DNA), HCV (RNA), HIV (RNA), and HTLV (RNA) were purified from 200 µl of positive clinical serum samples using the Viral Nucleic Acid Extraction Kit II. Real time qPCR and 1-step QRT PCR reactions were then conducted using the ABI 7300 Sequence Detection System (3 replications of each copy number). Serum samples containing various amounts of DNA/RNA viruses ranging from 10E1 to 10E6 copies/ml were successfully detected and identified. The low C_t values indicate a high number of target nucleic acid in the sample.

XI. Related Products:

Product Name	Cat. No.	Size
Magnetic Beads Viral DNA/RNA Kit	K1409-48, -96	-48, -96 rxns
Mag-Adenovirus Purification Kit	K1459-20, -100	20, 100 preps
Mag-Lentivirus and Retrovirus Purification Kit	K1458-20, -100	20, 100 preps
Viral DNA Extraction Kit	K1446-100	100 tests

XII. General Troubleshooting Guide:

Problem	Possible Reasons/Solution
Clogged Column	<ul style="list-style-type: none"> • Insufficient disruption and/or homogenization; Too much starting material • Centrifugation temperature was too low (should be 20°C-25°C)
Low RNA Yield	<ul style="list-style-type: none"> • Insufficient disruption and/or homogenization • DNA/RNA still bound to the VB Column membrane • Too much starting material; Ethanol carryover
RNA Degradation	<ul style="list-style-type: none"> • Harvested sample not immediately stabilized • Inappropriate handling of starting material; RNase contamination

FOR RESEARCH USE ONLY! Not to be used on humans.