

#### **Gentaur Europe BVBA** Voortstraat 49, 1910 Kampenhout BELGIUM Tel 0032 16 58 90 45 <u>info@gentaur.com</u>

01/19



# Fite's Staining Kit

(Cat# K1425-125; Store at RT)

#### I. Introduction:

The Fite's Staining Kit (for Leprosy and Nocardia) is intended for use in the histological visualization of mycobacterium Lepra bacillus (leprosy) and Nocardia. This kit may be used on formalin-fixed, paraffin-embedded or frozen sections.

Lepra bacillus: Red; Nocardia: Red; Background: Blue

#### II. Application:

- Histological applications
- For in vitro diagnostic use

# III. Sample Type:

- Formalin-fixed, paraffin-embedded (5 microns) or frozen sections.
- Control Tissue: Any well-fixed paraffin embedded Nocardia or Lepra bacillus infected tissue.

#### IV. Kit Contents:

••			
Components	K1425-125	Part Number	Storage Temperature
Xylene-Peanut Oil Solution	125 ml	K1425-125-1	RT
Carbol Fuchsin Solution	125 ml	K1425-125-2	RT
Acid Alcohol Solution (0.5%)	500 ml	K1425-125-3	RT
Methylene Blue Solution	125 ml	K1425-125-4	RT

# V. User Supplied Reagents and Equipment:

- · Distilled water
- · Coplin jars
- Xylene/Xylene Substitute
- Synthetic resin

# VI. Shipment and Storage:

All the reagents are shipped and stored at RT.

# VII. Reagent Preparation:

- Do not use if reagents become cloudy.
- Do not use past expiration date.
- Use caution when handling reagents.
- Non-Sterile.

# VIII. Lepra bacillus Procedure (Standard):

- 1. Deparaffinize sections in 2 changes of Xylene-Peanut Oil Solution for 12 minutes each.
- 2. Air dry slide for 15 min "without" removing oil film covering tissue section. Remaining film prevents de-staining of Lepra bacillus during differentiation.
- 3. Rinse slide in several changes of distilled water.
- 4. Incubate slide in Carbol Fuchsin Solution for 15 min.
- 5. Rinse slide in several changes of distilled water.
- 6. Differentiate section in Acid Alcohol Solution (1%) until background is pink.
- 7. Rinse slide in distilled water and check by microscope for correct differentiation.
- 8. Rinse in running tap water for 1 minute followed by 1 rinse in distilled water.
- 9. Dip slide 2-3 times in Methylene Blue Solution.
- 10. Dip slide quickly in distilled water and check by microscope for correct staining.
- 11. Air dry slide at room temperature.
- 12. Dip slide several times in Xylene or Xylene Substitute.
- 13. Mount in synthetic resin.

# Nocardia Procedure:

Preparation of Reagents Prior to Beginning:

Prepare Diluted Acid Alcohol Solution by mixing 25 ml of Acid Alcohol Solution (1%) with 25 ml of Distilled Water.

#### Procedure:

- 1. Deparaffinize sections in 2 changes of Xylene-Peanut Oil Solution for 12 minutes each.
- 2. Air dry slide for 15 minutes "without" removing oil film covering tissue section. Remaining film prevents de-staining of Lepra bacillus during differentiation.
- 3. Rinse slide in several changes of distilled water.
- 4. Incubate slide in Carbol Fuchsin Solution for 15 minutes.
- 5. Rinse slide in several changes of distilled water.
- 6. Dip slide 2-3 times in Diluted Acid Alcohol Solution.
- 7. Rinse slide in distilled water and check by microscope for correct differentiation. Avoid decolorizing the Nocardia organism.
- 8. Rinse in running tap water for 1 minute followed by 1 rinse in distilled water.
- 9. Dip slide 2-3 times in Methylene Blue Solution.
- 10. Dip slide quickly in distilled water and check by microscope for correct staining.
- 11. Air dry slide at room temperature.
- 12. Dip slide several times in Xylene or Xylene Substitute.
- 13. Mount in synthetic resin.



#### **Gentaur Europe BVBA** Voortstraat 49, 1910 Kampenhout BELGIUM Tel 0032 16 58 90 45 info@gentaur.com




# 01/19

# IX. Troubleshooting:

- Acid-fastness of the leprosy organisms is enhanced when the waxy capsule is protected by the mixture of xylene/peanut oil and the
  avoidance of dehydrating solutions.
- It is important to blot well; residual oil may produce staining artifact.
- If over-stained with methylene blue, organisms may be masked. Check microscopically before air drying. If over-stained, remove
  methylene blue with Acid Alcohol 1%; rinse thoroughly; repeat Methylene Blue Solution step with a shorter timing.
- If laboratory tap water is generally acidic, the methylene blue stain may be pale. Adjust staining time accordingly.
- A small percentage of Nocardia sp. organisms may resist taking the red stain and remain blue due to the growth phase of the individual organism.

# X. Data:



Staining using Fite's Staining Kit

FOR RESEARCH USE ONLY! Not to be used on humans.