




SARS-CoV-2 S1 Protein-ACE2 Binding Inhibitor Screening Kit

rev 10/20

(Catalog # K2050-100; 100 assays; Store at -20°C)

I. Introduction:

Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2), also known as the 2019 Novel Coronavirus (2019-nCoV) or human coronavirus 2019 (HCoV-19 or hCoV-19), is the cause of the Coronavirus Disease 2019 (COVID-19) pandemic. It is a RNA virus that causes severe respiratory diseases in humans. SARS-CoV-2 coronavirus contains four main structural proteins namely Spike (S), Membrane (M), Envelope (E), and Nucleocapsid (N) protein. Spike protein is located on the outer envelope of the virion and mediates the viral entry and thus, plays an important role in inducing neutralizing antibodies and protective immunity. S protein consists of S1 and S2 subunits. The S1 subunit contains a receptor-binding domain that can specifically bind to the host-receptor namely Angiotensin Converting Enzyme 2 (ACE2), which facilitates the entry of the virus into the target cells including respiratory or intestinal epithelial cells, endothelial cells, alveolar monocytes or macrophages. The receptor recognition step is an important determinant of the viral infectivity, pathogenesis & host range. Therefore, an intervention strategy that targets S1 protein and ACE2 interaction presents an important target for vaccination or antiviral strategies that includes small molecules and therapeutic antibodies. **BioVision's SARS-CoV-2 S1 Protein-ACE2 Binding Inhibitor Screening Kit** can be used to screen for potential inhibitors of S1 protein binding to human ACE2. In this assay, the binding of S1 protein to biotinylated human ACE2 is detected using Streptavidin-HRP. Subsequently, a TMB substrate is added to visualize the HRP enzymatic reaction thereby generating a blue colored product that changes to yellow once the stop solution is added. The density of the yellow color is proportional to the binding of S1 protein to Human ACE2. However, in the presence of potent inhibitor(s), the binding of S1 protein to Human ACE2 is suppressed thereby preventing the color generation. The assay kit is adapted to a 96-well format and provides a reliable test for high throughput screening of potential inhibitors of S1 protein binding to Human ACE2.



II. Application:

- Screening or characterizing inhibitors of SARS-CoV-2, S1 protein binding to human ACE2

III. Kit Contents:

Components	K2050-100	Cap Code	Part Number
S1 Protein coated Microplate	8 x 12 strips	---	K2050-100-1
Biotinylated Human ACE2	1 vial	Green	K2050-100-2
Streptavidin-HRP	25 µl	Amber	K2050-100-3
SARS-CoV-2 Inhibitor	70 µl	Red	K2050-100-4
TMB Substrate	20 ml	Amber/NM	K2050-100-5
Stop Solution	20 ml	NM	K2050-100-6
Wash Buffer (10X)	50 ml	NM	K2050-100-7
Assay Diluent	50 ml	Blue	K2050-100-8
Plate Sealers	4	---	K2050-100-9

IV. User Supplied Reagents and Equipment:

- Microplate reader capable of measuring absorbance at 650 nm and 450 nm
- Adjustable pipettes and pipette tips. Multichannel pipettes are recommended
- Deionized water
- Eppendorf tubes
- Absorbent paper

V. Storage Conditions and Reagent Preparation:

Store the kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening. Read the entire protocol before performing the assay. Upon opening, use within two months.

- S1 Protein coated Microplate:** Store at -20°C
- Biotinylated Human ACE2:** Reconstitute the vial in 60 µl Assay Diluent. Divide into aliquots & store at -20°C. Keep on ice while in use.
- Streptavidin-HRP:** After opening, store at 4°C, protected from light.
- TMB Substrate, Stop Solution and Assay Diluent:** After opening, store at 4°C. Bring to room temperature (RT) before use.
- SARS-CoV-2 Inhibitor:** Divide into aliquots and store at -20°C. Keep on ice while in use.
- Wash Buffer (10X):** Bring the bottle to RT. Prepare 1X Wash Buffer for the assay by diluting the 10X Wash Buffer with dH₂O. The 1X Wash Buffer can be stored at 4°C for one month.

VI. SARS-CoV-2 S1 Protein and Human ACE2 Binding Inhibitor Screening Protocol:

1. Screening Compounds, Inhibitor Control and Background Control Preparation:

Note: Wash the S1 Protein coated Microplate strip(s) one time with 1X Wash Buffer (250 µl/well) before adding the reagents.

Sample compound [S]: Dissolve the sample compound(s) at 50X or higher concentration in an appropriate solvent. Further dilute to 2X (the desired concentration) with Assay Diluent. Candidate antibodies can be prepared to proper dilution with Assay Diluent. Human serum can be diluted to at least 1:10 fold with Assay Diluent. Add 50 µl of diluted sample compound(s) or candidate samples into the designated wells of the S1 Protein pre-coated microplate(s).

Binding Control [Binding] (No Inhibitor): Add 50 µl of Assay Diluent to the designated well(s).

Incubate the plate by shaking gently at RT for 30 min, protected from light.

	[S]	[IC]	[Binding]	[BC]
Diluted Sample Compound(s)	50 µl	-	-	-
Diluted SARS-CoV-2 Inhibitor	-	50 µl	-	-
Assay Diluent	-	-	50 µl	100 µl

c)

Legend:
● Normal Human Serum
● COVID-19 Patient Serum

Y-axis: OD 450 nm (0 to 5)
X-axis: Serum Sample No. (0 to 40)

Serum Sample No.	Normal Human Serum (OD 450 nm)	COVID-19 Patient Serum (OD 450 nm)
1	2.8	1.1
2	2.8	0.6
3	2.7	0.4
4	2.7	1.2
5	2.6	0.1
6	2.7	0.1
7	2.7	0.4
8	2.8	0.1
9	2.8	0.4
10	2.8	0.4
11	2.8	0.1
12	2.8	0.4
13	2.8	0.1
14	2.8	0.1
15	2.8	0.1
16	2.8	0.1
17	2.8	0.1
18	2.8	0.1
19	2.8	0.1
20	2.8	0.1
21	2.8	0.1
22	2.8	0.1
23	2.8	0.1
24	2.8	0.1
25	2.8	0.1
26	2.8	0.1
27	2.8	0.1
28	2.8	0.1
29	2.8	0.1
30	2.8	0.1
31	2.8	0.1
32	2.8	0.1
33	2.8	0.1
34	2.8	0.1
35	2.8	0.1
36	2.8	0.1
37	2.8	0.1
38	2.8	0.1
39	2.8	0.1
40	2.8	0.1