



Human Platelet Isolation Kit

rev 01/21

(Catalog # K2071-25; 25 isolations; Store at -20 °C)

I. Introduction:

Platelets, also known as thrombocytes are small (2-4 μm in diameter), lens-shaped, anuclear cells found in the blood. They are produced in the bone marrow by the fragmentation of megakaryocytes. After fragmentation, they migrate into the circulation where they make up 4-7% of all blood cells. Like many other types of blood cells, platelets have multiple functions. Their primary role is in blood clotting. They also have an important role in innate immune response, inflammatory processes, and destruction of pathogens. Platelets can be found in two distinct forms namely unactivated and activated. The average lifespan of a platelet in blood is 3-10 days. **BioVision's human Platelet Isolation Kit** allows the *in vitro* isolation of intact, viable, unactivated platelets. The kit enables the high recovery of platelets of ($\geq 1 \times 10^7$ platelets/ml) yielding approximately 80% of the total platelets present in 1 ml of whole blood. The viability stain included in the kit is used to identify living platelets indicating that > 90% of the platelets in the isolated fraction is viable. The isolated fraction contains > 95% platelets and has minimal contamination of red blood cells and leukocytes.

II. Applications:

- Isolation of platelets from whole blood.
- Determination of platelet viability, purity, and quality.
- Studying platelet activation and characterization of platelet morphology.
- *In vitro* assays to evaluate primary platelet functions such as clotting, chemokine release, adhesion, chemotaxis etc.
- Studying platelet surface proteins.

III. Sample Types:

- Fresh whole blood collected < 8 hrs prior to platelet isolation, with EDTA or citrate anti-coagulants.

Notes:

1. For isolating healthy platelets, blood donor should not be taking aspirin for at least 48 hr or anti-platelet medications such as Plavix or Brilinta.
2. Platelets are fragile structures. Thus, it requires gentle handling and treatment to retain their *in vivo* properties and to prevent activation.

IV. Kit Contents:

Components	K2071-25	Cap Code	Part Number
Gradient Dilution Buffer	110 ml	NM	K2071-25-1
Density Gradient Media	25 ml	NM	K2071-25-2
Platelet Storage Media	50 ml x 2	NM	K2071-25-3
BSA Solution	5 ml	NM	K2071-25-4
Viability Stain	1 vial	Green	K2071-25-5

V. User Supplied Reagents and Equipment:

- Laminar flow hood to keep the reagents sterile
- DMSO
- PBS
- 15 ml conical tubes (polypropylene or polyethylene)
- Centrifuge with Swinging bucket Rotor
- Fluorescent Microscope with a dual FITC/TRITC Filter
- Hemocytometer
- Human blood 5 ml collected in EDTA or ACD anticoagulant
- Rocker platform if planning on storing platelets >24 hours
- Red Blood Cell Lysis Buffer (if desired)
- Multi-well plates, sterile

VI. Storage Conditions and Reagent Preparation:

Store kit at -20 °C, protected from light. Briefly centrifuge vials prior to opening. Read the entire protocol before performing the assay.

- **Gradient Dilution Buffer, Density Gradient Media, and Platelet Storage Media:** Use in a sterile environment to prevent contamination of isolated platelets. Bring to room temperature (RT) before use and mix well.
- **Viability Staining (lyophilized):** Reconstitute vial in 100 μl DMSO. **Light sensitive, do not expose to intense light.** Store at -20°C, protected from light.

VII. Platelet Isolation and Viability Assay Protocol (5 ml total volume of whole blood/isolation):

Work in a sterile environment. Use universal precautions when handling blood products and human body fluids. Handle blood and platelet samples carefully to avoid activation of the platelets.

This kit has sufficient reagents to isolate platelets from 5 ml samples.

1. Sample Preparation:

Prepare the following solutions to create the density gradient:

a. Density Barrier Solution (DBS, 1.072 g/ml). Add 5.0 volumes of Density Gradient Media (**DGM**) to 22 volumes of Gradient Dilution Buffer (**GDB**). For 5ml of DBS, add 0.925 ml DGM to 4.07 ml GDB. Mix well by inverting the mixture several times.

b. Blood: Gently invert 4-5 times to mix.

c. Platelet Storage Buffer: If you are using the buffer to store the isolated platelets, add 200 μl BSA to 0.8 ml of Platelet Storage Buffer.

Notes:

1. If either of the reagents is opened outside the hood, filter sterilization of the reagent is recommended.
2. If Platelet Storage Buffer is used for washing the platelets, BSA is not necessary.

