




# Elastin Assay Kit (Colorimetric)

02/21

(Catalog # K2077-100; 100 assays; Store at 4 °C)

## I. Introduction:

Elastin is a key protein of the extracellular matrix that provides resilience and elasticity to tissues and organs. It is highly elastic and is primarily present in extensible tissues including arteries, lung, ligaments and skin to resume their shape after stretching or contracting. Because of its critical role in normal development and function of vital organs, either an impairment of elastin synthesis or proteolytic degradation of the insoluble elastin fibers results in major clinical pathologies. Mutations in the Elastin gene may lead to diseases such as Williams-Beuren syndrome, Cutis laxa and Supravalvular aortic stenosis. Additionally, a number of skin diseases have been associated with abnormalities of elastin and marked fragmentation of elastin is observed in pulmonary emphysema. **BioVision's Elastin Assay Kit** is a quantitative dye-binding based method for the analysis of elastin levels in tissues and cultured cells. Due to its unique chemical composition and highly cross-linked nature, elastin is stable. But in this kit, elastin is denatured and solubilized under high heat and extreme pH conditions unlike other proteins. The isolated elastin then reacts with a dye to form a colored product, which is measured by absorbance at 513 nm. The colorimetric signal is directly proportional to the amount of elastin in samples. The assay is easy to perform, reliable and selective. It has a linear range of 5 µg-100 µg elastin per well.



## II. Application:

- Estimation of elastin levels in various sample types.

## III. Sample Types:

- Tissue lysate
- Cell lysate

## IV. Kit Contents:

Components	K2077-100	Cap Code	Part Number
Elastin Dye Buffer	35 ml	Amber/NM	K2077-100-1
Precipitating Reagent	3 ml	NM	K2077-100-2
Dye Dissociation Buffer	35 ml	Amber/WM	K2077-100-3
Oxalic Acid (1 M)	35 ml	NM	K2077-100-4
Elastin Standard (10 mg/ml)	2 x 1 ml	Yellow	K2077-100-5

## V. User Supplied Reagents and Equipment:

- Distilled Water
- PBS (optional)
- 96-well clear plate with flat bottom
- Multi-well spectrophotometer
- Mechanical shaker
- Metal heating block with the thermostat set between 95 °C and 100 °C

## VI. Storage Conditions and Reagent Preparation:

Upon arrival, store the kit at 4 °C, protected from light. Briefly centrifuge all small vials before opening. Read the entire protocol before performing the assay.

- **Elastin Dye Buffer and Dye Dissociation Buffer:** Store at 4 °C, protected from light. Bring to room temperature (RT) before use.
- **Precipitating Reagent:** Store at 4 °C. Keep on ice while in use.
- **Oxalic Acid (1 M):** Store at 4 °C. Bring to RT before use. Prepare 0.25 M Oxalic Acid for the assay by diluting 1 M Oxalic Acid with dH<sub>2</sub>O. **Note:** Oxalic Acid is corrosive, wear gloves when handling.
- **Elastin Standard (10 mg/ml):** Store at 4 °C. Use within two months.

## VII. Elastin Assay Protocol:

### 1. Sample Preparation:

- For tissues:** Rinse tissue samples with ice-cold dH<sub>2</sub>O or PBS to remove any residual blood and mince using clean scissors. Transfer the minced tissue (0.2-5 mg) to a microfuge tube and add 1 ml of ice-cold 0.25 M Oxalic Acid. Homogenize the tissue and incubate on a metal heating block at 100 °C for 1 hour. Vortex frequently and thoroughly. **For adherent cells:** Cultured cells can be removed using trypsin or any non-enzyme based cell dissociation solution. Wash the cells using PBS, centrifuge and aspirate the PBS. Resuspend the cell pellet in 0.5-1 ml of 0.25 M Oxalic acid per ~1 x 10<sup>6</sup> cells in a microfuge tube. Vortex thoroughly and incubate on a metal heating block at 100 °C for 1 hour. **Notes:** a) Do not tighten the tube caps while heating. b) The extract volumes of 0.25 M Oxalic Acid added need to be recorded for calculating the elastin content.
- Following 1 hour incubation, bring the tubes to RT. Centrifuge the homogenate at ≥ 10,000 x g and 4 °C for 15 min and transfer the acidic supernatant to a new microfuge tube. **The supernatant will be used for the assay.** Add 1 ml 0.25 M oxalic acid/or the same volume of 0.25 M oxalic acid added previously, to the **residual tissue/cell pellet** in the tubes and heat again at 100 °C for 1 hour and vortex thoroughly. Repeat the previous steps up to 2-3 times to completely solubilize the tissue elastin. **The extract from the last round should not contain any elastin.**
- Add 100 µl of supernatant into a separate 1.5 ml eppendorf tube labeled as **Sample**. Add 100 µl of 0.25 M Oxalic Acid into another 1.5 ml eppendorf tube labeled as **Background Control**.

