


# Nicotinamide N-methyltransferase Inhibitor Screening Kit II (F)

03/21

(Catalog # K2081-100; 100 assays; Store at -20 °C)

## I. Introduction:

Nicotinamide N-methyltransferase (NNMT; EC 2.1.1.1) is an important cytosolic enzyme. It catalyzes the N-methylation of nicotinamide by transferring a methyl group from S-adenosyl-L-methionine (SAM) to nicotinamide resulting in the formation of 1-methyl-nicotinamide (MNA) and S-adenosyl-L-homocysteine. NNMT plays a significant role in the regulation of metabolic pathways in tissues such as adipose and liver as well as in cancer cells. Abnormal NNMT expression and activity has been implicated in a number of chronic diseases including cardiovascular diseases, cancer, obesity, diabetes, osteoarthritis, Parkinson's disease etc. Therefore, small molecule inhibitors of NNMT could be beneficial in developing therapeutics for diseases characterized by abnormal NNMT activity. **BioVision's Nicotinamide N-methyltransferase inhibitor screening Kit II** is a non-coupled, plate-based fluorometric assay. The assay utilizes a selective NNMT substrate that is methylated by SAM to generate a fluorescent product, measured at Ex/Em = 320/420 nm. In the presence of 1-Methylnicotinamide, a NNMT inhibitor, the enzymatic activity is inhibited thereby resulting in a reduced fluorescent signal. The kit provides a one-step, rapid, reliable test for high-throughput screening of NNMT inhibitors.



## II. Application:

Screening/studying/characterizing potential inhibitors of NNMT

## III. Kit Contents:

Components	K2081-100	Cap Code	Part Number
NNMT Assay Buffer	25 ml	WM	K2081-100-1
DTT (1 M)	100 µl	Blue	K2081-100-2
NNMT Enzyme	50 µl	Green	K2081-100-3
S-Adenosylmethionine (SAM)	1 vial	Yellow	K2081-100-4
NNMT Substrate	35 µl	Amber	K2081-100-5
1-Methylnicotinamide	20 µl	Red	K2081-100-6

## IV. User Supplied Reagents and Equipment:

- 96-well white plate with flat bottom
- Multi-well fluorescence plate reader
- dH<sub>2</sub>O

## V. Storage Conditions and Reagent Preparation:

Store kit at -20 °C, protected from light. Briefly centrifuge all small vials prior to opening. Read the entire protocol before performing the assay.

- **NNMT Assay Buffer:** Store at 4 °C or -20 °C. Bring to room temperature (RT) before use.
- **DTT (1 M):** Store at -20 °C. **Prepare NNMT Assay Buffer containing 1 mM DTT immediately before use.** Add 2 µl of DTT (1 M) stock solution to 1998 µl of NNMT Assay Buffer.
- **NNMT Enzyme:** Divide into aliquots and store at -20 °C. Avoid repeated freeze/thaw cycles. Keep on ice while in use. Stable for two months.
- **S-Adenosylmethionine (SAM):** Reconstitute the vial in 40 µl dH<sub>2</sub>O. Pipette up and down to dissolve completely. Divide into aliquots and store at -20 °C. Avoid repeated freeze/thaw cycles. Keep on ice while in use. Use within two months.
- **NNMT Substrate:** Store at -20 °C, protected from light. Keep on ice while in use. Use within two months. **Note: NNMT Substrate has a strong odor. Wear gloves and mask when handling.**
- **1-Methylnicotinamide (MNA):** Store at -20 °C. Keep on ice while in use. Use within two months.

## VI. NNMT Inhibitor Screening Protocol:

### 1. Screening Compounds, Inhibitor Control & Blank Control Preparation:

**Note: Prepare NNMT Assay Buffer containing 1 mM DTT immediately before the assay.**

**Sample Compound:** Dissolve Test sample(s) in an appropriate solvent to make stock solution (see VI. 1. Note). Further dilute to 10 X using NNMT Assay Buffer containing 1 mM DTT. Add 10 µl of Diluted test sample(s) into wells of 96-well white plate designated as Sample [S].

**Enzyme Control (No Inhibitor) and Background Control:** Add 10 µl of NNMT Assay Buffer containing 1 mM DTT into designated well(s) of 96-well white plate designated as **Enzyme Control [EC]** and **Background Control [BC]** respectively.

**Inhibitor Control (1-Methylnicotinamide; MNA):** Prepare a 10-fold dilution of MNA by adding 2 µl of the MNA stock solution to 18 µl of NNMT Assay Buffer containing 1 mM DTT, mix well. Add 10 µl of diluted MNA into designated well(s) of 96-well white plate designated as **Inhibitor Control [IC]**.



**Note:** Various organic solvents may reduce the NNMT enzymatic activity (see Fig b). Prepare the stock samples in dH<sub>2</sub>O, if possible. If samples have to be prepared in an organic solvent, we recommend dissolving the test samples to 1000X or higher concentrations. Further dilute the test samples(s) with NNMT Assay Buffer containing 1 mM DTT to minimize the effect of organic solvent(s). Prepare parallel well(s) with the same final concentration of organic solvent as **Solvent Control [SC]** to test the effect of the solvent on NNMT activity. **It is strongly recommended that if you use any organic solvent for preparing the test inhibitors you have to use the Solvent Control well signal instead of the Enzyme Control well signal.**

	[S]	[IC]	[EC]	[BC]	[SC]
Test Sample	10 µl	-	-	-	-
Diluted MNA	-	10 µl	-	-	-
Assay Buffer containing 1 mM DTT	-	-	10 µl	10 µl	-
Solvent Control	-	-	-	-	10 µl

**2. NNMT Enzyme Preparation:** Prepare 60-fold dilution of NNMT Enzyme in NNMT Assay Buffer containing 1 mM DTT. For example, add 2 µl of NNMT Enzyme with 118 µl NNMT Assay Buffer /1 mM DTT), mix well. Mix enough reagents for the number of assays to be performed.

	<u>Enzyme Mix</u>	<u>Background Mix</u>
Diluted NNMT	30 µl	-
NNMT Assay Buffer containing 1 mM DTT	40 µl	70 µl

Add 70 µl of Enzyme Mix to test sample(s) [S], Inhibitor Control [IC], Enzyme Control [EC] and Solvent Control [SC] wells and mix well. Add 70 µl of Background Mix to the Background Control [BC] well. **Incubate the plate at 37 °C for 15~20 min, protected from light.**

**3. NNMT Reaction Mix:** Prepare 50-fold dilution of SAM. For example, dilute 2 µl of SAM stock solution with 98 µl NNMT Assay Buffer containing 1 mM DTT. Prepare 50-fold dilution of NNMT substrate. For example, dilute 2 µl of NNMT substrate stock solution with 98 µl NNMT Assay Buffer containing 1 mM DTT. Mix enough reagents for the number of assays to be performed. For each well, prepare a total 20 µl NNMT Reaction Mix as follows:

<u>NNMT Reaction Mix</u>	
Diluted NNMT Substrate	10 µl
Diluted SAM	10 µl


Mix well and add 20 µl NNMT Reaction to all wells including test sample(s), Inhibitor Control, [EC], [SC] and [BC] wells. Mix well. **The total final reaction volume for each well will be 100 µl.** **Note:** Prepare NNMT Reaction Mix immediately before adding to the wells.

**4. Measurement:** Immediately measure the fluorescence (Ex/Em = 320/420 nm) in kinetic mode at 37 °C for 30~60 min using a fluorometric microtiter plate reader. Choose two time points (t<sub>1</sub> & t<sub>2</sub>) in the linear range of the plot and obtain the corresponding RFU for Sample (R<sub>S1</sub> and R<sub>S2</sub>) and Background Control (R<sub>B1</sub> and R<sub>B2</sub>). **Note:** The Enzyme progressive curve is hyperbolic, with an initial linear portion followed by progressively slower reaction. Use the initial portion to check the linear range of the reaction.

**5. Calculation:** Calculate the slope for all test samples [S], Enzyme Control [EC], Solvent Control [SC] and Background Controls [BC] by dividing the net ΔRFU (R<sub>t2</sub>-R<sub>t1</sub>) values over reaction time Δt (t<sub>2</sub>-t<sub>1</sub>). Subtract the Slope of Background Control values from [S], [EC] and [SC]. **If [SC] slope is significantly different when compared to [EC], use [SC] values to determine effect of tested compound.**

$$\% \text{ Inhibition} = \frac{\Delta \text{RFU of EC} - \Delta \text{RFU of S}}{\Delta \text{RFU of EC}} \times 100$$

$$\% \text{ Relative Activity} = \frac{\Delta \text{RFU of S}}{\Delta \text{RFU of EC}} \times 100$$



**Figures.** a) NNMT enzyme activity in the presence or absence of MNA. b) NNMT enzyme activity in the presence or absence of DMSO. c) Inhibition of NNMT enzyme activity by MNA. IC<sub>50</sub> of MNA was calculated as 25.45±2.5 µM. Assays were performed following the kit protocol.

**VII. Related Products:**

NNMT Inhibitor Screening Assay Kit (K822)  
Nicotinamide N-Methyltransferase, Human Recombinant (7261)

S-Adenosyl-L-methionine disulfate tosylate (2077)

**FOR RESEARCH USE ONLY! Not to be used on humans.**