

2. Standard Curve Preparation: Prepare 1 mM Bile Acids Standard by adding 10 µl of 100 mM TBA Standard to 990 µl of ddH₂O. Further dilute to 12 µM by adding 12 µl of 1 mM Bile Acids Standard to 988 µl ddH₂O. Add 0, 2, 4, 6, 8, and 10 µl of 12 µM TBA Standard into a series of wells in a 96-well plate to generate 0, 24, 48, 72, 96 and 120 pmol of Bile Acids/well. Adjust the volume to 50 µl/well with ddH₂O.

Note:

a. The assay measures Enzymatic Activity Rates (Abs/min). For maximum accuracy, we recommend to carry out a Standard Curve at the same time samples are measured.

3. Probe Mix: Dilute TBA Probe 50-fold (i.e. 2 µl TBA Probe + 98 µl TBA Probe Buffer). Mix enough reagents for the total number of wells to be assayed. Mix & add 100 µl of Probe mix/well to Standards, sample background & sample wells. Incubate for 10 min at 37 °C.

4. Reaction Mix: Prepare 50 µl Reaction Mix for each well to be assayed as below and mix well.

	<u>Reaction Mix</u>	<u>*Background Control Mix</u>
TBA Cycling Assay Buffer	46 µl	48 µl
TBA Cycling Enzyme Mix	2 µl	----
NADH	2 µl	2 µl

Add 50 µl of Reaction Mix into Standard, and sample wells. Mix well. * For background correction, add Background Control Mix to sample background control well(s) and mix well.

5. Measurement: Measure absorbance at 405 nm in a kinetic mode at 37 °C for 60 min, protected from light. Choose two time points (T₁ & T₂) in the linear range to calculate the slope of every assayed well. Slopes for Standards, sample background control, and samples should be calculated using the same time points.

6. Calculation: Subtract 0 TBA Standard slope from all Standard readings. Plot the TBA Standard Curve. If sample background control slope is significant, then subtract sample background control slope from sample slope. Apply corrected slope to the Standard Curve to get B nmol TBA in the sample well.

$$\text{Sample TBA Concentration (C)} = B/V \times D \text{ nmol/}\mu\text{l or mM}$$

Where: **B** is amount of TBA in the sample well from the Standard Curve (nmol)

V is sample volume added into the reaction well (µl)


D is sample dilution factor

Note: For spiked samples, correct for any sample interference by using the following equation:

$$\text{TBA amount in spiked sample wells (B)} = \frac{\text{OD}_{\text{Sample slope (corrected)}}}{(\text{OD}_{\text{Spiked sample slope (corrected)}}) - (\text{TBA}_{\text{Sample slope (corrected)}})} \times \text{TBA spiked (nmol)}$$

1 mM ≡ 1000 µM

Bile Acids Molecular Weight: 521.69



Figures: (a) Total Bile Acids Standard Curve. **(b)** Estimation of Bile Acids concentration in human serum and urine. 30 µl of each undiluted sample was assayed following the kit protocol. Bile Acids concentrations are: Serum (in µM): A: 3.56 ± 0.41, B: 2.41 ± 0.27, C: 1.63 ± 0.17, D: 1.34 ± 0.25, Urine: 0.16 ± 0.02 µM/mM Creatinine.

VIII. Related Products:

- Alanine Aminotransferase (ALT or SGPT) Activity Colorimetric/Fluorometric Assay Kit (K752)
- Aspartate Aminotransferase (AST or SGOT) Activity Colorimetric Assay Kit (K753)
- Gamma Glutamyl Transferase (GGT) Activity Colorimetric Assay Kit (K784), Fluorometric Assay Kit (K785)
- Adenosine Deaminase Activity Assay Kit (Colorimetric) (K321), (Fluorometric) (K328)
- Adenosine Deaminase Activity Assay Kit Bilirubin (Total and Direct) Colorimetric Assay Kit (K553)
- Sodium deoxycholate (2830)
- Chenodeoxycholic acid (2831)
- Lithocholic acid (2187)

FOR RESEARCH USE ONLY! Not to be used on humans.