




Total Collagen Assay Kit (Colorimetric)

(Catalog # K218-100; 100 assays; Store at -20°C)

Rev 07/19

I. Introduction:

Collagen is the most abundant insoluble protein found in the extracellular matrix and connective tissues. It can be found in skin, tendons, bone, cartilage, muscle, vitreous humor and ligaments, among other tissues. There are more than sixteen - well characterized types of collagens, but types I, II and III collagen comprise more than 80% content in mammals. The triple-helical structure of collagen is quite unique: it consists of a repeating pattern of a basic trimer: Glycine-Proline-Hydroxyproline. In cells, collagens are secreted as procollagens and these chains are transported into the Endoplasmic Reticulum, where, numerous post-translational modifications lead to the formation of a triple helix with disulfide bonds. Excessive production of collagen is linked to pathological conditions including liver cirrhosis, lung fibrosis, and tumor growth. BioVision's Collagen Assay Kit is a simple and sensitive assay to detect small amounts of collagens in a variety of samples. The assay is based on the acid hydrolysis of samples to form hydrolysates and Hydroxyproline. This released Hydroxyproline gets oxidized to form a reaction intermediate, which further in the reaction, forms a chromophore (Abs 560 nm). The assay is simple, sensitive and specific for collagen and can detect as low as 0.5 µg of collagen in a variety of samples such as tissue homogenates, biological fluids and purified proteins.



II. Application:

- Measurement of collagen in various sample types.

III. Sample Types:

- Mammalian tissues
- Protein/peptide hydrolysates
- Serum
- Urine

IV. Kit Contents:

| Components | K218-100 | Cap Code | Part Number |
|--------------------------------------|----------|----------|-------------|
| Oxidation Buffer | 10 ml | WM | K218-100-1 |
| Chloramine T Concentrate | 0.6 ml | Red | K218-100-2 |
| Perchloric Acid/Isopropanol Solution | 5 ml | NM | K218-100-3 |
| DMAB Concentrate (in DMSO) | 5 ml | Amber | K218-100-4 |
| Collagen I Standard (2 mg/ml) | 0.1 ml | Yellow | K218-100-5 |

V. User Supplied Reagents and Equipment:

- 96-well clear plate with flat bottom
- Multi-well spectrophotometer
- 12 M Hydrochloric Acid (Concentrated HCl)
- For hydrolysis: Polypropylene Vials (BV Cat. No. M1352) and Screw Caps (BV Cat. No. M1353)

VI. Storage Conditions and Reagent Preparation:

Store the kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay.

- **Chloramine T Reagent:** For each well to be analyzed, add 6 µl of Chloramine T Concentrate to 94 µl of Oxidation Buffer and mix well.
- **DMAB Concentrate (in DMSO):** For each well to be analyzed, add 50 µl of DMAB Concentrate to 50 µl of Perchloric acid/Isopropanol Solution and mix well. Keep on ice, protected from light.

Note: The reagent concentrates are stable as supplied. Once the concentrates have been diluted to working concentration, they are only stable for 1-2 hr, so prepare fresh reagents as necessary for the number of Samples and Standards to be quantified.

VII. Collagen Assay Protocol:

1. Sample Preparation: Tissue or protein/peptide samples: tissue samples (*i.e.* lung) should be homogenized in ddH₂O, using 100 µl ddH₂O for every 10 mg of tissue. To a 100 µl of sample homogenate, add 100 µl concentrated HCl (~12 M, not provided) in a pressure-tight polypropylene screw-capped vial. Hydrolyze Samples at 120°C for 3 hours (**see note c**). **Urine:** hydrolyze Samples with equal volumes of 12 M concentrated HCl (*i.e.* 100 µl Urine + 100 µl HCl) in a pressure-tight polypropylene screw-capped vial. After homogenization, clarify Samples with activated charcoal by adding 4 mg of activated charcoal. Vortex and centrifuge at 10000 x *g* for 3 min to remove precipitate and activated charcoal. Repeat if needed. Transfer 2-30 µl of each hydrolyzed Sample to a 96-well plate and evaporate to dryness by heating the plate at 70°C on a hot plate/dry heat block or microplate incubator.

Notes:

- For Unknown Samples, we recommend performing a pilot experiment to ensure readings are within the standard curve range and adjusting the volume of Sample hydrolysate accordingly (2-30 µl of hydrolysate may be used) or diluting hydrolysate if necessary.
- For Samples with an extremely low collagen concentration, we recommend running two test Sample wells in parallel and spiking one with a known amount of Collagen I Standard (4.0 µg) to ensure accurate determination.
- For Sample hydrolysis, polypropylene vials with tight-fitting screw-on caps (without O-rings) yield best results. We recommend Biovision's Polypropylene Vials and Caps (Cat. No. M1353 and M1352).

