



**Note:**

a. The 100 µM 4-MU is stable for 2 months at -20 °C, protected from light.

- 3. Substrate Hydrolysis:** Mix 4 µl of the Lysozyme Substrate with 60 µl of Lysozyme Assay Buffer, vortex briefly and keep at RT. Add 10 µl of prepared Substrate to each well containing the Test Samples, Reagent Background Control and Lysozyme Positive Control. Mix well. Incubate the plate at 37 °C for 30-60 min, protected from light. After incubation time, add 50 µl Lysozyme Stop Buffer to each well containing Samples, Positive Control, Background Control and Standards. Mix well.

**Notes:**

a. For Samples exhibiting significant background, prepare a parallel Samples well(s) as Sample Background Control. Add 10 µl of Lysozyme Assay Buffer to the Sample Background Control well.

b. Incubation time depends on the Lysozyme enzymatic activity in Samples. Longer incubation time may be required for Samples having low Lysozyme activity.

- 4. Measurement:** Measure the fluorescence intensity at Ex/Em= 360/445 nm at 37 °C in end point setting using a fluorescence microtiter plate reader.

- 5. Calculation:** Subtract 0 Standard reading from all Standard readings. Plot the 4-MU Standard Curve. Subtract the Reagent Background Control from Sample reading(s). If the Sample Background Control reading is higher than the Reagent Background Control, subtract Sample Background Control instead. Apply Sample ΔRFU to 4-MU Standard Curve to obtain the corresponding pmol of product formed (**B**, in pmol). Calculate the activity of Lysozyme in the Sample as:

$$\text{Sample Lysozyme Activity} = \frac{B}{(\Delta T \times V)} \times D = \text{pmol/min/ml} = \mu\text{U/ml}$$

Where: **B** = 4-MU from the Standard Curve (pmol)


$\Delta T$  = Reaction time (min.)

**V** = Sample volume added into the reaction well (ml)

**D** = Dilution Factor

**Note:** For samples using lysozyme inhibitor (see VII. 1. Notes: b), subtract the Sample with inhibitor reading from Sample readings. Apply the Sample ΔRFU to 4-MU Standard Curve to obtain the corresponding pmol of product formed (**B**, in pmol), and calculate the activity of Lysozyme in the samples using the above mentioned formula. Lysozyme specific activity can be expressed as U/mg of protein

**Unit Definition:** One unit of lysozyme activity is the amount of enzyme that generates 1.0 µmol of 4-MU per min., at pH 5.0 at 37 °C.



**Figures:** (a) 4-Methylumbelliferon Standard Curve, results from three independent experiments. (b) and (c) Measurement of Lysozyme activity in human samples and cell cultured samples. Undiluted pooled serum (10 µl), pooled tears (15 µl), pooled saliva (30 µl) or J774 cells (15 µg), HepG2 cells (80 µg) were incubated with Lysozyme Substrate for 60 min. All assays were performed following kit protocols.

**VIII. RELATED PRODUCTS:**

Lysozyme Inhibitor Screening Kit (K220)

Lysozyme (M1237)

OryzaExp™ Lysozyme, Human Recombinant (P1400)

EZlys™ Lysozyme, Human (8005)

Lysozyme Antibody (A1110)

Lysozyme ELISA Kit (K4243)

**FOR RESEARCH USE ONLY! Not to be used on humans.**