



# ToxOut™ Protein G (Sephacrose) Spin Antibody Purification Kit 07/19

(Catalog # K2506-5 5 columns; Store at 4 °C)

## I. Introduction:

Protein G beads are widely used for IgG purification for their ability to bind selectively immunoglobulins. BioVision's Protein G (Cat. No. 6510) is a genetically engineered protein containing three Ig-binding regions of native Protein G. The cell wall binding region, albumin binding region and other non-specific regions have been eliminated from the recombinant Protein G to ensure the maximum specific IgG binding. Protein G-Sepharose beads (Cat# M1301) for IgG purification contain covalently coupled recombinant Protein G to 6% cross-linked Sepharose beads, the most popular resin for protein purification. BioVision's ToxOut™ Protein G (Sephacrose) Spin Antibody Purification Kit is a simple, ready to use kit for endotoxin-free IgG purification and immunoprecipitation. This kit can rapidly purify IgG from serum, ascites and cell culture media from various species such as human, mouse, rat, goat and rabbit.

## II. Applications:

Antibody purification for endotoxin-free samples

## III. Sample Types:

Serum, ascites and cell culture media

## IV. Kit Contents:

Components	K2506-5	Part Number
ToxOut™ Protein G (Sephacrose) Spin Column	5 columns	K2506-5-1
ToxOut™ Equilibration Buffer	1.8 ml	K2506-5-2
ToxOut™ Binding Buffer	10 ml x 2	K2506-5-3
ToxOut™ Elution Buffer	10 ml	K2506-5-4
ToxOut™ Neutralization Buffer	1.5 ml	K2506-5-5
ToxOut™ Collection Tube	50 tubes	K2506-5-6

## V. Specifications:

ToxOut™ Protein A (Sephacrose) Spin Column capacity is different in different species: ~2 mg/column

Species	IgG Binding Strength
Human	++++
Mouse	++++
Rat	+
Goat	+
Rabbit	++++

Note: Strong Binding ++++/ Weak Binding+

## VI. User Supplied Reagents and Equipment:

- Endotoxin-free 15 ml centrifuge tubes

## VII. Storage and Handling:

- Store at 4°C. All buffers expire 2 months after opening.
- Read the entire protocol before performing the assay.
- Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay.
- To prevent endotoxin contamination from dust, solution or dirty lab ware, only use endotoxin-free solutions and tubes and proceed with extra caution.
- Do not let the resin dry out anytime.

## VIII. Antibody Purification Kit Protocol:

1. **Sample Preparation:** Centrifuge samples at 10,000 x g and 4°C for 25 minutes and transfer supernatant to new tubes. Equilibrate samples by mixing with Equilibration Buffer at ratio of 9:1. (ex. mix 90 µl of sample with 10 µl Equilibration Buffer)

**Note:** Total IgG should be < 2 mg/column to ensure elution rate.

2. **Protein G Spin-Column Preparation:** Snap off the bottom plug from spin column by twisting it gently and save for later use. Put a ToxOut™ Collection Tube at the bottom to collect flow-through. Centrifuge the column at 700 x g for 2 min (use this step for all washes and eluates) to remove storage buffer. Discard flow-through. Wash and equilibrate the column 2 more times with 0.3 ml Binding Buffer each time.

3. **Sample Incubation:** Put the bottom plug back, load the equilibrated sample and plug the top cap back. Incubate the column for 1 hour at room temperature or overnight at 4°C by slowly inverting the column to achieve maximum binding.


4. **Washing:** Unplug both top cap and bottom plug and spin the column at 700 x g for 2 min to collect non-adsorbed material. Wash the column with 0.3 ml Binding Buffer and centrifuge at 700 x g for 2 min. Repeat this step 4 times.

**Note:** Reserve all flow through and washes in a 15 ml centrifuge tube until satisfactory enrichment of IgG elution is confirmed.




- 5. Elution:** Prepare 6 ToxOut™ collection tubes (labeled 1-6) by adding 20 µl Neutralization Buffer in each tube. Put the column inside tube #1 and add 0.1 ml Elution Buffer. Incubate the column for 1-2 min then centrifuge at 700 x g for 2 min. Mix the eluate with the Neutralization Buffer inside the tube immediately. Repeat elution step 3-5 times, each time in a new ToxOut™ collection tube.
- 6. Analyses:** Measure the IgG concentration by measuring OD absorbance at 280 nm. (1.4 OD<sub>280</sub> = 1 mg/ml IgG) Combine the eluted fractions containing the purified IgG.

a)



b)



**Figure: a) Binding capacity test:** IgG Elution curve using rabbit serum sample. Sample was incubated for 1 hour at room temperature and eluted with 0.1 ml Elution Buffer in each fraction. The IgG recovery rate is 90%. **b) ToxOut™ Protein G Sepharose endotoxin test:** Endotoxin level in concentration (right) and total amount (left) before and after purification of rabbit IgG with ToxOut™ Protein G (Sepharose) Spin Antibody Purification Kit. The recovered IgG eluate has > 98% reduction of endotoxin.

#### IX. RELATED PRODUCTS:

- ToxOut™ Rapid Endotoxin Removal Kit (K2501)
- ToxOut™ Protein A (Sepharose) Antibody Purification Kit (K2503)
- ToxOut™ Protein A (Sepharose) Column (M2503)
- ToxOut™ Endotoxin-Free Protein A Sepharose (M1300)
- ToxOut™ Endotoxin-Free Protein G Sepharose (M1301)
- Protein G (Cat# 6510)

**FOR RESEARCH USE ONLY! Not to be used on humans.**