



Immunoprecipitation (IP) Kit

rev. 4/18

(Catalog # K286-25; 25 assays; Store at 4 °C)

I. Introduction:

Immunoprecipitation (IP) is widely used in research and development, by which a protein can be selectively purified from samples. BioVision's immunoprecipitation kit provides optimized buffers for preparing cell/tissue extracts, antigen binding and washing steps. The protein A/G Sepharose beads provided in the kit has higher binding capacity with broader antibody isotype binding than traditional protein A or protein G resins. The kit can be used in variety of immunoprecipitation or Co-IP studies.

II. Application:

- Immunoprecipitation (IP) and Co-IP
- Functional study of Immunoprecipitated proteins/complexes
- SDS-PAGE or western blot analysis of Immunoprecipitated proteins/complexes

III. Sample Type:

- Tissue or cell extracts
- Biological samples

IV. Kit Contents:

Components	K286-25	Cap Code	Part Number
Non-Denaturing Lysis Buffer	40 ml	Clear NM	K286-25-1
RIPA Lysis Buffer	40 ml	Blue NM	K286-25-2
Protease Inhibitor Cocktail (lyophilized)	1 vial	Red	K286-25-3
10X Wash Buffer	20 ml	Clear WM	K286-25-4
Protein A/G Sepharose	1 ml	Orange	K286-25-5

V. User Supplied Reagents and Equipment:

- Primary antibody to the targeted protein
- Rotary mixer
- Phosphate Buffered Saline (PBS), 1M Tris/HCl pH 8.5

VI. Storage and Handling:

- Store kit at 4°C, protected from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the experiment.

VII. Reagent Preparation and Storage Conditions:

- **Non-Denaturing Lysis Buffer and RIPA Lysis Buffer:** Store buffer at 4°C once opened. Add 2 µl protease inhibitor cocktail per ml of Non-Denaturing Buffer just before use. Make fresh each time.
- **Protease Inhibitor Cocktail:** Re-suspend in 250 µl of DMSO. Aliquot and store at -20°C.
- **10X Wash Buffer:** To make 1X buffer add 1 ml of 10X buffer to 9 ml deionized water.
- **Protein A/G Sepharose:** Store at 4°C once the kit is opened.

VIII. Immunoprecipitation (IP) Protocol:

1. Sample Preparation:

Cell Extracts: For adherent cells, remove media and wash cells with PBS. Place culture plate on ice; add cold Lysis Buffer. Keep the plate on ice for one minute. Scrape the cells and gently transfer the disrupted cell suspension into a chilled microcentrifuge tube. Mix on a rotary mixer at 4°C for 30 min. Centrifuge at 10,000 g for 10 min. at 4°C & transfer the cell extract to chilled fresh tubes. For Suspension Cells, collect cells by centrifugation. Wash cells with PBS at room temperature and collect cells again by centrifugation. Drain the PBS carefully and prepare the cell lysates as described for adherent cells.

Notes:

- The no. of cells needed for optimal immunoprecipitation depends on the concentration of target antigen present in the sample and the affinity of the antibody to the antigen.
- Use Non-Denaturing Lysis Buffer for maintaining protein activity, studying protein-protein interaction & for antigens that are detergent soluble and can be recognized in the native form by the antibody. This buffer can be used for IP and Co-IP. The RIPA buffer, more denaturing than the Non-Denaturing Buffer, has 0.1% SDS, 1% NP40 and 0.5% Deoxycholate and can be used for an IP and may work for a Co-IP depending on how tight the complexes are.
- Lysis Buffer guidelines: 100-200 µl/well (24-well plate), 250-400 µl/well (6-well plate), 250-500 µl (100 x 60 mm dish) or 500-1000 µl (100 x 100 mm) dish.
- Butt end of a pipette tip can be used to scrape cells from wells.

Tissue Extracts: Snap Freeze dissected tissue and immediately grind to a fine powder using a mortar and pestle in a liquid nitrogen bath. Transfer the ground tissue to a pre-weighed chilled tube. Weigh the powder and store at -80°C until use. Add 300 µl Lysis Buffer with protease inhibitors per 5 mg of tissue powder. Mix on a rocker at 4°C for about an hour. Pass the lysate through a 25 gauge needle 3X. Collect the lysate and centrifuge at high speed (10,000g) at 4°C for 5 min. to remove cell debris. Transfer the tissue extract (supernatant) to a fresh tube.

