




Adenylate Kinase (AK) Activity Assay Kit (Colorimetric/Fluorometric)

(Catalog # K350-100; 100 assays; Store at -20°C)

2/17

I. Introduction:

Adenylate Kinase (AK) (EC 2.7.4.3) is an abundant enzyme involved in energy metabolism and homeostasis of cellular adenine nucleotide ratios in different intracellular compartments. The enzyme is found in the nucleus, cytosol, or mitochondria (intermembrane space or matrix) of various kinds of tissues. Adenylate kinase acts on two molecules of ADP to generate ATP and AMP. Nine isoforms of adenylate kinase have been identified. Erythrocyte adenylate kinase deficiency is associated with hemolytic anemia. Adenylate kinase also plays an important role in post-ischemic recovery and in apoptosis. BioVision's AK Activity Assay kit can kinetically measure Adenylate Kinase activity by detecting adenosine triphosphate (ATP) generated from adenosine diphosphate (ADP) as a substrate. ATP is detected via a multi-step reaction, resulting in the generation of an intermediate that reacts with the Adenylate Kinase Probe forming an end product that can be measured colorimetrically (OD 570 nm) or fluorometrically (Ex/Em = 535/587 nm).



II. Application:

- Detection of Adenylate Kinase activity

III. Sample Type:

- Purified Adenylate Kinase
- Cell and tissue lysate
- Mitochondrial lysate

IV. Kit Contents:

Components	K350-100	Cap Code	Part Number
AK Assay Buffer	25 ml	WM	K350-100-1
AK Probe	200 µl	Red	K350-100-2
ADP Substrate	200 µl	Brown	K350-100-3
AK Convertor	1 vial	Blue	K350-100-4
AK Developer	1 vial	Green	K350-100-5
Positive Control (AK Enzyme)	1 vial	Clear	K350-100-6
ATP Standard (1 µmol)	1 vial	Yellow	K350-100-7

V. User Supplied Reagents and Equipment:

- 96-well clear (colorimetric) or black plate (fluorometric) with flat bottom.
- Microplate reader capable of absorbance or fluorescence detection
- Protease Inhibitor Cocktail (Cat. # K271 or equivalent)

VI. Storage Conditions and Reagent Preparation:-

Store kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay.

- AK Assay Buffer:** Bring to room temperature before use. Store at -20°C or 4°C.
- AK Convertor and AK Developer:** Reconstitute each with 220 µl AK Assay Buffer and mix gently by pipetting. Briefly centrifuge to collect the contents at the bottom of the tube. Aliquot and store at -20°C. Avoid repeated freeze/thaw.
- Positive Control (AK Enzyme):** Reconstitute with 55 µl deionized water. Store at -20°C. Avoid repeated freeze/thaw. Use within two months.
- ATP Standard:** Dissolve in 100 µl of distilled water to generate a 10 mM stock solution. Keep on ice while in use. Store at -20°C. Avoid repeated freeze/thaw.

VII. Adenylate Kinase Activity Assay Protocol:

1. Sample Preparation: Rinse tissue and transfer ~50 mg of fresh or frozen tissue (stored at -80°C) to a prechilled tube. Add 150 µl cold AK Assay Buffer containing protease inhibitor cocktail (not provided) and thoroughly homogenize tissue on ice using an electrical homogenizer. Transfer the tissue homogenate to a cold microfuge tube.

To prepare cell extract, add 150 µl cold Homogenization Buffer containing protease inhibitor cocktail (not provided) to 1-5 x 10⁶ fresh or frozen cells and pipette several times to disrupt the cells. Transfer cell homogenate including cell debris to a cold microfuge tube and agitate on a rotary shaker at 4°C for at least 15 min. Centrifuge the tissue or cell homogenate at 16,000 X g, 4°C for 10 min. Transfer the clarified supernatant to a fresh pre-chilled tube & store on ice. Use lysates immediately to assay Adenylate Kinase activity. Mitochondria can be isolated using BioVision's Mitochondria Isolation Kit (K288) and solubilized in AK Assay Buffer for 10 min. on ice prior to use.

Add 2-50 µl of cell/tissue homogenate, mitochondrial lysate or purified protein into 96-well plate. For colorimetric assay, use 2-5 µl Positive Control. For fluorometric assay, dilute Positive Control 5x in AK Assay Buffer just before use. Add 2-5 µl of diluted Positive Control for the assay. Make up the volume of samples and Positive Control to 50 µl/well with AK Assay Buffer. Add 50 µl AK Assay Buffer to one well as reagent background control.

Notes:

- For unknown samples, we suggest testing several doses to ensure the readings are within the Standard Curve range.

