





- c. For samples having high protein content, we recommend deproteinizing the samples (tissue or cell lysate or serum, e.g.) using 10 kDa Spin Column (Cat. # 1997 or equivalent). Add sample to the spin column, centrifuge at 10,000 X g, 4°C for 10 min. Collect the filtrate.
- d. In cases where high background remains, sample may be treated with Sample Clean-Up Mix by adding 2 µl Sample Clean-Up Mix to 100 µl Sample and incubation at 37°C for 30 min. Boil sample at 90°C for ten minutes and filter through a 0.2 µm filter before proceeding with Assay.
- e. To ensure accurate determination of Histamine in the test samples or for samples having low concentrations of Histamine, we recommend spiking samples with a known amount of Histamine Standard (e.g. 100 pmol) and running them in parallel with unspiked samples.
- 2. Standard Curve Preparation:** Prepare 1 mM Histamine Standard by adding 10 µl of 50 mM Histamine Standard to 490 µl of ddH<sub>2</sub>O. Further dilute the 1 mM Histamine Standard by adding 10 µl to 190 µl Assay Buffer to generate a 50 µM Histamine Standard. Add 0, 2, 4, 6, 8, and 10 µl of 50 µM Histamine Standard into a series of wells in a 96-well plate to generate 0, 100, 200, 300, 400, and 500 pmol (0.5 nmol) of Histamine/well. Adjust the volume to 50 µl/well with Histamine Assay Buffer.
- 3. Reaction Mix:** Mix enough reagents for the total number of wells to be assayed. For each well, prepare 50 µl of Reaction Mix containing:

	<u>Reaction Mix</u>	<u>*Background Control Mix</u>
Histamine Assay Buffer	45.6 µl	47.6 µl
Histamine Enzyme Mix	2.0 µl	----
Histamine Developer	2.0 µl	2.0 µl
Histamine Probe	0.4 µl	0.4 µl

Mix well. Add 50 µl of Reaction Mix into Standard and sample wells. Mix.

\* For samples having background, add Background Control Mix to background control well(s) and mix.

- 4. Measurement:** Incubate plate at 37°C for 30 min. Measure fluorescence at 535 nm excitation/587 nm emission in end point mode.
- 5. Calculation:** Subtract 0 Histamine Standard reading from all readings. Plot the Histamine Standard Curve. If sample background control is significant, then subtract sample background control reading from sample reading to obtain corrected absorbance. Apply corrected absorbance to Standard Curve to get B nmol Histamine in the sample well.

$$\text{Sample Histamine Concentration (C)} = \frac{B}{V} \times D \text{ pmol/}\mu\text{l or } \mu\text{M}$$

Where: **B** is amount of Histamine in the sample well from Standard Curve (pmol)  
**V** is sample volume added into the reaction well (µl)  
**D** is sample dilution factor


**Note:** For spiked samples, correct for any sample interference by using the following equation:

$$\text{Histamine amount in spike sample well (B)} = \frac{\text{RFU}_{\text{sample (corrected)}}}{(\text{RFU}_{\text{sample+Histamine Std (corrected)}}) - (\text{RFU}_{\text{sample (corrected)}})} * \text{Histamine spike (nmol)}$$


Histamine molecular weight: 111.15 g/mol

1 µM ≡ 111.15 ppb

a)



b)



**Figure:** (a) Histamine Standard Curve. (b) Estimation of Histamine concentration in soy sauce. Soy sauce was prepared according to above protocol. Histamine concentrations are: Regular Soy Sauce: 720 ± 100 µM, Low Sodium Soy Sauce: 1010 ± 20 µM.

#### VIII. Related Products:

Nitric Oxide Colorimetric Assay Kit (K262)

Nitric Oxide Fluorometric Assay Kit (K252)

Histidine Decarboxylase Antibody (3691)

Histamine Colorimetric Assay Kit (K506)

Choline/Acetylcholine Quantification Assay Kit (K625)

α-Ketoglutarate Colorimetric/Fluorometric Assay Kit (K677)

Total Polyamine Assay Kit (Fluorometric) (K475)

Diamine Oxidase Activity Assay Kit (K496) Dounce Tissue

Homogenizer (1997)

**FOR RESEARCH USE ONLY! Not to be used on humans.**