



EZCell[™] Phagocytosis Assay Kit (Green Zymosan)

(Catalog # K397-100; 100 assays; Store at 4°C)

I. Introduction:

Phagocytosis in mammals serves as an important first line defense mechanism against invading pathogens. It is also essential for continuous clearance of dying cells, tissue remodeling, and acquisition of nutrients for some cells. Phagocytosis is a specific form of endocytosis initiated by recognition and binding of foreign particles by cell surface receptors, followed by their engulfment, and formation of phagosomes. Maturing phagosomes transform to phagolysosomes which destroy the pathogen through enzymes and toxic peroxides. Zymosan prepared from yeast cell wall (*Saccharomyces cerevisiae*), and consisting of protein-carbohydrate complexes is frequently used as a pathogen in phagocytosis assays. BioVision's EZCell[™] Phagocytosis Assay Kit (Green Zymosan) utilizes pre-labeled Zymosan particles as a tool for rapid and accurate detection and quantification of *in vitro* phagocytosis by fluorescent microscope, spectrophotometer or flow cytometry. Our kit provides a robust screening system for activators and/or inhibitors of phagocytosis and Toll-like receptors ligands (TLR).

II. Applications:

- Rapid detection, quantification and validation of phagocytosis in convenient 96-well format
- · Tracking ligand internalization and screening for effectors of phagocytosis

III. Sample Type:

• Phagocytic cell culture: adherent or suspension cells capable of phagocytosis

IV. Kit Contents:

Components	K397-100	Cap Code	Part Number
Phagocytosis Assay Buffer	2 X 100 ml	NM	K397-100-1
Buffer Additive	2 X 1 ml	Blue	K397-100-2
Green Zymosan	500 µl	Green	K397-100-3
10X Quenching Solution	500 µl	Yellow	K397-100-4

V. User Supplied Reagents and Equipment:

- A 6-, 12-, 24-, or 96-well clear plates should be used for culturing the cells only, and measurements of fluorescence should be
 performed in opaque white plates. Alternatively, opaque black or white plates with clear bottoms can be used for culturing and
 measurements
- Adherent or suspension cells capable of phagocytosis and appropriate media (e.g., JM774 or U937).
- Stock solutions of effectors of interest (for example, cytochalasin D, inhibitor of actin cytoskeletal rearrangement)
- Multi-well spectrophotometer measuring excitation and emission at 490 and 520 nm, respectively.
- Fluorescent microscope (optional) for observation or flow cytometer equipped with laser capable of excitation at 488 nm.

VI. Storage Conditions and Reagent Preparation:

- Store the entire kit at 4°C protected from light. Read the protocol before performing the assay.
- **Phagocytosis Assay Buffer**: Upon arrival, combine one entire vial of Buffer Additive with one Phagocytosis Assay Buffer, mix well. Use sterile pipetting technique throughout the assay.
- Green Zymosan Suspension: Equilibrate the suspension to room temperature and vortex gently for 5 seconds before each use.
- Quenching Solution: Dilute the content of the vial into 4.5 ml of 1X Phagocytosis Assay Buffer.

VII. Phagocytosis Assay Protocol:

- 1. Preparation of control and experimental wells: Subculture cells capable of phagocytosis (adherent or suspension) in appropriate medium. Day before the experiment, determine viability (BioVision Cat. K305, K303) and re-suspend to the concentration of 1 5 x 10⁶ viable cells/ml. Seed 100 µl of cell culture per well omitting the negative control wells and incubate overnight at 37 °C, 5 % CO₂. Next day, change the media and proceed to the phagocytosis effector assay. Your experiment should always consist of parallel negative, positive and experimental wells respectively.
- 2. Phagocytosis effector assay: Add 100 μl of cell culture media containing your effector of interest (not provided in the kit) at desired concentration (e.g. 20 μM Cytochalasin D) to each of the experimental wells. Aliquot 100 μl of media to each of the positive and 200 μl media to each of the negative control wells respectively. Incubate for 1 hour at 37 °C, 5 % CO₂.
- **3. Phagocytosis of Zymosan:** Add 5 µl of Zymosan slurry to all the wells. Immediately transfer the plate back to the incubator for 2 3 hours. The incubation time may be adjusted according to your normal protocol.
- **4. Green Zymosan Standard Curve:** Add 0, 1, 2, 3 and 4 μl of Green Zymosan slurry into a series of wells in 96-well plate. Adjust the volume to 100 μl with Phagocytosis Assay Buffer. Mix well. Immediately measure fluorescence using plate reader at Ex/Em=490/520 nm respectively. Subtract 0 Standard reading from all the readings and plot the Standard Curve.
- 5. Detection of internalized particles: Harvest the cells from each well into a separate tube and pellet for 5 minutes at 400 X g. Carefully aspirate off the media and gently re-suspend the cell pellets in 300 µl of ice cold Phagocytosis Assay Buffer containing the effector of interest at the same concentration as in the assay media. Centrifuge for 5 minutes at 400 X g and repeat the washing step 3 more times. Finally, suspend the cells in 200 µl of ice cold Phagocytosis Assay Buffer and proceed to the preferred method of detection.
- 6. Detection: Cells can be analyzed by scanning or imaging of all experimental and control wells in the plate reader or directly under the fluorescent microscope at Ex/Em at 490/520 nm, respectively.





Optional: For plate reader and microscope detection, re-suspend the cell pellets in 50 µl of the diluted Quenching Solution and incubate for two minutes at room temperature. Centrifuge the tube for 5 minutes at 400 X g and carefully remove the Quenching solution. Suspend in 200 µl of ice cold Phagocytosis Assay Buffer,

- a. For microscope and plate reader: Transfer 100 µl of each sample into a separate well and record the fluorescence.
- b. For flow cytometry: In the flow cytometry compatible vessel mix 100 µl of cell suspension and 900 µl of the Phagocytosis Assay Buffer. Analyze immediately in the FL1 channel of flow cytometer equipped with laser capable of excitation at 488 nm
- 7. Calculation: To calculate the net phagocytosis substract the average RFU of the no-cell negative-control wells from all positive control and experimental wells. The phagocytosis response to the experimental effector (% Effect) can be expressed as follows:

% Effect = $\frac{\text{Net experimental phagocytosis} \times 100\%}{\text{Net positive control phagocytosis}}$

The following figures demonstrate typical results with the EZCell[™] Phagocytosis Assay Kit (Green Zymosan). These data is for reference only and should not be used to interpret actual results. Your data will depend on the cell type and tested compound.



Figure: J774 macrophages were seeded overnight

at 5 x 10⁵ of viable cells/well. The next day the cells were pretreated with 20 µM Cytohalasin D for 1 h at 37°C prior to addition of 5 µl of Zymosan particles. Phagocytosis was conducted for 2 hours and the amount of engulfed Zymosan was determined as described in the Assay Protocol. (I) Inhibition of phagocytosis. Panel A and B: images of non-treated cells. Panel C and D: treatment with Cytohalasin D. (II) Flow cytometry plot. Red line: untreated control cells; green line: macrophages with engulfed Zymosan particles; blue line: inhibition of phagocytosis by Cytochalasin D. (III) Zymosan Standard curve.

VIII. RELATED PRODUCTS:

EZViable™ Calcein AM Cell Viability Assay Kit (Fluorometric) (K305) VisionBlue™ Quick Cell Viability Fluorometric Assay Kit (K303) Cell-Mediated Cytotoxicity Fluorometric Assay Kit (7-AAD/CFSE) Live/Dead Cell Viability Assay Kit (for Mammalian Cells) (K502) LDH-Cytotoxicity Colorimetric Assay Kit (K311) PicoProbe™ LDH-Cytotoxicity Fluorometric Assay Kit (K314)

FOR RESEARCH USE ONLY! Not to be used on humans.