



EZCell™ Phagocytosis Assay Kit (Red Zymosan)

4/17

(Catalog # K398-100; 100 assays; Store at 4°C)

I. Introduction:

Phagocytosis in mammals serves as an important first line defense mechanism against invading pathogens. It is also essential for continuous clearance of dying cells, tissue remodeling, and acquisition of nutrients for some cells. Phagocytosis is a specific form of endocytosis initiated by recognition and binding of foreign particles by cell surface receptors, followed by their engulfment, and formation of phagosomes. Maturing phagosomes transform to phagolysosomes which destroy the pathogen through enzymes and toxic peroxides. Zymosan prepared from yeast cell wall (*Saccharomyces cerevisiae*), and consisting of protein-carbohydrate complexes is frequently used as a pathogen in phagocytosis assays. BioVision's EZCell™ Phagocytosis Assay Kit (Red Zymosan) utilizes pre-labeled Zymosan particles as a tool for rapid and accurate detection and quantification of *in vitro* phagocytosis by fluorescent microscope, spectrophotometer or flow cytometry. Our kit provides a robust screening system for activators and/or inhibitors of phagocytosis and Toll-like receptors ligands (TLR).

II. Applications:

- Rapid detection, quantification and validation of phagocytosis in convenient 96-well format
- Tracking ligand internalization and screening for effectors of phagocytosis

III. Sample Type:

- Phagocytic cell culture: adherent or suspension cells capable of phagocytosis.

IV. Kit Contents:

Components	K398-100	Cap Code	Part Number
Phagocytosis Assay Buffer	2 X 100 ml	NM	K398-100-1
Buffer Additive	2 X 1 ml	Blue	K398-100-2
Red Zymosan	500 µl	Red	K398-100-3
10X Quenching Solution	500 µl	Yellow	K398-100-4

V. User Supplied Reagents and Equipment:

- A 6-, 12-, 24-, or 96-well clear plates should be used only for cell culturing. The measurement of fluorescence should be performed in opaque plates with clear bottoms. Alternatively, sterile opaque plates with clear bottoms can be used for both, culturing and measurements.
- Adherent or suspension cells capable of phagocytosis and appropriate media (e.g., JM774 or U937).
- Stock solutions of effectors of interest (for example, cytochalasin D, inhibitor of actin cytoskeletal rearrangement)
- Multi-well spectrophotometer measuring excitation and emission at 540 and 570 nm, respectively.
- Fluorescent microscope (optional) for observation or flow cytometer equipped with laser capable of excitation at 550 nm.

VI. Storage Conditions and Reagent Preparation:

Store the entire kit at 4°C protected from light. Read the protocol before performing the assay.

- **Phagocytosis Assay Buffer:** Upon arrival, combine **one** entire vial of Buffer Additive with **one** Phagocytosis Assay Buffer, mix well. Use sterile pipetting technique throughout the assay.
- **Red Zymosan:** Before each use, equilibrate the suspension to room temperature and vortex gently for 5 seconds.
- **Quenching Solution:** Dilute the content of the vial into 4.5 ml of 1X Phagocytosis Assay Buffer.

VII. Phagocytosis Assay Protocol:

- Preparation of control and experimental wells:** Subculture cells capable of phagocytosis in appropriate medium. The day prior to the experiment obtain a culture of $1 - 5 \times 10^6$ viable cells/ml. Aliquot 100 µl of the cell culture per well omitting the negative control wells and incubate the plate overnight at 37 °C, 5 % CO₂. Next day, change the media and proceed to the phagocytosis effector assay. Your experiment should always consist of parallel negative, positive and experimental wells respectively.
- Phagocytosis effector assay:** Add 100 µl of cell culture media containing your effector of interest (not provided in the kit) at desired concentration (e.g. 20 µM Cytochalasin D) to each of the experimental wells. Aliquot 100 µl of media to each of the positive and 200 µl media to each of the negative control wells respectively. Incubate for 1 hour at 37 °C, 5 % CO₂.
- Phagocytosis of Red Zymosan:** Add 5 µl of Zymosan slurry to all the wells. Immediately transfer the plate back to the incubator for 2 - 3 hours. The incubation time may be adjusted according to your protocol.
- Red Zymosan Standard Curve:** Add 0, 1, 2, 3 and 4 µl of Red Zymosan slurry into a series of wells in 96-well plate. Adjust the volume to 100 µl with Phagocytosis Assay Buffer. Mix well. Immediately measure fluorescence using plate reader at Ex/Em 540/570 nm respectively. Subtract 0 Standard reading from all the readings and plot the Standard Curve.
- Sample preparation:** Harvest the cells by centrifugation for 5 minutes at 400 X g. Carefully aspirate off the media and gently resuspend the cell pellets in 300 µl of ice cold Phagocytosis Assay Buffer containing the effector of interest at the same concentration as in the assay media. Centrifuge for 5 minutes at 400 X g and repeat the washing step 3 more times. Finally, suspend the cells in 200 µl of ice cold Phagocytosis Assay Buffer and proceed to the preferred method of detection.
- Detection:** Cells can be analyzed by FACS, fluorescent microscopy or by scanning of all experimental and control wells in the plate reader at Ex/Em at 540/570 nm, respectively.

