



Apo-BrdU-IHC™ *In Situ* DNA Fragmentation Assay Kit

(Catalog #K403-50; 50 assays; Store at -20°C)

I. INTRODUCTION:

Internucleosomal DNA fragmentation is a hallmark of apoptosis in mammalian cells. BioVision's **Apo-BrdU-IHC™ Kit** is a two-color TUNEL (Terminal deoxynucleotide transferase dUTP Nick End Labeling) assay for labeling DNA breaks to detect apoptotic cells by immunohistochemistry. The kit contains positive/negative control slides for assessing reagent performance; reaction and blocking buffers for processing individual steps in the assay; proteinase K; terminal deoxynucleotidyl transferase enzyme (TdT), bromodeoxyuridine triphosphate (Br-dUTP), biotin labeled antiBrdU antibody for labeling DNA breaks, horseradish peroxidase streptavidin conjugate, DAB, H₂O₂/Urea tablets for color generation and methyl green solution for counter staining the cells.

II. KIT CONTENTS:

| Component | Color Code | Volume | Store Temp. | Part Number |
|---|-------------|----------|-------------|-------------|
| Control Slides-pos/neg | natural box | 2 ea | -20° C | 403-50-1 |
| Blocking Buffer | white cap | 22 ml | -20° C | 403-50-2 |
| H ₂ O ₂ /Urea Tablets | amber vial | 6 ea | -20° C | 403-50-3 |
| Proteinase K | pink cap | 0.11 ml | -20° C | 403-50-4 |
| DAB Tablets | amber vial | 6 ea | -20° C | 403-50-5 |
| TdT Enzymes | yellow cap | 0.041 ml | -20° C | 403-50-6 |
| Br-dUTP | violet cap | 0.44 ml | -20° C | 403-50-7 |
| 200X Conjugate | clear cap | 0.035 ml | -20° C | 403-50-8 |
| 5X Reaction Buffer | green cap | 1.75 ml | +4° C | 403-50-9 |
| Anti-BrdU-Biotin Antibody | orange cap | 0.275 ml | +4° C | 403-50-10 |
| Methyl Green | natural cap | 6 ml | room temp. | 403-50-11 |

III. GENERAL CONSIDERATIONS:

- The components of this kit are for **Research Use Only**. To avoid reagent loss, centrifuge vials before using.
- After initial defrosting, the 5X Reaction Buffer and Anti-BrdU-Biotin Antibody should be stored at 4° C, and Methyl Green should be stored at room temperature. **Do Not Refreeze.**
- The control slides contain a mixture of apoptotic and non-apoptotic cells, allowing visualization of both positive & negative labeling within the same microscope field.
- Incubation time for proteinase K, DNase I (not provided), and the end labeling of the DNA may need to be empirically determined for your particular cell type and slide preparation. Use this protocol as a starting guideline.

IV. ASSAY PROTOCOLS: Staining of Paraffin Embedded Tissue (PET)

PET-A. Deparaffinization & Rehydration

- Immerse slides in xylene for 5 minutes at room temperature. Repeat using fresh xylene for second 5 minutes incubation.
- Immerse slides in 100% ethanol for 5 minutes at room temperature. Repeat using fresh 100% ethanol for second 5 minutes.
- Immerse slides in 90% ethanol for 3 min, then 80% ethanol for 3 min, and then 70% ethanol for 3 minutes at room temperature.
- Immerse slides briefly into 1X PBS and carefully dry the glass slide around the specimen.

If processing the kit's control slides simultaneously with unknown samples, please refer to the CFS protocol (CSF-A).

**** At this point it may be helpful to encircle the specimen using a waxed pen or a hydrophobic marker. ****

Do not let tissue specimen dry out at any step! If necessary, cover or immerse specimen in 1X PBS to keep hydrated!

PET-B. Permeabilization, Inactivation of Endogenous Peroxidase, & Equilibration

- Dilute only enough **Proteinase K (pink cap)** needed 1:100 in 10 mM Tris pH 8. Cover the entire specimen with 100 µl proteinase K. Incubate at room temperature for 20 minutes. **DO NOT OVER INCUBATE.**
- Rinse slide with 1X PBS. Gently tap off excess liquid and carefully dry the glass slide around the specimen.
- Dilute 30% H₂O₂ 1:10 in methanol. Cover the entire specimen with 100 µl of 3% H₂O₂. Incubate at room temperature for 5 minutes. **DO NOT OVER INCUBATE.**
- Rinse slide with 1X PBS. Gently tap off excess liquid and carefully dry the glass slide around the specimen.
- Dilute only enough **5X Reaction Buffer (green cap)** as needed 1:5 with dH₂O. Cover the entire specimen with 100 µl of the 1X Reaction Buffer. Incubate at room temperature for 10 to 30 minutes while preparing the labeling reaction mixture below.



For research use only

PET-C. End Labeling Reaction & Detection

10. Prepare the **Complete Labeling Reaction Mixture** as follows (Note: Mix only enough DNA Labeling Solution to complete the number of assays prepared per session. The DNA Labeling Solution is active for approximately 24 hours.)

| DNA Labeling Solution | 1 assay | 5 assays | 10 assays |
|--------------------------------|--------------|---------------|---------------|
| 5X Reaction Buffer (green cap) | 10 µl | 50 µl | 100 µl |
| TdT Enzyme (yellow cap) | 0.75 µl | 3.75 µl | 7.5 µl |
| Br-dUTP (violet cap) | 8 µl | 40 µl | 80 µl |
| Distilled H ₂ O | 32.25 µl | 161.25 µl | 322.5 µl |
| Total Volume | 51 µl | 255 µl | 510 µl |

11. Carefully blot the 1X Reaction Buffer from the specimen, taking care not to touch the specimen. Immediately apply **50 µl** of Complete Labeling Reaction Mixture (prepared above) onto each specimen except for the control slides which require only **25 µl each** (Note: The use of a cover slip at this point assures even distribution of the reaction mixture and prevents evaporation during incubation).

12. Cover the specimen with a piece of Parafilm cut slightly larger than the specimen (HINT: Folding up one corner of the Parafilm cover slip will aid in its application and removal). Place slides in a humid chamber and incubate at 37° C for 1 to 1.5 hours.

NOTE: The DNA End Labeling Reaction can also be carried out at 22-24° C overnight for the control slides. For samples other than the control slides provided in the kit, incubation times at 37° C may need to be adjusted to longer or shorter periods depending on the characteristics of the tissue you used.

13. Remove Parafilm cover slip and rinse slide with PBS. Gently tap off excess liquid and carefully dry the glass around the specimen.

14. Cover the entire specimen with 100 µl of **Blocking Buffer (white cap)**. Incubate at room temperature for 10 minutes. Carefully blot the Blocking Buffer from the specimen, taking care not to touch the specimen.

15. Immediately cover specimen with 100 µl of **Antibody Solution** (prepared as described below).

| Antibody Solution | 1 assay | 5 assays | 10 assays |
|--|---------------|---------------|----------------|
| Anti-BrdU-Biotin Antibody (orange cap) | 5 µl | 25 µl | 50 µl |
| Blocking Buffer (white cap) | 95 µl | 475 µl | 950 µl |
| Total Volume | 100 µl | 500 µl | 1000 µl |

16. Incubate with the **Antibody Solution** in the dark for 1-1.5 hours at room temperature (Hint: Cover slides with aluminum foil).

17. Rinse slide in PBS. Gently tap off excess liquid and carefully dry the glass around the specimen. Cover the entire specimen with 100 µl of **Blocking Buffer (white cap)**.

18. Dilute only enough of the **200X Conjugate (black cap)** needed 1:200 in Blocking Buffer (white cap). Prepared as described below.

| Conjugate Solution | 1 assay | 5 assays | 10 assays |
|-----------------------------|---------|----------|-----------|
| 200X Conjugate (black cap) | 0.5 µl | 2.5 µl | 5.0 µl |
| Blocking Buffer (white cap) | 100 µl | 500 µl | 1000 µl |

19. Carefully blot the Blocking Buffer from the specimen, taking care not to touch the specimen. Immediately apply 100 µl of diluted conjugate to the specimen. Incubate at room temperature for 30 minutes.

20. **Five minutes** before concluding incubation prepare DAB solution by dissolving one tablet of **DAB (amber vial)** and one tablet of **H₂O₂/Urea** (amber vial) in one ml of tap **H₂O**. This yields enough DAB solution for 10 specimens. (Note: Tap H₂O may contain metal ions that enhance the DAB reaction. *DAB is highly carcinogenic and care should be taken when handling.*)

21. Rinse slides with 1X PBS. Gently tap off excess liquid and carefully dry the glass slide around the specimen. Cover the entire specimen with 100 µl of DAB solution. Incubate at room temperature for 15 minutes. Rinse slides with H₂O and blot.

PET-D. Counterstain

22. Immediately cover the entire specimen with 100 µl of **Methyl Green Counterstain (natural cap)** solution. Incubate at room temperature for 3 minutes. Press edge of the slide against an absorbent towel to draw off most of the counterstain and place in a coplin jar slide holder.

22. Dip slides 2 times briefly into 100% ethanol. Blot slides briefly on an absorbent towel. Repeat step 22 using fresh 100% ethanol.

23. Blot slides briefly on an absorbent towel. Dip slides into xylene (or xylene substitute). Wipe excess xylene from back of slide and around specimen.10. Mount a glass cover slip using a mounting media such as permount (r) over the specimen.

24. Mount a glass cover slip using a mounting media such as permount (r) over the specimen.



V. ASSAY PROTOCOLS:

Staining of Cell Preparations Fixed on Slides (CFS)

The following protocol describes the method for measuring apoptosis in the positive and negative control slides that are provided in the APO-BRDU-IHC™ kit. The same procedure should be employed for measuring apoptosis in your own slide specimens.

Important points to remember before starting this assay:

The cells must be fixed prior to performing this assay. DO NOT LET THE CELLS DRY OUT BETWEEN OR DURING ANY STEPS!

To avoid loss of cells from glass slides during washing steps, it is recommended that slides be dipped into a beaker of 1X PBS rather than rinsed with a wash bottle.

CFS-A. Cell Fixation, Rehydration, & Permeabilization

1. Pellet cells at 300xg for 5 minutes at 4 ° C. Remove media. Add enough 1-4% formaldehyde (in PBS pH 7.4) to the pelleted cells to create a cell density of 1×10^6 cells/ml and incubate 15 minutes at room temperature.
2. Centrifuge at 300xg for 5 minutes at room temperature and resuspend at the same density in 70% ethanol. An aliquot of fixed cells (100-300 μ l) can then be adhered to glass slides by directly placing the suspension onto the slide or by using a Cytospin. Slides precoated with poly-L-lysine may enhance cell adherence.
3. Immerse slides in 1X PBS for 10 minutes at room temperature. Carefully dry the glass around the specimen.
**** At this point it may be helpful to encircle the specimen using a waxed pen or a hydrophobic slide marker ****
4. Dilute **Proteinase K (pink cap)** 1:100 in Tris pH 8. Cover the entire specimen with 50 - 100 μ l of the diluted proteinase K. Incubate at room temperature for **5 minutes**. DO NOT OVER INCUBATE.
5. Dip slide 2-3 times into a beaker of 1X PBS. Gently tap off excess liquid and carefully dry the glass slide around the specimen.

All the remaining steps of staining cells fixed on slides are identical to those steps outlined in the previous section for staining of paraffin-embedded tissue sections:

INACTIVATION OF ENDOGENOUS PEROXIDASES & EQUILIBRATION (PET-B, from Step 7)

END LABELING & DETECTION (PET-C)

COUNTERSTAIN (PET-D)

VI. ASSAY PROTOCOLS

Staining of Tissue Cryosections (TCS)

Important points to remember before starting this assay:

*Fixation of cryopreserved tissue is **required** prior to performing this assay. DO NOT LET THE TISSUE DRY OUT BETWEEN OR DURING ANY STEPS! if necessary cover or immerse the slide in 1X PBS to keep hydrated. To avoid loss of tissue from glass slides during washing steps, it is recommended that slides be dipped 2-3 times into a beaker of 1X PBS rather than rinsed with a wash bottle.*

TCS-A. Tissue Fixation, Rehydration, & Permeabilization

1. Immerse slides in 4% formaldehyde (in PBS, pH 7.4) for 15 minutes at room temperature. Gently tap off excess liquid and carefully dry the glass slide around the specimen.
2. Immerse slides in 1X PBS for 15 minutes at room temperature. Carefully dry the glass slide around the specimen.
****At this point it may be helpful to encircle the specimen using a waxed pen or hydrophobic slide marker****
3. Dilute **proteinase K (pink cap)** 1:100 in 10 mM Tris, pH 8. Cover the entire specimen with 50-100 μ l of the diluted proteinase K solution. Incubate at room temperature for **10 minutes**. DO NOT OVER INCUBATE.
4. Dip slide 2 - 3 times into a beaker of 1X PBS. Gently tap off excess liquid and carefully dry the glass slide around the specimen.

All the remaining steps of staining tissue cryosections on slides are identical to those steps outlined for staining of paraffin-embedded tissue sections:

INACTIVATION OF ENDOGENOUS PEROXIDASES & EQUILIBRATION (PET-B, from Step 7)

END LABELING & DETECTION (PET-C)

COUNTERSTAIN (PET-D)

VII. Technical Tips and Frequently Asked Questions About the APO-BRDU-IHC™ Assay

1. High background on slides all the cells are brown. There are many possible explanations for this:
-Did the sample dry out at any time during the staining?



- Did the H₂O₂/metalloid solution evaporate on the sample in the inactivation of endogenous peroxidases step (PET-C)?
 - Is this a false positive due to improper fixation of the tissue?
 - Was there a period of time between removal of the tissue and fixation when apoptosis (or necrosis) could have occurred?
 - Is their endogenous peroxidase activity?
 - Or non-specific binding of strep-avidin HRP? Try staining a control with no TDT.
 - Some people suggest increasing the H₂O₂ concentration to 5%, but over incubation with H₂O₂ or proteinase K can also damage DNA and create a false brown background.
2. Cell fixation using a DNA crossing linking chemical fixative is an important step in analyzing apoptosis. Unfixed cells may lose smaller fragments of DNA that are not chemically fixed in place inside the cell during washing steps. The researcher may have to explore alternative fixation and permeabilization methods to fully exploit their systems.
 3. Wash, don't squirt your slides.

RELATED PRODUCTS:

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents

Cell Fractionation System

- Mitochondria/Cytosol Fractionation Kit
- Nuclear/Cytosol Fractionation Kit
- Membrane Protein Extraction Kit
- Cytosol/Particulate Rapid Separation Kit
- Mammalian Cell Extraction Kit
- FractionPREP Fractionation System

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- High Throughput Apoptosis/Cell Viability Assay Kits
- LDH-Cytotoxicity Assay Kit

Cell Damage & Repair

- HDAC Fluorometric & Colorimetric Assays & Drug Discovery Kits
- HAT Colorimetric Assay Kit & Reagents
- DNA Damage Quantification Kit
- Glutathione Fluorometric & Colorimetric Assay Kits
- GST Activity Assay Kit
- Nitric Oxide Fluorometric & Colorimetric Assay Kits

Signal Transduction

- Camp & cGMP Assay Kits
- Akt & JNK Activity Assay Kits
- Beta-Secretase Activity Assay Kit

Adipocyte & Lipid Transfer

- Recombinant Adiponectin, Survivin, & Leptin
- CETP Activity Assay & Drug Discovery Kits
- Total Cholesterol Quantification Kit

Molecular Biology & Reporter Assays

- siRNA Vectors
- Cloning Insert Quick Screening Kit
- Mitochondrial & Genomic DNA Isolation Kits
- 5 Minutes DNA Ligation Kit
- 20 Minutes Gel Staining/Destaining Kit

Antibodies & Recombinant Proteins (many)