



# Apo-BrdU-Red™ *In Situ* DNA Fragmentation Assay Kit

(Catalog #K404-60; 60 assays; Store Components Separately-see below)

## I. INTRODUCTION:

Internucleosomal DNA fragmentation is a hallmark of apoptosis in mammalian cells. BioVision's **Apo-BrdU-Red *In Situ* DNA Fragmentation Assay Kit** provides complete components including positive and negative control cells for conveniently detecting DNA fragmentation by fluorescence microscopy or flow cytometry. The kit utilizes Br-dUTP (brominated deoxyuridine triphosphate nucleotides) which is more readily incorporated into DNA strand breaks than other larger ligands (e.g., fluorescein, biotin or digoxigenin). The greater incorporation gives rise to brighter signal when the Br-dUTP sites are identified by a Red fluorescence labeled anti-BrdU monoclonal antibody. The assay is suitable for studying apoptosis with GFP transfected cells.

## II. KIT CONTENTS:

Components	Cap Color	Volume	Store Temp.	Part#
Positive Control Cells	brown	5 ml	-20°C	K404-60-1
Negative Control Cells	natural	5 ml	-20°C	K404-60-2
Wash Buffer	blue	120 ml	+4°C	K404-60-3
Reaction Buffer	green	0.6 ml	+4°C	K404-60-4
TdT Enzymes	yellow	45 µl	-20°C	K404-60-5
Br-dUTP	violet	0.48 ml	-20°C	K404-60-6
Rinse Buffer	red	120 ml	+4°C	K404-60-7
Anti-BrdU-Red Antibody	orange	0.3 ml	+4°C	K404-60-8
7-AAD/RNase Staining Buffer	amber bottle	30 ml	+4°C	K404-60-9

## III. APOBRDU ASSAY PROTOCOL FOR CULTURED CELLS:

### A. Cell Fixation

1. Induce apoptosis in cells by desired method. Concurrently incubate a control culture without induction.
2. Pellet  $1-5 \times 10^6$  cells and resuspend in 0.5 ml of PBS.
3. Fix the cells by adding 5 ml of 4% (w/v) formaldehyde in PBS and place on ice for 15 minutes.
4. Centrifuge the cells for 5 min at  $300 \times g$  and discard the supernatant.
5. Wash cells in 5 ml of PBS and pellet the cells by centrifugation. Repeat one time the wash and centrifugation step.
6. Resuspend the cells in 0.5 ml of PBS.
7. Add the cells to 5 ml of ice-cold 70% (v/v) ethanol. Let cells stand for a minimum of 30 min (or overnight if you prefer) on ice or in the freezer.
8. Store the cells in 70% (v/v) ethanol at  $-20^\circ\text{C}$  until use. Cells can be stored at  $-20^\circ\text{C}$  for several days before use.

### B. Detection by Flow Cytometry:

The procedures can be used for both control cells and your testing cells.

1. Resuspend the fixed cells by swirling the vials. Remove 1 ml aliquots of the cell suspension ( $\sim 1 \times 10^6$  cells per ml) and place in 12 x 75 mm tubes. Centrifuge ( $300 \times g$ ) for 5 min and carefully remove the ethanol by aspiration.
2. Resuspend the cells with 1 ml of **Wash Buffer** (blue cap). Centrifuge as before and remove supernatant carefully by aspiration.
3. Repeat one time the washing step (step 2).
4. Resuspend in 50 µl of the **DNA Labeling Solution** prepared as below:


DNA Labeling Solution	1 assay	10 assays
TdT Reaction Buffer (green cap)	10 µl	100 µl
TdT Enzyme (yellow cap)	0.75 µl	7.5 µl
Br-dUTP (violet cap)	8 µl	80 µl
ddH <sub>2</sub> O	32.25 µl	322.5 µl
Total Volume	51 µl	510 µl

5. Incubate the cells in the **DNA Labeling Solution** for 60 min at  $37^\circ\text{C}$ . Shake cells every 15 min to resuspend.
6. Add 1 ml of **Rinse Buffer** (red cap) to each tube and centrifuge for 5 min. Remove supernatant by aspiration.
7. Repeat one time the rinsing step (step 6).
8. Resuspend cells in 0.1 ml of the **Antibody Solution** prepared as below:

Antibody Solution	1 assay	10 assays
Anti-BrdU-Red Antibody (orange cap)	5 µl	50 µl
Rinse Buffer (red cap)	95 µl	950 µl



9. Incubate the cells with the **Antibody Solution** in the dark for 30 min at room temperature.
10. Add 0.5 ml of **7-AAD/RNase A Solution** (amber bottle).  
Note: If DNA cell cycle information is not necessary, add 0.5 ml PBS instead and continue with step 12.
11. Incubate the cells in the dark for 30 min at room temperature.
12. Analyze the cells by flow cytometry (Ex/Em = 488/576 nm for BrdU-Red and 488/655 nm for 7-AAD). Cells should be analyzed within 3 hours of staining.



**RELATED PRODUCTS:**

Apoptosis Detection Kits & Reagents

- Annexin V Kits & Bulk Reagents
- Caspase Assay Kits & Reagents
- Mitochondrial Apoptosis Kits & Reagents
- Nuclear Apoptosis Kits & Reagents
- Apoptosis Inducers and Set
- Apoptotic Cell Isolation Kit
- Apoptosis siRNA Expression Vectors

Cell Proliferation & Senescence

- Quick Cell Proliferation Assay Kit
- Senescence Detection Kit
- LDH-Cytotoxicity Assay Kit
- Bioluminescence Cytotoxicity Assay Kit
- Live/Dead Cell Staining Kit

Metabolism Assays

- Ascorbic acid, Lactate, Pyruvate Assay Kits
- Glucose, sucrose, Galactose, Lactose, Maltose Assay Kits
- Choline, Ethanol, Uric acid, L-Amino acid assay Kits
- ADP/ATP, NAD/NADH, NADP/NADPH Quantitation Kits