



Alpha Galactosidase (α-Gal) Activity Assay Kit (Fluorometric)

04/18

(Catalog # K407-100; 100 assays; Store at -20°C)

I. Introduction:

Alpha-Galactosidase (α -Gal; EC 3.2.1.22) hydrolyzes alpha-galactosyl moieties found in glycolipids and glycoproteins. In mammals, α -Gal hydrolyzes poly- and oligosaccharides commonly found in dietary sources that are difficult to digest. Therefore, α -Gal is used in dietary supplements that help to reduce the production of intestinal gases due to consumption of certain foods. It is known total α -Gal activity is due to two major isozymes with unique, yet different thermostability profiles. Alpha-Galactosidase A, is thermolabile and represents approximately 90% of total α -Gal activity. Fabry Disease, a lysosomal disease disorder, is characterized by mutations in alpha-Galactosidase A. These mutations cause abnormal accumulation of glycosphingolipids in lysosomes. BioVision's Alpha Galactosidase A ctivity Assay Kit provides a simple, rapid way to monitor total α -Gal activity in wide variety of biological samples. In this kit, α -Gal cleaves a synthetic specific substrate releasing a fluorophore, which can be easily quantified (Ex/Em= 360/445 nm). The assay is specific, sensitive and can detect as low as 0.1 μ U of α -Galactosidase activity.

α-Gal Substrate Alpha Galactosidase

Cleaved Substrate + Fluorescent Product (Ex/Em= 360/445 nm)

II. Applications:

• Measurement of α-Galactosidase activity in various samples

III. Sample Type:

- Tissue Homogenates: kidney, etc.
- Cell Lysates: U937, etc.
- · Biological fluids: Saliva, etc.

IV. Kit Contents:

Components	K407-100	Cap Code	Part Number
α-Gal Assay Buffer	25 ml	NM	K407-100-1
α-Gal Stop Buffer	25 ml	WM	K407-100-2
α-Gal Substrate	220 µl	Blue	K407-100-3
4-Methylumbelliferone Standard	35 µl	Yellow	K407-100-4
α-Gal Positive Control	1 vial	Green	K407-100-5

V. User Supplied Reagents and Equipment:

- Multi-well spectrophotometer (ELISA reader)
- 96-well white plate with flat bottom is preferred for this assay. 96-well clear plate can also be used.
- Dounce Tissue Homogenizer (Cat. #1998)

VI. Storage Conditions and Reagent Preparation:

Store kit at -20°C, protect from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay. Upon opening, use within two months.

- α-Gal Assay Buffer and Stop Buffer: Store at 4 °C or -20 °C. Bring to 37 °C before use.
- α-Gal Substrate: Light sensitive. Thaw at room temperature. Store at -20 °C.
- 4-Methylumbelliferone Standard (5 mM): Light sensitive. Thaw at room temperature. Store at -20 °C.
- α-Gal Positive Control: Reconstitute with 20 µl α-Gal Assay Buffer and mix thoroughly. Store at -20 °C. Keep on ice while in use. Use within two months.

VII. α-Gal Activity Assay Protocol:

 Sample Preparation: For tissue and cells: Homogenize tissue (10 mg) or pelleted cells (~5 x 10⁵) with 100 µl ice-cold α-Gal Assay Buffer and keep on ice for 10 min. Centrifuge samples at 12,000 x g at 4 °C for 10 min. and collect the supernatant. Dilute the supernatant 10-20 fold in α-Gal Assay Buffer. Add 2-10 µl of diluted samples into a 96-well plate that will be designated as Sample(s).
 For biological fluids: Undiluted fluids can be added directly to the well. Add 2-10 µl of samples into well(s) in a 96-well plate that will be designated as Samples. For Reagent Background Control: add same volume of α-Gal Assay Buffer in parallel well(s). For Positive Control: dilute reconstituted α-Gal Positive Control 1:10 fold with α-Gal Assay Buffer prior to the assay and add 2-6 µl of diluted α-Gal Positive Control into desired wells(s). Adjust the volume of Positive Control, Sample(s), and Reagent Background Control to 40 µl/well with α-Gal Assay Buffer.

Note:

- a. We suggest using 3-5 different volumes of the samples per well to ensure the readings are within the standard curve range and the progress curve rates are within the linear range.
- b. Do not store unused diluted α-Gal Positive Control.
- 2. Standard Curve Preparation: Prepare a 100 μM 4-Methylumbelliferone (4-MU) Standard by adding 10 μl of 5 mM 4-MU to 490 μl α-Gal Assay Buffer in amber tube. Further dilute the 100 μM Standard solution 5-fold by adding 20 μl of 100 μM 4-MU to 80 μl α-Gal Assay Buffer to generate 20 μM 4-MU Standard. Add 0, 2, 4, 6, 8, 10 μl of 20 μM 4-MU standard into a series of wells to generate 0, 40, 80, 120, 160, 200 pmol/well of 4-MU Standard respectively. Adjust the volume to 60 μl/well with α-Gal Assay Buffer.





Note: Equilibrate the α-Gal Assay Buffer to 37 °C prior to the assay.

- 3. Substrate Hydrolysis: Prepare sufficient volume of 10-fold dilution of α-Gal Substrate (i.e. Dilute 4 µl of α-Gal stock Substrate with 36 µl of α-Gal Assay Buffer), vortex briefly. Add 20 µl of diluted α-Gal Substrate to each well containing the test Sample(s), Positive Control and Reagent Background Control. *The total volume in each well (i.e. Samples, Positive Control and Reagent Background Control)* should be 60 µl). Mix well and incubate at 37 °C for 2 hours, avoid light. After incubation, add 200 µl of α-Gal Stop Buffer to each well containing Sample(s), Positive Control, Reagent Background Control and Standards. Mix well.
 Note:
 - a. Equilibrate the α-Gal Stop Buffer to 37 °C prior to the assay.
- b. Standards can be prepared at the end of the incubation time, and measured in end-point mode.
- 4. Measurement: Measure fluorescence intensity (Ex/Em= 360/445 nm) at 37°C using an end-point setting.
- **5.** Calculation: Subtract 0 Standard reading from all Standard readings. Plot the 4-MU Standard Curve; subtract the Reagent Background Control reading from all Sample readings. Apply sample ΔRFU to 4-MU Standard Curve to obtain the corresponding pmol of product formed (**B**, in pmol) and calculate the activity of α-Galactosidase activity in the sample as:

Specific Sample α -Galactosidase Activity = B/ (2 x V x P) x D (pmol/h/mg \equiv 0.0167 μ U/mg)

- Where: **B** = 4-MU amount in sample well from Standard Curve (pmol)
 - **2** = Reaction time (hour)
 - V = Sample volume added into the reaction well (ml)
 - **P** = Initial Sample Concentration in mg-protein/ml (mgP/ml)
 - **D** = Sample Dilution Factor

1 pmol/h= 0.0167 pmol/min \equiv 0.0167 μ U

Unit Definition: One unit of α -Galactosidase activity is the amount of enzyme that generates 1.0 µmol of 4-Methylumbelliferone per min., at pH 4.5 at 37 °C.



Figure: (a) 4-Methylumbelliferon Standard Curve. Results are from multiple experiments. (b) α -Galactosidase Activity in Mouse Kidney Tissue Extracts (1 µg protein) and U937 Cell Lysates (0.2 µg protein). (c) Measurement of α -Galactosidase Activity in undiluted Human Pooled Saliva (5 µl). All assays were performed following kit protocol.

VIII. RELATED PRODUCTS:

β-Galactosidase Activity Assay Kit (K821)

β-Galactosidase Staining Kit (K802)

β-Galactosidase Inhibitor Screening Kit (K827)

α-L-Fucosidase Activity Assay Kit (Fluorometric) (K542)

α-L-Fucosidase (FUCA1) Assay Kit (Colorimetric)

β-glucuronidase Activity Assay Kit (K514)

β-N-Acetylglucosaminidase Activity Assay Kit (Colorimetric) (K733)

EZClick™ O-GlcNAc Modified Glycoprotein Assay Kit (FACS/Microscopy, Green Fluorescence) (K714)

EZClick™ Sialic Acid (ManAz) Modified Glycoprotein Assay Kit (FACS/Microscopy, Green Fluorescence) (K441)

EZClick™ O-GalNAc Modified Glycoprotein Assay Kit (FACS/Microscopy, Green Fluorescence)

Dounce Tissue Homogenizer (1998)

FOR RESEARCH USE ONLY! Not to be used on humans

.. ..

CL 05005 TTCL IT (400) 403 1000 T (400) 403 1001 I