


Fatty Acid Uptake Assay Kit (Cell-Based)

(Catalog # K408-100; 100 Reactions; Store at -20°C)

rev 06/19

I. Introduction:

Fatty acids are essential dietary nutrients consisting of a carboxylic acid moiety linked to a long aliphatic hydrocarbon chain. They are the fundamental building blocks of more complex lipids such as those that form cellular membranes. In the blood, long-chain fatty acids travel bound to serum albumin and are transported into cells by members of the SLC27 family of transmembrane transporter proteins. Uptake of long-chain free fatty acids into adipocytes is essential for regulating circulating lipid levels and cellular lipid metabolism. Abnormalities in adipocyte fatty acid uptake and lipid metabolism are major contributing factors to obesity, diabetes/metabolic syndrome, cardiovascular disorders and nonalcoholic fatty liver disease (NAFLD). As diet and exercise consistently fail to achieve any lasting results in combating obesity, pharmacological targeting of fatty acid uptake is being investigated as a potential treatment strategy. BioVision's Fatty Acid Uptake Assay Kit enables rapid determination of real-time fatty acid uptake kinetics in live cells expressing fatty acid transporter proteins. The assay uses a highly fluorescent long-chain fatty acid analogue (Ex/Em = 488/523 nm) that acts as a substrate for transmembrane fatty acid transporters (FATPs) and accumulates in intracellular lipid droplets. A proprietary non-toxic membrane-impermeable quenching agent eliminates any fluorescence arising from the extracellular space, ensuring specific measurement of intracellular fatty acid accumulation without requiring any tedious wash steps. The assay is highly sensitive, has a simple "mix-and-read" protocol and is high-throughput adaptable for screening of fatty acid uptake-modulating test compounds or treatment conditions. The kit contains a complete set of reagents sufficient for performing 100 reactions in a 96-well plate format.



II. Applications:

- Characterization of fatty acid uptake kinetics in live cells.
- Screening and characterization of drugs and test compounds for modulation of fatty acid uptake.

III. Sample Type:

- Adherent cells expressing transmembrane fatty acid transporters (e.g. differentiated 3T3-L1 or 3T3-F442A adipocytes)

IV. Kit Contents:

Components	K408-100	Cap Code	Part Number
Uptake Assay Buffer	25 ml	WM	K408-100-1
Extracellular Quenching Solution (100X)	200 µl	Amber	K408-100-2
Fluorescent Fatty Acid Probe (200X)	100 µl	Red	K408-100-3

V. User Supplied Reagents and Equipment:

- Cell line for testing: cells that express membrane fatty acid transporter proteins (differentiated 3T3-L1 adipocytes or heterologous cells stably transfected with desired FATP(s))
- Serum-free/phenol red-free cell culture medium and 5% CO₂ cell culture incubator
- Multiwell fluorescence microplate reader (capable of bottom read)
- Black-walled 96-well tissue-culture plates with clear flat bottom wells (e.g. Corning #3720 or equivalent)

VI. Storage Conditions and Reagent Preparation:

Store kit at -20°C and protect from light. Briefly centrifuge all small vials prior to opening. Open all of the reagents under sterile conditions (e.g. a cell culture hood) only. Read entire protocol before performing the assay procedure.

- **Uptake Assay Buffer:** Allow to thaw to room temperature and open under sterile conditions. Store at 4°C.
- **Extracellular Quenching Solution (100X):** Aliquot the stock solution as desired and store aliquots at -20°C, **protected from light.** Stable for at least 4 freeze/thaw cycles.
- **Fluorescent Fatty Acid Probe (200X):** Aliquot the stock solution as desired and store aliquots at -20°C, **protected from light.** Avoid repeated freeze/thaw cycles.

VII. Fatty Acid Uptake Assay Protocol:

The procedure described below employs differentiated 3T3-L1 adipocytes as a model cell line for measuring fatty acid uptake. Undifferentiated 3T3-L1 cells (ATCC CL-173) display a fibroblast-like morphology, but become rounded and accumulate lipid droplets several days after the initiation of differentiation. Primary adipocytes derived from tissues or other cell lines known to import long-chain fatty acids (such as Caco-2 or HepG2 cells) may also be used if desired.

1. 3T3-L1 Adipocyte Differentiation and Seeding:

- a. Grow 3T3-L1 fibroblasts in preadipocyte growth medium (DMEM with 10% BCS) until cells are maximally confluent (>95%). At two days post-confluence, replace preadipocyte medium with differentiation medium (DMEM/F12 with 10% FBS) containing a cocktail of induction agents (BioVision Cat. # K579). Incubate cells for 3 days at 37°C in a humidified cell culture incubator with a 5% CO₂ atmosphere, then replace differentiation medium with maintenance medium (DMEM/F12 with 10% FBS, supplemented with 1.5 µg/ml insulin). Change maintenance medium every 2-3 days, until at least 90% of the cells exhibit an adipocyte-like morphology, with intracellular lipid droplets visible by brightfield microscopy (typically 7-10 days after initiation of differentiation). Once 3T3-L1 adipocytes are fully differentiated, the insulin-supplemented medium should be replaced with standard DMEM/F12 with 10% FBS.

