



Lipid Peroxidation Colorimetric Assay Kit

03/19

(Catalog # K454-100; 100 assays; Store at 4°C)

I. Introduction:

Lipid peroxidation is a type of oxidative damage that affects cellular membranes, lipoproteins, and other lipid-containing molecules. This series of chain reactions occur between highly reactive oxygen species and unsaturated fatty acids in cell membranes, which leads to cell damage. Lipid peroxides, derived from polyunsaturated fatty acids are unstable and decompose forming malondialdehyde (MDA) and 4-hydroxyalkenals. The measurement of MDA is widely used and is considered as an indicator of lipid peroxidation. Traditionally, 2-Thiobarbituric acid (TBA) method has been widely used for the colorimetric detection and measurement of MDA in samples. However, this reaction is relatively nonspecific since both free and protein-bound MDA react with TBA. BioVision's Lipid Peroxidation Assay Kit uses a proprietary set of reagents that minimizes the interference from other lipid peroxidation products including 4-hydroxyalkenals. The principle is based on the reaction between a specific chromogenic reagent and MDA at 45 °C. The assay is simple to perform, and yields a strong, yet stable colorimetric signal (OD 586 nm) that is proportional to the amount of MDA in samples. Our assay can detect as low as 3 nmol MDA in samples.



II. Applications:

- Measurement of lipid peroxidation in tissues and serum samples

III. Sample Type:

- Serum
- Plasma
- Tissues (e.g. liver)

IV. Kit Contents:

Components	K454-100	Cap Code	Part Number
Chromogenic Reagent	6 ml	NM	K454-100-1
Diluent A	5 ml	Red NM	K454-100-2
Diluent B	5 ml	Blue NM	K454-100-3
BHT	1 ml	Purple	K454-100-4
MDA Standard	200 μ l	Yellow	K454-100-5

V. User Supplied Reagents and Equipment:

- 96-well clear plate with flat bottom
- Multi-well spectrophotometer
- Concentrated HCl (12 N)
- Acetonitrile
- Isopropyl alcohol

VI. Storage Conditions and Reagent Preparation:

Upon arrival, store kit at 4°C, protected from light. Briefly centrifuge small vials prior to opening. Read the entire protocol before performing the assay.

- **Chromogenic Reagent:** Prepare working solution (500 μ l/assay) immediately before use. Dilute 10-fold with acetonitrile, (i.e. 1 ml Chromogenic Reagent + 9 ml Acetonitrile). Dilute enough reagents for the number of assays to be performed. Discard unused diluted reagents
- **Diluent A:** Prepare working solution (~200 μ l/assay) immediately before use. Dilute Diluent A 62.5-fold with **isopropanol** (i.e. 64 μ l Diluent A + 3936 μ l Isopropyl alcohol). Dilute enough reagents for the number of assays to be performed. Discard unused diluted reagents.
- **Diluent B:** Prepare working solution (500 μ l/100 mg of tissue). Dilute Diluent B 10-fold with ddH₂O (i.e. 500 μ l Diluent B + 4.5 ml of ddH₂O). Diluted solution can be stored at 4 °C and should be used within 2 months.
- **Developer Mix:** Prepare Developer Mix (650 μ l/assay) as follows: Mix diluted Chromogenic Reagent and diluted Diluent A at a ratio of 3:1 (i.e. 495 μ l diluted Chromogenic Reagent + 165 μ l diluted Diluent A). Discard unused diluted reagents.
- **MDA Standard:** Ready-to-use as supplied.

VII. Lipid Peroxidation Assay Protocol:

1. Sample Preparation:

For tissues: Prepare homogenizing buffer (1 ml of diluted Diluent B + 3 μ l BHT). Rapidly homogenize Samples (~200 mg) in 1 ml ice-cold homogenizing buffer. Keep on ice for 10 min. Centrifuge at 3000 x g and 4°C for 10 min. Transfer supernatant to a fresh tube. Bring pH of the tissue supernatant to 1.5 using concentrated HCl.

For serum: Add 3 μ l BHT to 1 ml of Serum Sample and adjust pH to 1.5.

Sample Hydrolysis: For all Samples: Aliquot 200 μ l of acidified supernatants (pH 1.5), and hydrolyze at 60 °C for 80 minutes.

Notes: a) Approximately 10 μ l of concentrated HCl is required to bring 1 ml of Serum Sample to pH 1.5. Users need to verify the amount of HCl needed for their respective Samples.

b) Make sure to bring the heat block to 60°C before starting the hydrolysis step.

