




# Actin Polymerization/Depolymerization Assay Kit (Fluorometric)

rev 11/20

(Catalog # K457-100; 100 assays; Store kit at -20 °C)

## I. Introduction:

Actin, a highly conserved and abundant protein in eukaryotic cells, is one of the major components of cytoskeleton. It can be found as monomeric globular protein, called G actin or it can polymerize into filamentous actin, named F actin. Actin plays major roles in cell division, cell motility, cell signaling, organelle movement, etc. Mammals have 6 isoforms of Actin, which can be divided into 3 classes,  $\alpha$ ,  $\beta$  and  $\gamma$ . Muscle Actin is  $\alpha$  class and all other non-muscle actins belong to  $\beta$  and  $\gamma$ -classes. Understanding the effect of different drugs, proteins, etc. on Actin Polymerization and Depolymerization is very important for understanding cellular machinery. More importantly because Cytoskeleton is a very important target for cancer therapy. **BioVision's Actin Polymerization/Depolymerization Assay Kit** can be used to study the effect of different compounds, proteins and tissue extracts on Actin polymerization and depolymerization. The kit utilizes a proprietary Pyrene-labeled Actin molecule that develops a higher fluorescent signal if it undergoes polymerization. The signal can be easily detected using a fluorescence microplate reader. The assay is simple, high-throughput compatible and can be completed in less than three hours.



## II. Applications:

- Study and quantitate the effect of different compounds, proteins and tissue extracts on Actin polymerization and/or depolymerization.
- Evaluation of critical concentrations of actin polymerization in different conditions.

#### III. Sample Type

• Protein, Tissue Extracts, Compounds/Chemotherapeutic Agents.

## IV. Kit Contents:

Components	K457-100	Cap Code	Part Number
Buffer G	20 ml	WM	K457-100-1
Buffer P (10X)	1.5 ml	Clear	K457-100-2
Labeled Rabbit Muscle Actin	4 vials	Blue	K457-100-3
ATP (100 mM)	2 X 100 µl	Yellow	K457-100-4

# V. User Supplied Reagents and Equipment:

- 96-well plate with flat bottom. Black plates are required for this assay.
- Multi-well fluorescence microplate reader.
- DT1

# VI. Storage Conditions and Reagent Preparation:

Store kit at -20 °C, protected from light. Briefly centrifuge small vials at low speed prior to opening. Read the entire protocol before performing the assay.

- Buffer G: Store the buffer at -20 °C. Supplement Buffer G with 0.2 mM ATP and 0.5 mM DTT (For example: add 2 μl of 100 mM ATP, and 5 μl of 100 mM DTT to 993 μl of Buffer G). Prepare as needed. (DTT is not provided).
- Buffer P (10X): Store the buffer at -20 °C. Supplement Buffer P with 10 mM ATP (for example: add 10 μl of 100 mM ATP per 90 μl of 10X Buffer P). Avoid multiple freeze thaw cycles. Prepare as needed.
- Labeled Rabbit Muscle Actin: Store at -20 °C. Keep the tubes in dark to avoid photobleaching. Reconstitute vial as needed. Before experiment, reconstitute the contents of one vial with 500 µl of supplemented Buffer G. After reconstitution, *keep the tube on ice for 1 hour.* Once re-constituted, Actin can be flash frozen and saved at -80C up to 1 week. Stored actin may lose activity by 30%. Use lyophilized Actin within three months. Avoid multiple freeze thaws.
- ATP (100 mM): Ready to use. Store at -20 °C. Thaw and aliquot before use. Avoid multiple freeze thaw cycles.

## VII. Actin Polymerization/Depolymerization Assay Protocol:

Note: Avoid exposing Labeled Actin to light for extended periods of time. Protect labeled actin from light.

Actin Polymerization/Depolymerization experiments use Buffer G supplemented with ATP and Buffer P supplemented with ATP (See Section VI; Reagent Preparation). For brevity, these buffers will be referred as Supp. Buffer G and Supp. Buffer P respectively.

1. Actin Polymerization Assay: Prepare sample background, positive and sample on a black 96 microplate following the table below:

	Sample Background Control	Positive Control	Sample
Supp. Buffer G	70 µl	70 µl	60 µl
Actin	20 μΙ	20 μΙ	20 µl
Test Sample	-	-	10 µl



#### **Gentaur Europe BVBA** Voortstraat 49, 1910 Kampenhout BELGIUM Tel 0032 16 58 90 45 <u>info@gentaur.com</u>



Mix well. Incubate microplate for 15 min, or preferred incubation time based on your protocols at room temperature (RT). After incubation, **For Background Control:** add 10 µl of Supp. Buffer G; **For Positive Control and Sample Test:** Add 10 µl of Supp. Buffer P (10X) to each well containing samples and positive control. Mix and then start data acquisition (see step 3).

Note: If the initial signal is too high, incubate on ice in dark for 1 hr.

2. Actin Depolymerization Assay: First, to make polymerized Actin (F Actin), incubate Actin, Supp. Buffer P, Supp. Buffer G, test sample(s) based on the following table:

	Negative Control	Sample
Supp. 10X Buffer G	60 µl	60 µl
Supp. Buffer P (10X)	10 µl	10 µl
Actin	20 µl	20 µl


Incubate the plate at RT for one hr to polymerize Actin protected from light. To make sure that the polymerization is complete, you can take a measurement after 1 hour. **For Negative Control**: add 10 µl of the solvent of test sample/Supp. Buffer G. **For Sample**: add 10 µl of test sample, start data acquisition (see step 3).


- 3. Measurement: Measure Fluorescence Ex/Em = 365/410 nm in kinetic mode for 1 hr. at RT. Choose two time points (t<sub>FINAL</sub> & t<sub>INITIAL</sub>) in the linear range of the plot and obtain the corresponding values for the fluorescence (RFU<sub>FINAL</sub> and RFU<sub>INITIAL</sub>). For Actin Polymerization assay calculate (RFU<sub>FINAL</sub> RFU<sub>INITIAL</sub>)/Δt and for Actin Depolymerization assay calculate (RFU<sub>INITIAL</sub> RFU<sub>FINAL</sub>)/Δt.
- **4. Calculation:** To calculate the effect of test sample on Actin polymerization and/or Actin depolymerization, calculate ΔRFU<sub>P</sub>, ΔRFU<sub>G</sub> and ΔRFU<sub>S</sub> as indicated in the following equations:


ΔRFU<sub>G</sub> = Generated fluorescence of Actin in presence of Buffer G (monomeric actin)

 $\Delta RFU_p = Generated fluorescence of Actin in presence of Buffer P (polymeric actin)$ 

ΔRFU<sub>s</sub> = Generated fluorescence of Actin with test sample







**Figures:** (a) **Actin Polymerization:** Actin Polymerization is induced by Buffer P. The process is inhibited by Latrunculin A (23 μΜ) (b) **Actin Depolymerization:** Polymerized Actin is depolymerized by Latrunculin A (23 μΜ). Assays were performed following the kit protocol. **Note:** Latrunculin A is Actin polymerization inhibitor *in vitro* and *in vivo* by the formation of a 1:1 complex with monomeric G-actin. Latrunculin A acts a depolymerization agent acting on Actin filaments (Factin).

## VIII. Related Products:

Latrunculin A (2181)

Latrunculin B (2182)

Actin, Rabbit Muscle

FOR RESEARCH USE ONLY! Not to be used on humans.