

- d. Optional:** For samples having background, prepare parallel sample well(s) as sample background control. Use same amount of tissue/cell homogenate or purified enzyme as in the sample well. Adjust the final volume to 100 μ l with HDAC6 Assay Buffer.
- 2. Standard Curve Preparation:** Dilute AFC Standard to 20 μ M by adding 10 μ l of 1 mM AFC Standard to 0.5 ml of HDAC6 Assay Buffer. Add 0, 5, 10, 15, 20, and 25 μ l of diluted 20 μ M AFC Standard into a series of wells in a 96-well white plate and adjust the final volume to 100 μ l/well with HDAC6 Assay Buffer to generate 0, 100, 200, 300, 400, and 500 pmol/well of AFC Standard respectively. Mix well.

- 3. HDAC6 Substrate Mix Preparation:** Prepare 50 μ l of HDAC6 Substrate Mix per well as given below:

48 μ l HDAC6 Assay Buffer
2 μ l HDAC6 Substrate

Dissolve the Substrate Mix by vigorous vortexing. Add 50 μ l of Substrate Mix solution into each Sample, and Positive Control well.

Note: Do not add Substrate Mix to the sample Background Control and Standard wells.

- 4. Incubation and Measurement:** Mix well, cover the plate and incubate at 37°C for 30 min. To stop the reaction, add 10 μ l of Developer to each well. Mix well, cover the plate and incubate at 37°C for 10 min to generate fluorescence. The fluorescence signal is stable for at least 30 min after the addition of Developer. Measure the fluorescence at Ex/Em 380/490 nm in an end point mode at 37°C.

5. Calculations:

a. AFC Standard Curve: Obtain change in the RFU (Δ RFU) by subtracting fluorescence of the 0 Standard Controls from those containing all standards. Plot the Δ RFU against pmol of AFC. The plot should be linear; determine the slope **A** (Δ RFU/pmol) of the curve.

b. Samples: Calculate Δ RFU (**B**) for each Sample by subtracting fluorescence of corresponding Inhibitor Control. Using this value, calculate Sample HDAC6 activity using following equation.

$$\text{Sample HDAC6 Activity} \left(X, \frac{\text{U}}{\text{ml}} \right) = \frac{B \times 1000}{A \times T \times V} \times \text{Dilution Factor}$$

$$\text{Sample HDAC6 Activity} \left(\frac{\text{U}}{\text{mg}} \right) = \frac{X}{P}$$

Where, **B** = Sample HDAC6 Activity as calculated (Δ RFU)


A = Slope of the AFC standard curve (Δ RFU/pmol)

T = 30 = Assay Time (min)

V = μ l of Sample used in the assay

P = Protein concentration in the lysate (mg/ml)

Unit Definition: 1 U is the amount of HDAC6 required to deacetylate 1 pmol of HDAC6 Substrate per min under the assay conditions.



Figures: AFC Standard Curve (A), HDAC6 activity with different amounts of Positive Control (B), and in rat liver, kidney lysates along with HeLa Nuclear Extract (Cat. No. 1641), uninduced (Cat. No. 1106) and Camptothecin-induced Jurkat Cell Extracts (Cat. No. 1107) (C) are shown in the figure (n = 3). The assays were performed according to the kit protocol.

VIII. RELATED PRODUCTS:

- | | |
|--|--|
| HDAC Activity Colorimetric Assay Kit (K331) | HDAC Activity Fluorometric Assay Kit (K330) |
| HDAC Inhibitor Drug Screening Kit (Fluorometric) (K340) | HDAC1 Immunoprecipitation (IP) & Activity Assay Kit (K342) |
| HDAC2 Immunoprecipitation (IP) & Activity Assay Kit (K341) | HDAC3 Activity Fluorometric Assay Kit (K343) |
| HDAC3 Immunoprecipitation (IP) & Activity Assay Kit (K344) | HDAC3 Inhibitor Screening Kit (Fluorometric) (K363) |
| HDAC8 Activity Fluorometric Assay Kit (K348) | HDAC8 Inhibitor Screening Kit (Fluorometric) (K368) |
| HeLa Nuclear Extract (1641) | Jurkat Cell Extract (Uninduced) (1106) |
| Jurkat Cell Extract (Induced) (1107) | |

FOR RESEARCH USE ONLY! Not to be used on humans.