




Phosphatidylethanolamine Assay Kit (Fluorometric)

rev 06/20

(Catalog # K499-100; 100 assays; Store at -20°C)

I. Introduction:

Phosphatidylethanolamine (PE), also known as *cephalin* is the second most abundant phospholipid in animal and plant tissues, and is present on the cytoplasmic side of the plasma membrane. Phosphatidylethanolamine is a neutral phospholipid consisting of a phosphatidyl group ester linked to an ethanolamine molecule. Its functions include membrane fission/fusion, maintenance of membrane curvature and stabilization of membrane proteins, since it can form hydrogen bonds with proteins through an ionizable amine group. It acts as a chaperone during assembly of membrane proteins and aids in their translocation from the cytoplasm to the membrane. It is also involved in secretion of very low density lipoproteins in the liver. **BioVision's Phosphatidylethanolamine Assay Kit** is a microplate based enzymatic assay for the quantitation of PE in cells and tissues. PE Converter hydrolyses PE to an intermediate, which converts a colorless probe to a fluorescent product via enzymatic reaction (Ex/Em = 535/587 nm). The intermediate formed through PE converter hydrolysis is specific to phosphatidylethanolamine. Thus no other phospholipids (*i.e.* phosphatidylcholine, phosphatidylinositol or phosphatidic acid) will be detected, making the kit highly specific. This assay kit can detect as low as 0.2 nmol per well.



II. Applications:

- Measurement of Phosphatidylethanolamine content in lipid extracts from cells and tissues

III. Sample Type:

- Cell lipid extract
- Tissue lipid extract

IV. Kit Contents:

Components	K499-100	Cap Code	Part Number
PE Assay Buffer	25 ml	WM	K499-100-1
PE Converter	1 vial	White	K499-100-2
PE Developer	600 µl	Orange	K499-100-3
PE Enzyme Mix	1 vial	Green	K499-100-4
PE Probe	200 µl	Red	K499-100-5
PE Standard (1 mM)	100 µl	Yellow	K499-100-6

V. User Supplied Reagents and Equipment:

- 96-well flat bottom clear plate
- Multi-well spectrophotometer
- Deionized water
- Triton X-100 (peroxide free)
- Water bath/heating plate

VI. Storage Conditions and Reagent Preparation:

Upon arrival, store the kit at -20°C, protected from light. Centrifuge vials prior to opening. Read the protocol before performing the assay.

- **PE Assay Buffer:** Warm to room temperature (RT) before use.
- **PE Converter and PE Enzyme Mix:** Store at -20°C. Lyophilized vials are stable for at least 6 months. Reconstitute each vial in 220 µl assay buffer before use. Aliquot remaining components. Store at -20°C. Reconstituted vials are stable for at least two months.
- **PE Developer:** Store -20°C. Thaw on ice before use. Divide into aliquots and store the remaining at -20°C. Product is stable for at least six months.
- **PE Probe:** Thaw at RT before use.
- **PE Standard:** Thaw in a water bath at 45°C for 15-20 min. The solution should look clear. Divide into aliquots and store the remaining at -20°C. Thaw in a water bath at 45°C before next use.

VII. Phosphatidylethanolamine Assay Protocol:

1. Sample Preparation: Prepare a 5% (v/v) solution of peroxide free triton X-100 in water. Homogenize tissue (~100 mg; non-perfused) or cells (~1 million) in 1 ml solution containing 5% Triton X-100 in water. Protein content in the sample may be determined at this stage if desired. We recommend *BCA protein assay kit (BV# K813-2500)*. Heat the samples to 80°C in a water bath for 5 - 10 minutes or until the solution becomes cloudy, then cool down to RT. Repeat the heating step once more to solubilize all lipids and allow the solution to cool to RT again. Centrifuge (10000 X g, 10 min, 4°C) and collect the supernatant, which contains solubilized lipids. If not used immediately, store the supernatant in -80°C. Add 2 to 10 µl of samples into wells of a 96-well clear plate. For each sample prepare two wells; "Sample Background Control" and "Sample". Bring the volume in "Sample" wells to 50 µl and in "Sample Background Control" to 70 µl using PE Assay Buffer respectively.

Note:

- a. Different dilutions of sample should be tested to make sure that Phosphatidylethanolamine concentration falls in the linear range of the assay.



- b. Samples should be diluted using PE Assay Buffer.
- 2. Phosphatidylethanolamine Standard Curve:** Dilute the 1 mM Phosphatidylethanolamine 1:10 in PE Assay Buffer to obtain 100 μ M PE solution. For example, to prepare 200 μ l of 100 μ M PE solution mix 20 μ l of 1 mM PE solution in 180 μ l of PE Assay Buffer and incubate at 45°C for 30 minutes. Add 0, 5, 10, 20, 30, and 40 μ l of the 100 μ M standard to wells of the 96 well plate to obtain 0, 0.5, 1, 2, 3 and 4 nmol, of Phosphatidylethanolamine per well. Bring up the total volume in these wells to 50 μ l with PE Assay buffer.
- 3. Converter Mix Preparation:** Mix enough reagents for the number of assays to be performed. For each Sample and Standard well, prepare 20 μ l Converter Mix as mentioned below:

<u>Converter Mix</u>	
PE Assay Buffer	18 μ l
PE Converter	2 μ l

Add the Converter Mix to wells containing Samples and Standards. Mix well. *Do not add the convertor mix to "Sample Background Control" wells.* Incubate at 45°C for 1 hr.

- 4. Reaction Mix:** Mix enough reagents for the number of assays to be performed. For each well, prepare 30 μ l:

<u>Reaction Mix</u>	
PE Assay Buffer	20 μ l
PE Developer	6 μ l
PE Enzyme Mix	2 μ l
PE Probe	2 μ l

Add the Reaction Mix to all wells. Mix well. Incubate at 40°C for 3 hr.

- 5. Measurement:** Record fluorescence in end point mode at Ex/Em = 535/587 nm.
- 6. Calculations:** Subtract 0 Standard PE reading from all PE Standard readings. Plot the Phosphatidylethanolamine Standard Curve. Subtract the Sample background control readings from Sample readings. If 0 PE Standard readings are higher than the Sample Background Control readings, subtract those from Sample readings instead. Apply corrected RFU to Standard Curve to get B nmol PE in the Sample well.

$$\text{PE concentration in sample: } C = (B / V) \times D \text{ (nmol/ml or } \mu\text{M)}$$


Where **B** = Amount of Phosphatidylethanolamine in the Sample well from the Standard Curve (nmol)

V = Volume of sample added into the well (ml)

D = Dilution factor

PE molecular weight: 726 g/mol

PE concentrations can also be expressed as nmol PE per mg protein or nmol PE per mg tissue weight.



Figures: (a) Phosphatidylethanolamine Standard Curve. (b) Phosphatidylethanolamine content in rat liver (100 μ g wet tissue) and HEK 293 cells (25 μ g protein). Sample preparation and assay was carried out according to kit protocol.

VIII. RELATED PRODUCTS

Phospholipid Assay Kit (K351)

Cardiolipin Assay Kit (K944)

Phosphatidylcholine Colorimetric/Fluorometric Assay Kit (K576)

FOR RESEARCH USE ONLY! Not to be used on humans.