

- c. Add 2-20 μl of neutralized sample extract to desired well(s) in a white 96-well plate. For each test sample, prepare *two parallel sample wells*, with one well serving as a sample background control. Adjust the volume of all wells to 80 μl /well with Collagen Assay Buffer.

Notes:

- It is important to keep samples chilled during the homogenization procedure. Heat generated by homogenization can cause denaturation and crosslinking of soluble tropocollagen fibrils, rendering them insoluble in acid solutions.
- Soluble collagen levels and the effectiveness of acid solubilization can vary tremendously between tissues. Collagen present in "tough" tissues (such as cartilage or connective tissue) is highly crosslinked and tends to be resistant to acid solubilization.

Cell Culture Medium (Secreted Soluble Collagen):

- d. Collagen secreted into cell culture medium may be assayed directly, without the need for acid solubilization. Remove a sample of culture medium without detaching cells, centrifuge at 10000 x g for 15 min at 4°C to pellet any debris and transfer the clarified supernatant to a new microfuge tube. Add 10-20 μl of clarified medium to two parallel wells (one well will serve as a sample background control) and adjust the volume to 80 μl /well with Collagen Assay Buffer.
2. **Standard Curve Preparation:** Prepare a 0.2 mg/ml collagen solution by adding 20 μl of the 2 mg/ml Collagen I Standard to 180 μl of ddH₂O. Add 0, 2, 4, 6, 8, and 10 μl of the 0.2 mg/ml working solution into a series of wells, generating 0, 0.4, 0.8, 1.2, 1.6 and 2 μg of collagen/well. Adjust the volume of all standard wells (including the 0 μg /well reagent blank) to 80 μl /well with Collagen Assay Buffer.

3. Reaction Preparation:

- a. Prepare a working solution of collagenase by diluting the reconstituted Collagenase Enzyme Mix stock solution at 1:10 ratio with Collagen Assay Buffer. Prepare enough of the working solution to add 20 μl to each reaction well (for 10 reactions, mix 20 μl of Collagenase Enzyme Mix stock and 180 μl of Collagen Assay Buffer). Add 20 μl of collagenase working solution to test sample and standard wells. For sample background control wells, add 20 μl of Collagen Assay Buffer only. Incubate plate at 37°C for 60 min.
- b. Following incubation, prepare detection reaction solution by diluting the reconstituted Peptide Labelling Reagent stock in Detection Reaction Buffer at a 1:30 ratio. Prepare enough of the working solution to add 75 μl to each reaction well (for 10 reactions, mix 25 μl of Peptide Labelling Reagent stock and 725 μl of Detection Reaction Buffer). Add 75 μl of detection reaction solution to all test sample and standard wells (including sample background control wells) and incubate plate (protected from light) at 37°C for 5 min.
- c. Prepare developer working solution by diluting the Developer Solution stock in ddH₂O at a 1:10 ratio. Prepare enough diluted developer solution to add 25 μl to each reaction well (for 10 reactions, mix 25 μl of Developer Solution stock and 225 μl of ddH₂O). Add 25 μl of detection reaction solution to all test sample and standard wells (including sample background control wells). Incubate the plate (protected from light) at 37°C for 15 min with gentle orbital shaking to ensure well contents are effectively mixed.

4. **Measurement:** Measure the fluorescence (Ex/Em = 376/468 nm) of all test sample and standard curve wells in endpoint mode.

5. **Calculations:** For the collagen standard curve, subtract the reagent blank (0 μg /well collagen standard) fluorescence (RFU) value from all standard readings, plot the background-subtracted values and calculate the slope of the standard curve. For test samples, subtract the corresponding sample background control well RFU value from the sample reading ($F = RFU_{\text{Sample}} - RFU_{\text{BC}}$) and apply the background-subtracted fluorescence (F) to the standard curve to get B μg of soluble collagen in the well.


$$\text{Sample Soluble Collagen Concentration} = \frac{B}{V} \times D = \mu\text{g}/\mu\text{l}$$

Where: B is the amount of collagen, calculated from the standard curve (in μg)

V is the volume of sample added to the well (in μl)

D is the sample dilution factor (if applicable, $D=1$ for undiluted samples)

Note: The calculation above gives the amount of collagen in the sample added to the well. The dilution factor D is only needed if the sample is diluted *after* the neutralization step. When calculating the amount of collagen in the original sample homogenate, remember to account for the 2-fold dilution that occurs during neutralization of the acidic homogenate.



Figures: (a) Collagen I Standard curve. (b) Estimation of acid-soluble collagen content in rat tissues. Rat leg muscle, lung and heart samples were homogenized in 0.5 M acetic acid, incubated overnight at 4°C to solubilize collagen and neutralized with 0.5 M NaOH. Collagen levels (calculated as $\mu\text{g collagen}/\text{mg wet tissue}$) for the samples were: $3.83 \pm 0.47 \mu\text{g}/\text{mg}$ for muscle, $1.47 \pm 0.22 \mu\text{g}/\text{mg}$ for lung and $2.85 \pm 0.08 \mu\text{g}/\text{mg}$ for heart. (c) Estimation of secreted soluble collagen in cultured cell growth medium (DMEM/F12 medium). 3T3-L1 fibroblasts were grown to ~80% confluence in T-75 flasks and then treated with either vehicle or TGF- β (Cat. # 6479; 100 ng/ml), a cytokine known to stimulate synthesis and secretion of tropocollagen fibrils. After 48 hours of treatment, the culture medium was removed, centrifuged to remove debris and was assayed directly (each 10 μl clarified medium, undiluted). TGF- β treatment resulted in a roughly 1.5-fold increase in collagen secretion. Data are mean \pm SEM of 3 replicates, assayed according to the kit protocol.

VIII. RELATED PRODUCTS:

- | | | |
|--|--|--|
| Total Collagen Assay Kit (K218 and K406) | Hydroxyproline Assay Kit (K555 and K226) | Glycine Assay Kit (K589) |
| Collagenase Activity Assay Kit (K792) | Collagenase Inhibitor Screening Kit (K833) | Collagen-Fluorescein Conjugate (M1304) |

FOR RESEARCH USE ONLY! Not to be used on humans.