



# EZClick™ Total Phospholipid Assay Kit (Cell-Based)

rev 0720

(Catalog # K537-100; 100 assays; Store at -20°C)

## I. Introduction:

Phospholipids are major component of the bilayers of all plasma membranes. A single phospholipid molecule consists of a phosphate group on one ("head"), and two side-by-side chains of fatty acids that make up the "tails". The phosphate head groups can be modified with organic molecules such as Choline (Cho). Cho-containing phospholipids (Phosphatidylcholines; PC) are critical for structural membrane integrity, cellular metabolism and signaling either as individual molecules or precursors of secondary messengers. Changes in global synthesis of Cho-containing phospholipids are an essential parameter in analysis of cellular response to both, physiological and pathological conditions, environmental stress, or drug treatment. To date, phospholipids biochemistry, cell biology and metabolism remain obscure, due to limited methods for their direct cellular visualization. **BioVision's EZClick™ Global Phospholipid Synthesis Assay Kit** offers a simple and robust method to label and visualize newly synthesized phospholipids *in vivo*. Based on the metabolic incorporation of the choline analogs directly into their structure, modified phospholipid molecules can be detected with high sensitivity and spatial resolution by click chemistry with azide-containing dyes (Ex/Em= 494/521 nm). This kit enables quantitative analyses of global biosynthesis/turnover of Cho-containing phospholipids in cells. Cells show strong incorporation of Cho analogs into all classes of phospholipids that can be assayed by microplate reader and fluorescence microscope. The kit provides sufficient materials for 100 assays.

## II. Applications:

- Detection and quantification of biosynthesis, subcellular localization and turnover of phospholipids
- Evaluating effects of anti-cancer drugs and genotoxic agents on phospholipids
- Screening for genotoxic compounds and effectors of phospholipid biosynthesis in proliferating cells

## III. Sample Type:

- Suspension or adherent cell cultures

## IV. Kit Contents:

Components	K537-100	Cap Code	Part Number
EZClick™ Wash Buffer (10X)	25 ml	NM	K537-100-1
Fixative Solution	10 ml	WM	K537-100-2
Permeabilization Buffer (10X)	25 ml	NM/Blue	K537-100-3
EZClick™ Phospholipid Label (1000X)	10 µl	Clear	K537-100-4
Copper Reagent (100X)	100 µl	Blue	K537-100-5
EZClick™ Fluorescent Azide (100X)	100 µl	Green/Amber Vial	K537-100-6
Reducing Agent (20X)	500 µl	Yellow	K537-100-7
EZClick™ Total DNA Stain (1000X)	20 µl	Blue/Amber Vial	K537-100-8

## V. User Supplied Reagents and Equipment:

- Tissue culture vessels and appropriate culturing media
- A 6-, 12-, 24-, or 96-well clear plates should be used only for cell culturing. The measurement of fluorescence should be performed in black opaque plates. Alternatively, sterile opaque plates can be used for both, culturing and measurements
- Phosphate Buffered Saline (PBS, pH 7.4)
- Multi-well spectrophotometer and Fluorescence microscope (optionally) capable of measuring Ex/Em= 494/521 nm spectra

## VI. Storage Conditions and Reagent Preparation:

Upon arrival, store the entire kit at -20°C protected from light. Briefly centrifuge small vials prior to opening. Read the entire protocol before performing the assay.

- **10X Wash Buffer and 10X Permeabilization Buffers: Thaw at 37°C to dissolve completely.** Dilute the 10X stocks 1:10 in sterile water, mix well. Store at 4°C.
- **Fixative Solution: Divide into aliquots and store at -20°C, protected from light.**
- **Remaining components:** Store at -20°C protected from light. While in use, keep on ice and minimize light exposure.

## VII. Total Phospholipid Assay Protocol:

### Notes:

This assay was developed with Jurkat (suspension) and HeLa (adherent) cells and can be modified for any cell line. The following protocol has been optimized for a 96-well opaque plate at  $1 \times 10^6$  cells per well, using fewer cells per well or clear plate will result in decreased signal. We suggest testing growth conditions, cell number per well and several concentrations of the EZClick™ Phospholipid Label, to find the best experimental design for your cell type. The assay volume is 100 µl; adjust volumes accordingly for other plate formats. Avoid stressing the cells by washes or temperature changes prior to incubation with EZClick™ Phospholipid Label. All steps should be carried out at room temperature (RT) unless otherwise specified; equilibrate all buffers to RT prior to the experiment.

### 1. Labeling of control and experimental cells:

- Plate suspension or adherent cells at a desired density and allow for overnight recovery before treatment. Ensure that adherent cells are sub-confluent. Include appropriate controls and account for cell loss during the processing. **Negative control** (cells not exposed to the Phospholipid Label or treatment), **Background control** (cells treated with EZClick cocktail only), **Positive control** (cells incubated with 1X Phospholipid Label only).
- Next day, remove the media and treat the cells with appropriate effectors according to your protocol; do not add treatment to the positive and negative control cells. For **suspension cells**: Centrifuge the plate at 500 x g (or the lowest centrifuge setting) for 5 minutes at RT to pellet the cells. Tilt the plate and **gently** remove the media with a pipette tip. Avoid excessive centrifugation speeds, which can damage the cells. *Use these centrifugation settings throughout the entire protocol!*
- Dilute EZClick™ Phospholipid Label (1000X) to 1X final concentration with culture medium and add into the **Experimental** and **Positive control** cells respectively. Incubate the cells for additional 24 hours, or time required by your experimental protocol in a 37°C incubator. Do not add the EZClick™ Phospholipid Label into the **Negative control** cells. Do not remove the drug-containing media during incubation with 1X Phospholipid Label to avoid potential reversibility of drug action on label incorporation. Harvest the

