

3. DNA Probe: Prepare enough DNA probe by diluting the 200X DNA Probe stock (e.g. add 10 μ l of 200X DNA Probe stock into 1990 μ l of 1X DNA buffer). Add 100 μ l of the 1X DNA Probe into each well containing DNA standards and DNA samples. Mix well, cover the plate and incubate for 5 min at 25 °C and protected from light.

Note: Working solutions should be prepared prior to running your experiments, protected from light and *should be used within 4 hours*.

4. Measurement: Measure fluorescence (Ex/Em = 492/528nm) in a microplate reader in endpoint mode.

5. Calculation: Subtract 0 standard readings from all standard readings. Plot the wide range and low DNA concentration standard curves separately. Apply fluorescence from DNA samples to the DNA standard curves to get **B** ng of DNA in the sample well.

$$\text{Sample DNA Concentration} = \frac{B}{V} \times D \text{ (ng/ } \mu\text{l)}$$

Where: **V** is the volume of sample added to the well (in μ l)

D is the sample dilution factor (if applicable, D=1 for undiluted samples)




Figure: (a) λ DNA standard curve; (b) Specificity of the assay: the probe specifically recognizes dsDNA only. ssDNA and RNA are not detected. (c) 150 ng/well of λ DNA mixed with different concentrations of inorganic and organic contaminants commonly found in DNA samples ([contaminants] based on a 10 μ l of DNA samples/well). Assays were performed following the kit protocols.

VIII. RELATED PRODUCTS:

- Bacterial Genomic DNA Isolation Kit (K309)
- DNA Damage Quantification Colorimetric Kit (K253)
- DNA Cleanup Maxi Kit (K1369)
- Genomic DNA Isolation Kit (K281)
- Mitochondrial DNA Isolation Kit (K280)

FOR RESEARCH USE ONLY! Not to be used on humans.