

- a. Prepare a 1:10 dilution of viability stain with cell suspension by adding 2 µl of Viability Stain to 18 µl of washed PBMC suspension in a 1.5 ml centrifuge tube. Carefully resuspend cell pellet. Inoculate hemocytometer with 10 µl of stained PBMC suspension. Determine and record the total cell count with a Bright-field microscope.
- b. With the same Region of Interest (ROI) in view, illuminate the slide with light from the fluorescent lamp with the FITC/TRITC filter combination to count the fluorescent cells. *If a small amount of incidental white light illuminates the hemocytometer, the grid will be visible allowing the viewer to see the same ROI as visible with fluorescent light.* Live cells will fluoresce green. Dead cells will fluoresce red. Tally the number of green and red cells to complete the calculations using the equations below.

3. Measurement:

- a. Number of RBCs = Total Cell Count (hemocytometer) – Total Number of Fluorescent Cells
- b. Live Cells Fraction = Number of Green Fluorescent Cells/Total Number of Fluorescent Cells × 100
- c. Dead Cells Fraction = Number of Red Cells/Total Number of Fluorescent Cells × 100
- d. Percentage RBCs = Number of RBCs/Total Cell Count (hemocytometer) × 100

Note: Contamination of PBMCs with RBCs may affect downstream applications including, but not limited to flow cytometry or T Cell Killing Assays. For that reason, we recommend repeating the separation with remaining density gradient media if RBCs are >10% of total cells count.




Figure A. Illustration of conical tube held at appropriate angle while blood layered on top of DGM. **Figure B.1.** Layers of DGM and whole blood prior to and after centrifugation showing the separation of layers in the conical tube. **Figure B.2** illustrates separation of four layers (plasma, PBMCs, DGM and RBCs). **Figure C.** Brightfield image of hemocytometer showing total cells(left); image from Fluorescent microscope with Rhodamine/FITC filters of same ROI showing live (green) and dead (red) cells (middle); merge of two panels (right).

VIII. Related Products

- 1X Red Blood Cell Lysis Buffer (5830)
- Propidium Iodide (1056)
- 10X Red Blood Cell Lysis Buffer (5831)

FOR RESEARCH USE ONLY! Not to be used on humans.