




Arginase I (ARG1) Inhibitor Screening Kit (Colorimetric)

2/18

(Catalog #K567-100; 100 assays; Store at -20°C)

I. Introduction:

Arginase (EC 3.5.3.1) is a manganese-containing enzyme which catalyzes the conversion of arginine into urea and ornithine, which is the final reaction in the urea cycle. Arginase I (ARG1) is the liver isoform of arginase. Recent studies showed that ARG1 expression by mature myeloid cells in tumor environment as demonstrated in a 3LL murine lung carcinoma model causes L-Arginine depletion by tumor-associated myeloid cells (TAMC). L-arginine depletion suppresses immune-response against tumor cells due to inhibition to T-cell proliferation. In addition, the depletion of arginine increases the reactive nitrogen species (NOS) and reactive oxygen species (ROS), which, in consequence, induce T-cell apoptosis and supports antigenic cell proliferation. BioVision's Arginase I (ARG1) Inhibitor Screening Kit is designed for screening ARG1 inhibitors. Two substituted-2-amino-6-hexanoic acids have been studied as arginase inhibitor and in this kit, Amino-2-Borono-6-Hexanoic Acid (ABH) is provided as a positive control. The ARG1 activity is monitored by the increase in absorbance readings (OD 450 nm), while potential inhibitors will cause a decrease of absorbance. The assay kit is simple, quick and can be used to identify and characterize ARG1 inhibitors in a high-throughput format.



II. Applications:

- Screening for inhibitors of human arginase I (ARG1)

III. Kit Contents:

Components	K567-100	Cap Code	Part Number
Assay Buffer	25 ml	WM	K567-100-1
ARG1 Substrate	1 vial	White	K567-100-2
ARG1 Probe Mix A	12 ml	NM/Blue	K567-100-3
ARG1 Probe Mix B	12 ml	NM/Brown	K567-100-4
Human ARG1	1 vial	Green	K567-100-5
ABH (in DMSO)	20 μ l	Purple	K567-100-6

IV. User Supplied Reagents and Equipment:

- 96-well clear plate with flat bottom
- Multi-well spectrophotometer (ELISA reader)

V. Storage Conditions and Reagent Preparation:

Store kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay. Use within two months of opening.

- **Assay Buffer:** Warm to room temperature before use. Store at 4°C or -20°C.
- **ARG1 Substrate:** Reconstitute with 250 μ l dH₂O. Pipette up and down to dissolve. Store at -20°C.
- **ARG1 Probe Mix A:** Ready to use as supplied. Warm to room temperature before use. Store at 4°C or -20°C. Keep away from light.
- **ARG1 Probe Mix B:** Ready to use as supplied. Warm to room temperature before use. Store at 4°C or -20°C. Keep away from light.
- **Human ARG1:** Reconstitute with 220 μ l ARG1 Assay Buffer. Pipette up and down to dissolve. Aliquot and store at -20°C. Avoid repeated freeze/thaw cycles. Keep on ice while in use.
- **ABH (in DMSO):** Ready to use as supplied. Warm to room temperature before use

VI. ARG1 Inhibitor Screening Assay Protocol:

1. Test compounds, Inhibitor Control, Enzyme Control & Background Control Preparations:

Dissolve candidate inhibitors at 1000X highest final test concentration into an appropriate solvent. Dilute to 5X the desired test concentration with ARG1 Assay Buffer. Add 10 μ l diluted test inhibitor or Assay buffer into designated wells as sample screen [S]. Add 10 μ l of Assay Buffer to a well designated as Enzyme Control [EC] (no inhibitor) respectively. **For ABH control:** Dilute ARG1 inhibitor by adding 2 μ l of the stock solution into 18 μ l of ARG1 Assay Buffer. Add 10 μ l of the diluted ABH inhibitor into one wells labeled as Inhibitor Control [IC]. **For Background Control [BC]:** Add 10 μ l of the diluted ABH inhibitor and 30 μ l of Assay Buffer in a well designated as Background Control [BC]. If you are screening test compounds that have significant absorbance (OD 450 m) at the 5X final concentration prepare background controls as described above.

2. Enzyme Solution Preparation: Mix enough reagents for the number of assays to be performed. For each well, prepare 30 μ l ARG1 Enzyme Solution:

Assay Buffer	28 μ l
ARG1 Enzyme	2 μ l

Mix and add 30 μ l of the ARG1 enzyme solution into all wells except Background Control well(s). Add 30 μ l of Assay Buffer into Background Control well(s). Mix well, and incubate the plate for 5 min at 37 °C.

