




Transglutaminase Activity Assay Kit (Colorimetric)

2/15

(Catalog # K571-100; 100 assays; Store at -20°C)

I. Introduction:

Transglutaminases (EC 2.3.2.13) are calcium dependent enzymes that catalyze the post-translational modification of proteins by formation of isopeptide bonds. This occurs either through protein cross-linking via formation of γ -glutamyl- ϵ -lysine bonds or through incorporation of primary amines at selected peptide-bound glutamine residues. The transglutaminase enzyme family comprises the intracellular forms (TG1, TG3 and TG5) expressed mostly in the epithelial tissue; TG2 which is both intracellular and extracellular and expressed in various tissue types; TG4 which is expressed in the prostate gland; factor XIII which is expressed in blood; TG6 and TG7, whose tissue distribution is unknown and band 4.2 (lacking enzymatic activity) which is present on erythrocyte membranes. Transglutaminases also exhibit GTPase, phosphodiesterase and protein kinase activity. Transglutaminases are associated with certain neurological and autoimmune disorders and also cancer. BioVision's Transglutaminase Activity Assay kit utilizes the deamidation reaction of the transglutaminase enzyme with a donor and acceptor substrate resulting in the formation of a hydroxamate product. The hydroxamate product reacts with the Stop Solution forming a purple complex that can be measured colorimetrically at 525 nm. The limit of quantification of this assay is ~10 μ J or 80 ng of recombinant hTG2 enzyme.



II. Applications:

Quantification of Transglutaminase enzyme activity

III. Sample Type

- Cell lysate
- Tissue lysate
- Recombinant Transglutaminase

IV. Kit Contents:

Components	K571-100	Cap Code	Part Number
TG Assay Buffer	12 ml	WM	K571-100-1
Homogenization Buffer (10x)	10 ml	NM Clear	K571-100-2
1M DTT	125 μ l	Blue	K571-100-3
Donor Substrate	1 Vial	NM Brown	K571-100-4
Acceptor Substrate	2 Vials	Red	K571-100-5
Hydroxamate Standard	1 Vial	Yellow	K571-100-6
Stop Solution	5 ml	NM Red	K571-100-7
Positive Control	1 Vial	Brown	K571-100-8
Plate Sealer	1	-	K571-100-9

V. User Supplied Reagents and Equipment:

- 96-well clear plate with flat bottom.
- Multi-well spectrophotometer capable of absorbance detection
- Protease Inhibitor Cocktail (Cat. # K272 or equivalent)
- Homogenizer
- 20% glycerol in water.
- Phosphate Buffered Saline (PBS)
- Deionized water

VI. Reagent Preparation and Storage Conditions:

Store kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay.

- **TG Assay Buffer:** Warm to 37°C before use.
- **Homogenization Buffer:** Thaw and make 1x buffer by diluting 1 volume 10x Homogenization Buffer with 9 volumes deionized water. Keep on ice while in use.
- **1 M DTT:** Store at -20°C. Thaw and keep on ice while in use. Use within two months.
- **Donor Substrate:** Reconstitute with 1.1 ml deionized water. Aliquot and store at -80°C. Avoid repeated freeze/thaw. Use within two months.
- **Acceptor Substrate:** Reconstitute each vial with 550 μ l deionized water as needed. Aliquot and store at -80°C. Use within two months.
- **Hydroxamate Standard:** Reconstitute with 330 μ l deionized water to make 10 mM stock Standard solution.
- **Stop Solution:** Store at 4°C or -20°C. Thaw before use. Keep on ice while in use.
- **Positive Control:** Reconstitute with 30 μ l 20% glycerol in water (not provided). Mix well, aliquot and store at -80°C. Avoid repeated freeze/thaw. Thaw and mix gently before use.

VII. Transglutaminase Activity Assay Protocol:

1. **Sample Preparation:** Add DTT to 1x Homogenization Buffer at a final concentration of 0.2 mM. Make fresh as needed. Rinse tissue with PBS and transfer ~100 mg of fresh or frozen tissue (stored at -80°C) to a prechilled homogenizer. Add 500 μ l cold Homogenization Buffer (with DTT) containing protease inhibitor cocktail (not provided) and thoroughly homogenize tissue on ice using a dounce or electrical homogenizer. To prepare cell extract, add 150-300 μ l cold Homogenization Buffer (with DTT) containing protease inhibitor

