



	Reaction Mix	*Background Control Mix
Lactate Assay Buffer	46 µl	48 µl
Lactate Enzyme Mix	2 µl	--
Probe	2 µl	2 µl

Mix well. Add 50 µl of the Reaction Mix to each well containing the Lactate Standards & test samples and mix well.

Note:

* For samples having high background, add 50 µl of Background Control Mix to sample background control well(s).

a. The fluorometric assay is ~10 times more sensitive than the colorimetric assay. Use 0.4 µl of the probe per reaction to decrease the background reading.

4. Measurement: Incubate the reaction for 30 min. at room temperature, protected from light. Measure absorbance (OD 570 nm) or fluorescence (Ex/Em = 535/590 nm) in a microplate reader.

5. Calculation: Subtract 0 Standard reading from all readings. If sample background control reading is significant then subtract the sample background control reading from sample reading. Plot the Lactate Standard Curve. For unspiked samples, apply the corrected OD to the Lactate Standard Curve to get B nmol of Lactate in the sample well.

$$\text{Sample Lactate concentration (C)} = \frac{B}{V} \times D \text{ nmol/}\mu\text{l or mM}$$

Where: **B** is the amount of Lactate in the sample well (nmol)

V is the sample volume added into the reaction well (µl)

D is the sample dilution factor

Note: For spiked samples, correct for any sample interference by subtracting the sample reading from spiked sample reading.

$$\text{For spiked samples, Lactate amount in sample well (B)} = \left(\frac{(\text{OD}_{\text{sample (corrected)}})}{(\text{OD}_{\text{sample + Lactate Std (corrected)}}) - (\text{OD}_{\text{sample (corrected)}})} \right) * \text{Lactate Spike (nmol)}$$

Lactic acid molecular weight: 90.08.




Figure: Lactate Standard Curve: colorimetric (a), fluorometric (b). c.) Quantitation of lactate in human saliva & serum. Saliva was centrifuged at 10000 x g for 10 min. at 4°C. 3 µl of supernatant was spiked with a known amount of lactate (4 nmol) as internal Standard. Serum was diluted 20-fold and 3 µl was assayed. Calculated concentrations (mg/dl): Saliva: 0.42 ± 0.05; Serum: 12.9 ± 0.5. Assays were performed according to the kit protocol.

IX. RELATED PRODUCTS:

Lactate Colorimetric Assay Kit II (K627)

PicoProbe™ Lactate Fluorometric Assay Kit (K638)

Lactate Dehydrogenase (LDH) Activity Assay Kit (K726)

D-Lactate Colorimetric Assay Kit (K667)

PicoProbe™ D-Lactate Fluorometric Assay Kit (K668)

LDH Antibody (3842)

FOR RESEARCH USE ONLY! Not to be used on humans.