

## L-Amino Acid Quantitation Colorimetric/Fluorometric Kit

(Catalog #K639-100; 100 assays; Store kit at -20°C)

### I. Introduction:

L-Amino acids are the most essential elements in biology. Accurately quantitating L-amino acids in body fluids or purified samples may provide valuable information for diagnostic or basic research studies. BioVision's L-Amino Acid Assay Kit provides a convenient means for directly detecting L-amino acids in biological samples. There is no requirement for sample pretreatment or purification when using this kit. The L-amino acid levels can be quantified using fluorometric (at Ex/Em = 535/587 nm) or colorimetric (at  $\lambda = 570$  nm) methods in 96-well plates.

### II. Kit Contents:

Components	K639-100	Cap Code	Parts Number
L-Amino Acid Assay Buffer	25 ml	WM	K639-100-1
L-Amino Acid Probe	0.2 ml	Red	K639-100-2A
L-Amino Acid Enzyme Mix	1 Vial	Green	K639-100-4
L-Amino Acid Standard (4nmol/ $\mu$ l)*	300 $\mu$ l	Yellow	K639-100-5

\*Mixture of all amino acids at equal molar ratio. Glycine is not detectable in this assay.

### III. Reagent Preparation and Storage Conditions:

**Probe:** Ready to use as supplied. Warm to room temperature to thaw the DMSO solution before use. Store at -20°C, protect from light. Use within two months.

**Enzyme Mix:** Dissolve in 220  $\mu$ l L-Amino Acid Assay Buffer. Pipette up and down to complete dissolve the content. Store at -20°C. Use within two months.

### IV. Assay Protocol:

- Standard Curve Preparations:** For the colorimetric assay, add 0, 2, 4, 6, 8, 10  $\mu$ l L-Amino Acid Standard into a series of wells of a 96-well plate to generate 0, 8, 16, 24, 32, 40 nmol/well of L-Amino Acid Standard. Adjust volume to 50  $\mu$ l/well with L-Amino Acid Assay Buffer.

For the fluorometric assay, dilute the L-Amino Acid to 0.4 nmol/ $\mu$ l by adding 10  $\mu$ l of the L-Amino Acid to 90  $\mu$ l of L-Amino Acid Assay Buffer, mix well. Add 0, 2, 4, 6, 8, 10  $\mu$ l into each well individually to generate 0, 0.8, 1.6, 2.4, 3.2, 4.0 nmol/well of the L-Amino Acid Standard. Adjust volume to 50  $\mu$ l/well with L-Amino Acid Assay Buffer.

- Sample Preparations:** Prepare test samples in a final volume of 50  $\mu$ l/well with L-Amino Acid Assay Buffer in the 96-well plate. We suggest using several doses of your sample to ensure the readings are within the standard curve range.
- Reaction Mix Preparation:** Mix enough reagents for the number of assays performed: For each well, prepare a total 50  $\mu$ l Reaction Mix containing:
  - 46  $\mu$ l L-Amino Acid Assay Buffer
  - 2  $\mu$ l L-Amino Acid Probe
  - 2  $\mu$ l L-Amino Acid Enzyme Mix,
- Add 50  $\mu$ l of the Reaction Mix to each well containing the L-Amino Acid standard or test samples.
- Incubate the reaction for 30 min at 37°C, protect from light.


at 570nm for the colorimetric assay or fluorescence at Ex/Em = 535/590 nm in a micro-plate reader.

- Correct background by subtracting the value of the 0 L-amino acid control from all samples (The background reading can be significant and must be subtracted from sample readings). Then apply the sample readings to the L- amino acid standard curve to obtain the total amino acid amount.

$$\text{L-Amino Acid Concentration} = A/Sv \text{ (nmol}/\mu\text{l or mM)}$$

Where A: L-Amino acid amount (nmol) from the standard curve based on Absorbance OD 570 or fluorescence of your samples.

Sv: Sample volume ( $\mu$ l) you added into the sample wells.



**L-Amino Acid Standard Curve.** Different amounts of L-Amino Acids were measured according to the kit procedure.

### V. Related Products:

Cholesterol/Cholesteryl Esters Quantitation Kit  
HDL, LDL/VLDL cholesterol assay Kit  
Glucose Assay Kit  
Lactate Assay Kit  
Glutathione Assay Kit  
NAD/NADH Assay Kit  
LDH-cytotoxicity Assay Kit  
Uric Acid Assay Kit  
Apoptosis Detection Kits



Problems	Cause	Solution
Assay not working	<ul style="list-style-type: none"> <li>• Use of ice-cold assay buffer</li> <li>• Omission of a step in the protocol</li> <li>• Plate read at incorrect wavelength</li> <li>• Use of a different 96-well plate</li> </ul>	<ul style="list-style-type: none"> <li>• Assay buffer must be at room temperature</li> <li>• Refer and follow the data sheet precisely</li> <li>• Check the wavelength in the data sheet and the filter settings of the instrument</li> <li>• Fluorescence: Black plates (clear bottoms) ; Luminescence: White plates ; Colorimeters: Clear plates</li> </ul>
Samples with erratic readings	<ul style="list-style-type: none"> <li>• Use of an incompatible sample type</li> <li>• Samples prepared in a different buffer</li> <li>• Samples were not deproteinized (if indicated in datasheet)</li> <li>• Cell/ tissue samples were not completely homogenized</li> <li>• Samples used after multiple free-thaw cycles</li> <li>• Presence of interfering substance in the sample</li> <li>• Use of old or inappropriately stored samples</li> </ul>	<ul style="list-style-type: none"> <li>• Refer data sheet for details about incompatible samples</li> <li>• Use the assay buffer provided in the kit or refer data sheet for instructions</li> <li>• Use the 10 kDa spin cut-off filter or PCA precipitation as indicated</li> <li>• Use Dounce homogenizer (increase the number of strokes); observe for lysis under microscope</li> <li>• Aliquot and freeze samples if needed to use multiple times</li> <li>• Troubleshoot if needed, deproteinize samples</li> <li>• Use fresh samples or store at correct temperatures till use</li> </ul>
Lower/ Higher readings in Samples and Standards	<ul style="list-style-type: none"> <li>• Improperly thawed components</li> <li>• Use of expired kit or improperly stored reagents</li> <li>• Allowing the reagents to sit for extended times on ice</li> <li>• Incorrect incubation times or temperatures</li> <li>• Incorrect volumes used</li> </ul>	<ul style="list-style-type: none"> <li>• Thaw all components completely and mix gently before use</li> <li>• Always check the expiry date and store the components appropriately</li> <li>• Always thaw and prepare fresh reaction mix before use</li> <li>• Refer datasheet &amp; verify correct incubation times and temperatures</li> <li>• Use calibrated pipettes and aliquot correctly</li> </ul>
Readings do not follow a linear pattern for Standard curve	<ul style="list-style-type: none"> <li>• Use of partially thawed components</li> <li>• Pipetting errors in the standard</li> <li>• Pipetting errors in the reaction mix</li> <li>• Air bubbles formed in well</li> <li>• Standard stock is at an incorrect concentration</li> <li>• Calculation errors</li> <li>• Substituting reagents from older kits/ lots</li> </ul>	<ul style="list-style-type: none"> <li>• Thaw and resuspend all components before preparing the reaction mix</li> <li>• Avoid pipetting small volumes</li> <li>• Prepare a master reaction mix whenever possible</li> <li>• Pipette gently against the wall of the tubes</li> <li>• Always refer the dilutions in the data sheet</li> <li>• Recheck calculations after referring the data sheet</li> <li>• Use fresh components from the same kit</li> </ul>
Unanticipated results	<ul style="list-style-type: none"> <li>• Measured at incorrect wavelength</li> <li>• Samples contain interfering substances</li> <li>• Use of incompatible sample type</li> <li>• Sample readings above/below the linear range</li> </ul>	<ul style="list-style-type: none"> <li>• Check the equipment and the filter setting</li> <li>• Troubleshoot if it interferes with the kit</li> <li>• Refer data sheet to check if sample is compatible with the kit or optimization is needed</li> <li>• Concentrate/ Dilute sample so as to be in the linear range</li> </ul>
<p><b>Note#</b> The most probable list of causes is under each problem section. Causes/ Solutions may overlap with other problems.</p>		