

- a. The next day, apply desired treatments to the experimental wells omitting the negative and positive control wells. To use Tamoxifen as positive control, dilute Tamoxifen (1000X) directly into the culture medium of positive control wells to obtain 1:1,000 - 1:10,000 dilution. Incubate the plate for the period of time required by your experimental protocol.
- b. Upon completion, gently aspirate off the culture medium from all wells and rinse cells briefly with 200 μ l of 1X PBS. For suspension cells: centrifuge the plate at 200 g (or the lowest centrifuge setting) for 3 minutes to gently deposit the cells onto the surface. Tilt the plate and gently remove the media by aspirating with a pipette tip. Rinse cells briefly with 200 μ l of PBS and spin again. It is important to avoid excessive centrifugation speeds, which can damage the cells. Make note of the place that is used, and perform subsequent aspirations from the same place.

2. Permeabilization and Blocking:

- a. Remove PBS and incubate the cells with 100 μ l of Fixative Solution for 15 min in the dark. Remove Fixative Solution by gentle aspiration for adherent cells, or centrifuge the plate at 200 g for 3 minutes followed by gentle aspiration for suspension cells.
- b. Wash cells three times with 100 μ l of Wash Buffer 5 min each. Remove the Wash Buffer. For suspension cells: centrifuge the plate between each wash at 200 g for 3 minutes.
- c. Incubate cells with 100 μ l of Blocking Buffer for 30 minutes. Remove the Blocking buffer by aspiration after centrifugation at 200 g for 3 minutes. While blocking, prepare the primary antibody and proceed to Immunofluorescence Staining.

3. Immunofluorescence Staining:

Notes:

The recommended dilution for primary and secondary antibodies is 1:100 but it may vary for different cell lines. To prevent cells from drying and photobleaching, the plate should be always covered and protected from light during the incubation periods.

- a. **Primary Antibody Incubation:** Dilute the Phospho-ERK1/2 Primary Antibody 1:100 in Wash Buffer. Add 100 μ l of antibody dilution into each well. Incubate the plate for 2 hours at room temperature, or for best results overnight at 4°C. Remove the antibody by aspiration or centrifugation for suspension cells. Rinse cells briefly three times with 100 μ l of Wash Buffer and remove the washes.
- b. **Secondary Antibody Incubation:** Dilute the Secondary Antibody 1:100 in Wash Buffer. Add 100 μ l of antibody dilution into each well. Incubate the plate for 2 hours at room temperature in the dark, or overnight at 4°C fridge protected from light. Repeat the removal of antibody and wash from step 3a.
- c. **DAPI Staining:** Dilute DAPI stain 1:1000 in Wash Buffer, aliquot 100 μ l to each well and incubate for 10 minutes in the dark. Remove the stain and rinse wells with 100 μ l of Wash Buffer. For removal of DAPI stain and wash follow the steps from paragraph 3a. Add 100 μ l of PBS into each well. Cells are ready to be imaged. For later analysis, store the plate at 4°C in the dark.
- d. Examine the staining under fluorescence microscope with 570 nm excitation and UV laser for Phospho-ERK1/2 Secondary Antibody and DAPI respectively.

Figure 1. 10% FBS Media





Figure 2. 1% FBS Media



Figures: Tamoxifen-induced translocation of phosphorylated ERK1/2 in MCF-7 cells. MCF-7 cells (1×10^5 cells per well) were grown, fixed and stained according to the included protocol. **Figure 1:** Cells grown in media supplemented with 10% FBS and treated with a vehicle (A) or 1X Tamoxifen for 20 min (B). Immunofluorescent staining revealed translocation of phosphorylated ERK1/2 from the cytoplasm (A) to nuclei (B). **Figure 2:** Cells grown in presence of 1% FBS in absence (A) or presence (B) of 1X Tamoxifen for 20 min exhibited translocation of phosphorylated ERK1/2 from nuclei (A) to cytoplasm (B). Bottom panels in Figures 1 and 2 show nuclear staining with DAPI.

VIII. RELATED PRODUCTS:

Phospho-p38 MAPK (Thr180+Tyr182) Translocation Assay Kit (Cell-Based) (K965-100)
Phospho-p53 (Ser15) Translocation Assay Kit (Cell-Based) (K966-100)
Tamoxifen Citrate (1551-1000)

Erk2 Antibody (3442-100)
p53 Nuclear Translocation Assay Kit (Cell-Based) (K961-100)
p53 (human) ELISA Kit (K4829-100)
Phospho-Mek1/2 Antibody (3519-100)
p44/42 MAPK (Erk1/2) Antibody (3085R-100)