



Notes: a) Reducing small molecules in some samples may interfere with PGI assay. We recommend removing the small molecules by ammonium sulfate precipitation method. Ammonium sulfate precipitation: Aliquot 10-100 μ l (~300-500 μ g) of sample to a clean centrifuge tube & add saturated ammonium sulfate (~4.1 M) to final concentration of 3.2 M. Incubate on ice for 20 min. Centrifuge at 14,000 rpm for 5 mins. Suspend the pellet to the original 10-100 μ l volume. b) For unknown samples, we suggest testing several doses to ensure the readings are within the Standard Curve range.

- PGI Positive Control:** Make fresh dilution of Positive Control by adding 2 μ l Positive control to 998 μ l dH₂O. Use 1-10 μ l of diluted Positive Control into the desired well(s) & adjust final volume to 50 μ l with Assay Buffer.
- Reaction Mix:** Mix enough reagents for the number of assays to be performed. For each well, prepare 50 μ l mix containing:

	Reaction Mix	Background Control Mix*
PGI Assay Buffer	44 μ l	46 μ l
PGI Enzyme Mix	2 μ l	2 μ l
PGI Developer	2 μ l	2 μ l
PGI Substrate	2 μ l	---

Add 50 μ l of the reaction mix to each well containing the Standard, Positive Control and test samples and 50 μ l of background control mix to each well containing the background control sample. Mix well.

***Note:** Background control is recommended for samples having high NADH background.

- Measurement:** Incubate for 20-60 min at room temperature and measure OD_{450nm}. **Note:** Incubation time depends on the phosphoglucose isomerase activity in the samples. We recommend measuring OD in a kinetic mode, and choose two time points (T₁ & T₂) in the linear range to calculate the PGI activity of the samples. The NADH Standard Curve can read in endpoint mode (i.e., at the end of incubation time).
- Calculation:** Subtract 0 Standard reading from all standard readings. Plot the NADH Standard Curve. Correct sample background by subtracting the value derived from the background control from all sample readings. Calculate the phosphoglucose isomerase activity of the test sample: $\Delta OD = A_2 - A_1$. Apply the ΔOD to the NADH Standard Curve to get B nmol of NADH generated during the reaction time ($\Delta T = T_2 - T_1$).

$$\text{Sample Phosphoglucose Isomerase activity} = \frac{B}{(\Delta T \times V)} \times \text{Dilution Factor} = \text{nmol/min/ml} = \text{mU/ml}$$

Where: **B** is the NADH amount from Standard Curve (nmol).

V is the sample volume added into the reaction well (ml).

ΔT is the reaction time (min).

Unit Definition: One unit of phosphoglucose isomerase is the amount of enzyme that generates 1.0 μ mol of NADH/min at pH 8.0 at 25°C.




Figure: NADH Standard Curve (a). Phosphoglucose Isomerase activity in bovine liver, HeLa cell lysate & Positive Control (b). Assays were performed following kit protocol.

IX. RELATED PRODUCTS:

Glucose Assay kit

Glucose Dehydrogenase Activity Assay Kit

Glucose-6-Phosphate Dehydrogenase Assay Kit

Glucose Uptake Colorimetric Assay Kit

Glycogen Assay Kit

Maltose and Glucose Assay Kit

NADP/NADPH Quantification Kit

Starch Assay Kit

Glucose and Sucrose Assay Kit

Glucose-1-Phosphate Colorimetric Assay Kit

PicoProbe™ Glucose-6-Phosphate Assay Kit

Glucose Uptake Fluorometric Assay Kit

Hexokinase Colorimetric Assay Kit

NAD/NADH Quantification Kit

Phosphoglucomutase Colorimetric Assay Kit

FOR RESEARCH USE ONLY! Not to be used on humans.