



Beta Galactosidase (β-Gal) Inhibitor Screening Kit (Fluorometric)

6/15

(Catalog # K827-100; 100 assays; Store at -20°C)

I. Introduction:

Beta Galactosidase (β -Gal, EC: 3.2.1.23) is an enzyme which hydrolyzes the β -galactosides into monosacchanides. β -Gal is widely used as a reporter gene in the field of molecular biology. Senescence Associated β -Gal (SA- β -Gal) is an isoform of β -Gal which has the optimal activity at pH 6.0, and is mostly used as a biomarker for senescent cells (K802). β -Gal is an essential enzyme in humans and its deficiency results in Morquio's Syndrome, a severe birth defect. β -Gal can also be used as a tool to study protein-protein interaction. In BioVision's Beta-Galactocidase Inhibitor Screening Kit, β -Gal converts β -Gal substrate to give an intensely fluorescent product (Ex/Em = 480/520 nm). In the presence of a β -Gal inhibitor, the reaction is impeded/abolished resulting in decrease or total loss of fluorescence. This assay kit can be used to screen/study/characterize the potential inhibitors of Beta Galactosidase. The assay is simple, high-throughput adaptable and can be performed within 30 min.

 $\beta \text{-Gal Substrate-Fluorescein} \quad \xrightarrow{\beta \text{-Galactosidase}} \quad \text{Galactose} + \quad \text{Fluorescence} \; (\text{Ex/Em} = 480/520 \; \text{nm})$

II. Application:

• Screening/characterizing/studying potential inhibitors of Beta-Galactosidase.

III. Kit Contents:

Components	K827-100	Cap Color	Part Number
β-Gal Assay Buffer	25 ml	WM	K827-100-1
β-Gal Substrate (in DMSO)	200 µl	Blue	K827-100-2
β-Galactosidase	1 vial	Purple	K827-100-3
β-Gal Inhibitor Control	1 vial	Orange	K827-100-4

IV. User Supplied Reagents and Equipment:

- 96-well white plate with flat bottom
- Multi-well spectrophotometer (fluorescent plate reader)

V. Storage Conditions and Reagent Preparation:

Store kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening. Read entire protocol before performing the assay.

- β-Gal Assay Buffer: Bring to room temperature before use. Store at -20°C or 4°C.
- β-Gal Substrate: Thaw at room temperature. Aliquot and store at -20°C.
- β-Galactosidase: Reconstitute with 550 μl β-Gal Assay Buffer. Aliquot and store at -20°C. Avoid freeze/thaw. Keep on ice while in use. Use within two months.
- β-Gal Inhibitor Control: Reconstitute with 200 μl dH₂O. Aliquot and store at -20°C. Avoid freeze/thaw. Keep on ice while in use. Use within two months

VI. β-Gal Inhibitor Screening Protocol:

1. Screen Compounds, Inhibitor Control, and Enzyme Control Preparation: Dissolve candidate inhibitors into an appropriate solvent at 100X the final concentration to be tested. Dilute to 2X desired test concentration with β-Gal Assay Buffer. Add 50 μl diluted candidate inhibitor or β-Gal Assay Buffer into desired wells, as Sample [S], or Enzyme Control [EC] (no inhibitor). For Inhibitor Control (IC), dilute Inhibitor Control 5 times by adding 20 μl Inhibitor Control to 80 μl β-Gal Assay Buffer. Add 50 μl of diluted Inhibitor Control into desired well(s).

Note: Solvents used to solubilize the inhibitors might affect the enzymatic activity. If solvent effect on enzymatic activity is a concern, prepare a solvent control well(s) with the same final concentration of the solvent(s) as in the inhibitor sample(s) as solvent controls (SC).

- 2. β-Gal Enzyme: Add 5 μl β-Gal Enzyme into Sample, Enzyme Control, and Inhibitor Control wells (if necessary, in Solvent Control wells). Incubate for 5 min. at 25°C. Add 55 μl of Assay Buffer into separate well designated as BC (Background Control).
- **3. Substrate Solution Preparation**: Make enough reagents for the number of assays to be performed. For each well, prepare 45 μl of Substrate solution containing:

β-Gal Assay Buffer 43 μl β-Gal Substrate 2 μl

Mix and add 45 µl of Substrate solution into each well. Mix well with gentle shaking, protected from light.

- **4. Measurement:** Measure fluorescence (Ex/Em = 480/520 nm) in kinetic mode for 5-30 min. at 37°C. Choose two time points (T₁ & T₂) in the linear range of the plot and obtain the corresponding values for the fluorescence (RFU₁ & RFU₂).
- 5. Calculations: Subtract from all samples the background (ΔBC = BC₂-BC₁). Calculate the slope for all Samples (S), including Enzyme Control (EC), by dividing the corrected ΔRFU (RFU₂-RFU₁) values with the time ΔT (T₂-T₁).

% Relative Inhibiton =
$$\frac{\text{(Slope of EC - Slope of S)}}{\text{Slope of EC}} \times 100$$



Gentaur Europe BVBA Voortstraat 49, 1910 Kampenhout BELGIUM Tel 0032 16 58 90 45 <u>info@gentaur.com</u>



Notes:

- a. Irreversible inhibitors that inhibit the β -Gal activity completely at the tested concentration will have Δ RFU = 0 and thus the % Relative Inhibition will be 100%.
- **b.** In case Solvent Control(s) has substantially different slope(s) than the EC, use SC slope(s) instead of Slope of EC in the equation above.




Figure: Inhibition of Beta-Galactosidase activity by β-Gal Inhibitor (b-D-Galactopyranosyl Amine). $IC_{50} = 25.30 \mu M$. Assay was performed following the kit protocol.

VII. Related Products:

Beta-Galactosidase Staining Kit (K802)
Factor Xa Inhibitor Screening Kit (K362)
Myeloperoxidase (MPO) Inhibitor Screening Kit (K746)
pCAF Inhibitor Screening Assay (K345)
HMG-CoA Reductase Activity/Inhibitor Screening kit (K588)
SIRT2 Inhibitor Screening Assay Kit (K322)
Beta-Galactosidase Activity Assay Kit (K821-100)

AHCY Inhibitor Screening Kit (K326) Human Calpain 1 Inhibitor Screening Kit (K244) p300 Inhibitor Screening Kit (K346) TACE Inhibitor Screening Assay Kit (K366) SIRT1 Inhibitor/Activator Screening Kit (K325) SIRT6 Inhibitor Screening Assay Kit (K323)

FOR RESEARCH USE ONLY! Not to be used on humans.