


Cytochrome P450 2C19 (CYP2C19) Activity Assay Kit (Fluorometric)

rev 06/19

(Catalog # K848-100; 100 Reactions; Store at -20°C)

I. Introduction:

Cytochrome P450 2C19 (CYP2C19, EC 1.14.14.1) is a member of the cytochrome P450 monooxidase (CYP) family of microsomal xenobiotic metabolism enzymes. CYPs are membrane-bound heme proteins responsible for Phase I biotransformation reactions, in which lipophilic drugs and other xenobiotic compounds are converted to more hydrophilic products to facilitate excretion from the body. CYP2C19 is primarily expressed in liver and intestinal tissue and catalyzes oxidation of neutral or weakly basic lipophilic molecules with 2-3 hydrogen bond donor/acceptor sites. Isoforms of the CYP2C subfamily are responsible for metabolism of nearly 20% of all small molecule drugs commonly used by humans. Polymorphisms in the human CYP2C19 gene have been implicated in clinical drug/drug interactions involving widely-prescribed drugs, including proton pump inhibitors, antiplatelet agents and anticonvulsants. In addition, for drugs whose pharmacological activity requires metabolism from a pro-drug form, CYP2C19 inhibition or allelic deficiency can lead to decreased drug efficacy. BioVision's CYP2C19 Activity Assay Kit enables rapid measurement of native or recombinant CYP2C19 activity in biological samples such as liver microsomes. The assay utilizes a non-fluorescent CYP2C19 substrate that is converted into a highly fluorescent metabolite detected in the visible range (Ex/Em = 406/468 nm), ensuring a high signal-to-background ratio with little interference by autofluorescence. A highly-selective CYP2C19 inhibitor is provided for determination of CYP2C19 activity in heterogeneous biological samples, where other CYP isozymes may contribute to substrate metabolism. The inhibitor displays greater than 100-fold selectivity for CYP2C19 over other CYPs, ensuring targeted inhibition. CYP2C19 specific activity is calculated by running parallel reactions in the presence and absence of the selective inhibitor and subtracting any residual activity detected with the inhibitor present. The kit contains a complete set of reagents sufficient for performing 50 sets of paired reactions (in the presence and absence of inhibitor).



II. Applications:

- Rapid assessment of native/recombinant CYP2C19 activity in fractions prepared from tissues and cells.
- Screening of drugs and novel ligands for interaction with native/recombinant CYP2C19.

III. Sample Type:

- Human liver microsomes and liver S9 fractions
- Lysates of tissues and cultured cells, primary hepatocytes
- Heterologously expressed recombinant CYP2C19 preparations

IV. Kit Contents:

Components	K848-100	Cap Code	Part Number
CYP2C19 Assay Buffer	100 ml	NM	K848-100-1
3-CHC Standard	1 vial	Yellow	K848-100-2
CYP2C19 Inhibitor ((+)-N-3-benzylnirvanol)	1 vial	Blue	K848-100-3
NADPH Generating System (100X)	1 vial	Green	K848-100-4
β-NADP ⁺ Stock (100X)	1 vial	Amber	K848-100-5
CYP2C19 Substrate	1 vial	Red	K848-100-6
Recombinant Human CYP2C19	1 vial	Violet	K848-100-7

V. User Supplied Reagents and Equipment:

- Multiwell fluorescence microplate reader
- Precision multi-channel pipette and reagent reservoir
- Anhydrous (reagent grade) acetonitrile and DMSO
- Black 96-well plates with flat bottom

VI. Storage Conditions and Reagent Preparation:

Store kit at -20°C and protect from light. Briefly centrifuge all small vials prior to opening. Allow the CYP2C19 Assay Buffer to warm to room temperature prior to use. Read entire protocol before performing the assay procedure.

- **3-CHC Standard:** Reconstitute in 110 μl of DMSO and vortex until fully dissolved to yield a 5 mM stock solution. The 3-CHC stock solution should be stored at -20°C and is stable for at least 3 freeze/thaw cycles.
- **CYP2C19 Inhibitor ((+)-N-3-benzylnirvanol):** Reconstitute in 55 μl of acetonitrile and vortex until fully dissolved to yield a 6 mM stock solution. The stock solution is stable for 2 months at -20°C. To obtain a 150 μM working solution of (+)-N-3-benzylnirvanol (5X final concentration), add 25 μl of the 6 mM stock solution to 975 μl of CYP2C19 Assay Buffer. The 150 μM working solution should be stored at -20°C and used within one month (stable for at least 3 freeze/thaw cycles).

