

Note: It is preferable to centrifuge tubes in a centrifuge with swing-bucket rotor and no brake to prevent disturbance of pellet.

- f. Remove conical tube and observe 4 layers (top to bottom): monocytes, Density Barrier Solution, plasma and RBCs. Collect monocytes that float to the top of the solution. There will be a slight gradient (or change in opacity) in the upper 2-3 ml of the top layer because of the monocytes distributed at the top (**Figure C**). Remove the top 2 ml (containing the monocytes) from each of the three tubes and suspend in a clean 50 ml conical tube. If cells are from the same patient, they can be pooled during the washes.
- g. Dilute the collected monocytes with 18 ml of MIB. Centrifuge cells at 700 x g for 10 minutes, 4°C. Repeat twice and resuspend the pellet in 5 ml of complete media.

2. Determination of Live Cell Count:

- a. Prepare a 1:10 dilution of Viability Stain with cell suspension by adding 18 µl of washed monocyte suspension and 2 µl of Viability Stain to a 1.5 ml centrifuge tube. Inoculate hemocytometer with 10 µl of stained monocyte suspension. Determine and record the total cell count with a Bright-field microscope.
- b. With the same Region of Interest (ROI) in view, reduce white light, open fluorescent lamp shutter and view cells with a FITC/TRITC filter to count the fluorescent cells. *If a small amount of incidental white light illuminates the hemocytometer, the grid will be visible allowing the viewer to see the same ROI as was visible with fluorescent light.* Live cells will fluoresce green. Dead cells will fluoresce red. Tally the number of green and red cells to complete the calculations using the equations below.

3. Evaluation of Purity in Cell Suspension:

- a. Remove 25 µl of monocyte suspension and add to a clean glass slide. Allow the solution on the slide to completely dry (15 min) or heat slide at 37°C to expedite the drying process. Add 50 µl of Giemsa Stain to dried cells and incubate with stain for 5 minutes. Wash slide until runoff is clear. Allow slide to dry. Cells and morphology are best viewed with oil immersion objective lens (**Figure D**). Perform differential cell count using the formulas described in **4. Measurement** (below).

4. Measurement:

a. RBCs, % = $\{(Total\ Cells\ (hemocytometer) - Number\ of\ Fluorescent\ Cells\ (red\ \&\ green))/Total\ Cells\ (hemocytometer)\} \times 100$

b. Live Cells, % = $Number\ of\ Green\ Fluorescent\ Cells / Total\ Number\ of\ Fluorescent\ Cells \times 100$

c. Dead Cells, % = $Number\ of\ Red\ Cells / Total\ Number\ of\ Fluorescent\ Cells \times 100$

Note: Contamination of monocytes with RBCs may affect downstream applications including, but not limited to flow cytometry, chemotaxis assays or cell differentiation protocols. For that reason, we recommend repeating the separation with remaining MDGM if RBCs are >10% of total PBMC count. Alternatively, RBC lysis buffer (Catalog #5830) can be used.

- d. Differential Cell Count: with a 40X objective, count a total of 200 cells. Observe the nucleus of each cell, and score those with unilobar nucleus and vacuolated cytoplasm as a monocyte. Multinucleate cells or those with lobed nuclei are scored as "not monocytes." Monocytes, % = $(Number\ of\ large\ unilobar\ nucleus) / (Total\ Number\ of\ Cells\ (200)) \times 100$.

Note: platelets are a small, anucleate blood component and not included in the differential cell counts.




Figure A. Illustration of conical tube held at 45° angle while DBS is layered on top of Dilute Blood. **B.** Layers of DBS and whole blood prior to and after centrifugation showing the separation of layers in the conical tube. **C.** Separation of four layers (monocytes, DBS, plasma, and RBCs). **D.** Bright-field image of monocytes stained with Giemsa. The large, "fluffy" cytoplasm is approximately two-fold larger in size than the amoeboid nucleus (M, monocyte). **E.** Bright-field image of cells on hemocytometer to determine the Total Cell Count (Left), Image from Fluorescent microscope with FITC/TRITC filters of same ROI showing live (green) and dead (red) cells (middle); merge of two panels (right).

VIII. RELATED PRODUCTS:

Catalog # 5830 1X Red Blood Cell Lysis Buffer

Catalog # 1056 Propidium Iodide

FOR RESEARCH USE ONLY! Not to be used on humans.