




# RNase Activity Detection/Quantification Assay Kit (Fluorometric)

02/19

(Catalog # K934-100; 100 assays; Store at -20°C)

## I. Introduction:

Many molecular biology applications require RNase-free environment due to the degradation of RNA samples by RNase. Obtaining full length, high quality RNA samples are essential, yet challenging in genomics settings. BioVision's RNase Activity Detection Assay Kit enables researchers to measure RNase activity in buffers, reagents, and other components as well as quantitatively evaluate RNase activity of recombinant enzymes in real time. The assay uses a highly sensitive, specific probe that releases a fluorescent product in the presence of active RNase. The limit of detection is 0.4 pg RNase/well and limit of quantification is 1.2 pg RNase/well.



## II. Applications:

- Measurement of RNase contamination in buffers and samples
- Quantitative analysis of RNase activity of purified enzymes

## III. Sample Type:

- Buffers, assay reagents, samples
- Purified enzyme

## IV. Kit Contents:

Components	K934-100	Cap Code	Part Number
10X RNase Buffer	1 ml	Red	K934-100-1
RNA Probe	1 vial	White	K934-100-2
Molecular Biology Grade Water	25 ml	NM	K934-100-3
Half-area Plate	1 plate	--	K934-100-4
RNase Positive Control	250 $\mu$ l	Blue	K934-100-5
Fluorescence Standard (10 mM in DMSO)	100 $\mu$ l	Yellow	K934-100-6

## V. User Supplied Reagents and Equipment:

- Spectrophotometer
- RNase-free barrier pipette tips
- Certified RNase-free reagents, buffers

## VI. Storage Conditions and Reagent Preparation:

**NOTE:** For accurate estimations, it is crucial to use molecular biology grade reagents (RNase free) during sample preparation and RNase barrier filter tips for sample pipetting at all times in order to avoid RNase contamination.

Store kit at -20°C, protected from light. Briefly centrifuge small vials prior to opening. Keep all components on ice while in use. Read entire protocol before performing the assay.

- **10X RNase Buffer, Molecular Biology Grade Water:** Warm to room temperature. Ready to use. Store at RT.
- **RNase Probe:** Reconstitute with 110  $\mu$ l Molecular Biology Grade Water. Aliquot and store at -20°C. Avoid multiple freeze-thaw cycles. Keep away from light.
- **Half-area Plate:** Ready to use. Store at RT.
- **RNase Positive Control:** Ready to use. Aliquot and store at -20°C.
- **Fluorescence Standard (10 mM):** Warm to room temperature. Ready to use. Store at -20°C. Keep away from light.

## VII. RNase Activity Detection Assay Protocol:

**1. Sample Preparation: For solutions with suspected RNase contamination:** Add 2-44  $\mu$ l of solution into a half-area plate. **For surfaces with suspected RNase contaminations:** Put 100  $\mu$ l of the Molecular Biology Grade Water onto the surface. Pipette 44  $\mu$ l of this liquid sample with a pipette and transfer it to a well in the half-area plate. **For purified RNase:** Prepare several dilutions with Molecular Biology Grade Water (provided) to make sure that the kinetic curve fall within the Standard Curve range. Add 2-44  $\mu$ l of sample to each well in the Half-area Plate. **For Positive Control:** Make a 100-fold dilution of the Positive Control by adding 5  $\mu$ l of the RNase Positive Control into 445  $\mu$ l of Molecular Biology Grade Water and 50  $\mu$ l of 10X RNase Buffer. Make a 10,000-fold dilution of the Positive Control by adding 5  $\mu$ l of the 100-fold dilution into 495  $\mu$ l of Molecular Biology Grade Water. Add 2-10  $\mu$ l of the 10,000-fold dilution into a well of the Half-area Plate. **For all Samples:** Bring volume of all Sample wells to 44  $\mu$ l with Molecular Biology Grade Water. **Negative Control (NC):** Aliquot 44  $\mu$ l Molecular Biology Grade Water to a well. Add 6  $\mu$ l 10 X RNase Buffer into all wells containing Samples, Positive Control and Negative Control. *Partial volume for all Sample wells should be 50  $\mu$ l.*

**2. RNA Probe Standard Curve:** Make a 100  $\mu$ M Standard solution by adding 10  $\mu$ l of 10 mM RNase Standard into 990  $\mu$ l Molecular Biology Grade Water. Make a 5  $\mu$ M Standard solution by adding 10  $\mu$ l of 100  $\mu$ M Standard solution into 190  $\mu$ l of Molecular Biology Grade Water. Add 0, 2, 4, 6, 8, 10  $\mu$ l of the 5  $\mu$ M Standard solution into a series of wells in the half-area plate, resulting in 10, 20, 30, 40, 50 pmol RNase Standard per well. Make up the wells to 54  $\mu$ l with Molecular Biology Grade Water. Add 6  $\mu$ l 10 X RNase Buffer into all Standard wells. *Volume for all Standard wells should be 60  $\mu$ l.*

**3. Reaction Mix:** Mix enough reagents for the number of assays to be performed. For each well containing Sample(s), Positive Control, and Negative Control, prepare 10  $\mu$ l Mix containing:

