




α-Glucosidase Inhibitor Screening Kit (Colorimetric)

06/18

(Catalog # K938-100; 100 assays; Store at -20°C)

I. Introduction:

α-Glucosidase (EC 3.2.1.20) is localized in the brush border of the small intestine and is responsible for the enzymatic hydrolysis of 1,4-linked polysaccharides, producing glucose as one of the main products. Due to the vital role of glucose as one of the main sources of energy in eukaryotes, α-Glucosidase is a target for the modulation of postprandial hyperglycemia. α-Glucosidase Inhibitors (AGIs) such as Acarbose, Miglitol and Voglibose are anti-diabetic medicines that help to reduce post-meal blood glucose levels by arresting glucose absorption in the gastrointestinal tract. In addition, recent research is also focused on the discovery of natural products that could act as α-Glucosidase Inhibitors. BioVision's α-Glucosidase Inhibitor Screening Kit can be used to screen potential inhibitors of this enzyme. It utilizes the ability of an active α-Glucosidase to cleave a synthetic substrate thus, releasing a chromophore (OD: 410 nm). In the presence of an α-Glucosidase specific inhibitor, the enzymatic activity is greatly reduced which is detected by a decrease of absorbance readings. The assay kit provides a rapid, simple and reliable test for high-throughput screening of α-Glucosidase inhibitors.



II. Applications:

- Screening/characterizing α-Glucosidase inhibitors

III. Kit Contents:

Components	K938-100	Cap Code	Part Number
α-Glucosidase Assay Buffer	25 ml	WM	K938-100-1
α-Glucosidase Substrate Mix	300 µl	Amber	K938-100-2
α-Glucosidase	1 vial	Blue	K938-100-3
Acarbose	140 µl	Red	K938-100-4

IV. User Supplied Reagents and Equipment:

- 96-well clear plate with flat bottom
- Temperature-controlled plate reader

V. Storage Conditions and Reagent Preparation:

Store kit at -20 °C, protect from light. Briefly centrifuge small vials prior to opening.

- **α-Glucosidase Assay Buffer:** Warm to room temperature before use. Store at 4°C or -20°C.
- **α-Glucosidase Substrate Mix:** Ready to use as supplied. If precipitate is observed, briefly sonicate contents. Store at -20°C.
- **α-Glucosidase:** Reconstitute with 100 µl dH₂O to prepare stock solution. Aliquot Stock Solution in 10 µl aliquots and store at -20 °C. Use aliquot only once. Once aliquoted use within two months.
- **Acarbose:** Ready to use. Keep on ice while in use. Use within two months.

VI. α-Glucosidase Inhibitor Screening Protocol:

1. Screening Compounds, Inhibitor Control & Background Control preparations: Samples [S] and Inhibitor Control [IC]: Dissolve test samples to 100X in a proper solvent. Further dilute to 10X using α-Glucosidase Assay Buffer. Add 10 µl of Diluted test compound, 10 µl of Acarbose into wells of 96-well clear plate designated as test samples [S] or Inhibitor Control [IC], respectively. **Enzyme Control [EC] and Background Control [BC]:** Add 10 and 20 µl of α-Glucosidase Assay Buffer into designated well(s) of 96-well clear plate, respectively. **IC₅₀ estimation (Optional):** prepare several dilutions of candidate(s) in α-Glucosidase Assay Buffer. Add 10 µl of each dilution into designated wells.

Note: Various organic solvents may reduce the α-Glucosidase enzymatic activity. Prepare parallel well(s) as Solvent Control [SC] to test the effect of the solvent on α-Glucosidase activity. If [SC] slope is significantly different when compared to EC, use [SC] values to determine effect of the respective tested compound (see Step 5).

	[S]	[IC]	[EC]	[BC]	[SC]
Test Sample	10 µl	-	-	-	-
Acarbose	-	10 µl	-	-	-
α-Glucosidase Assay Buffer	-	-	10 µl	20 µl	-
Solvent Control	-	-	-	-	10 µl

2. α-Glucosidase Enzyme Solution Preparation: Prepare a 20-fold dilution of α-Glucosidase (i.e. Dilute of 2 µl of α-Glucosidase with 38 µl of α-Glucosidase Assay Buffer), mix thoroughly and keep on ice. Add 10 µl of **Diluted α-Glucosidase Enzyme Solution** to each well containing Test Sample(s) [S], Inhibitor Control [IC], Enzyme Control [EC] and Solvent Control [SC]. Adjust the volume of each well to **80 µl/well** with α-Glucosidase Assay Buffer. **Mix well and incubate at room temperature for 15-20 min. Protect from light.**

Note: Do not store Diluted α-Glucosidase Enzyme Solution. Discard unused solution.

