

2. Standard Curve Preparation:

Colorimetric Assay: Dilute the Pyruvate Standard (100 mM) stock solution to 0.5 nmol/μl by mixing 5 μl of the Standard with 995 μl of ALT Assay Buffer. Add 0, 1, 2, 3, 4, 5 μl into a series of standard wells on a 384-well plate. Adjust the volume to 5 μl/well with ALT Assay Buffer to generate 0, 0.5, 1.0, 1.5, 2, and 2.5 nmol/well of Pyruvate Standard for the colorimetric assay.

Fluorometric Assay: Dilute the Pyruvate Standard (100 mM) stock solution to 1 nmol/μl by mixing 10 μl of the Standard with 990 μl of ALT Assay Buffer. Further dilute the Standard another 10-fold to 0.1 nmol/μl by mixing 10 μl of the 1 nmol/μl solution with 90 μl of ALT Assay Buffer. Add 0, 0.5, 1, 1.5, 2, 2.5 μl into a series of standard wells on a 384-well plate. Adjust the volume to 5 μl/well with ALT Assay Buffer to generate 0, 0.05, 0.1, 0.15, 0.2, and 0.25 nmol/well of Pyruvate Standard for the fluorometric assay.

3. Reaction Mix: Prepare enough Reaction Mix for the number of assays to be performed (including Pyruvate Standard Curve and Positive Control wells). For each well, prepare 25.0 μl Reaction Mix containing:

	<u>Colorimetric</u>	<u>Background</u>	<u>Fluorometric</u>	<u>Background</u>
ALT Assay Buffer	21.5 μl	24.0 μl	21.9 μl	24.4 μl
OxiRed™ Probe	0.5 μl	0.5 μl	0.1 μl	0.1 μl
ALT Enzyme Mix	0.5 μl	0.5 μl	0.5 μl	0.5 μl
ALT Substrate	2.5 μl	---	2.5 μl	---

Add 25.0 μl of the Reaction Mix to each well containing the test samples, Pyruvate Standards, or ALT Positive Control. Add 25.0 μl of the Background Mix to each well containing the background test samples. *The final volume will be 30 μl per well.*

Note: The fluorometric assay is ~10 times more sensitive than the colorimetric assay. Using 0.1 μl of the probe per reaction decreases the background reading and increases detection sensitivity significantly.

4. Measurement: Measure the absorbance (OD_{570 nm}) or fluorescence (Ex/Em = 535/587 nm) in kinetic mode for 60 min or longer at 37 °C. **Notes:** **a.** While the assay can be performed in either end-point or kinetic mode, we strongly recommend reading in kinetic mode in order to ensure that the measurements recorded are within the linear range of the reaction. **b.** Microplate reader settings may need to be adjusted according to the chosen 384-well plate. The dimensions of the used 384-well plate may be available in the manual provided by the plate manufacturer.

5. Calculation: Subtract the 0 Standard reading from all Standard readings and plot the Pyruvate Standard Curve. For each sample type, choose any two time points (t_1 and t_2) in the linear phase of the reaction curve and obtain the corresponding absorbance (A_1 and A_2) or fluorescence (RFU_1 and RFU_2) values at those time points and determine the change in absorbance or fluorescence signal over the time interval: $\Delta OD_{570\text{ nm}} = A_2 - A_1$ or $\Delta F = RFU_2 - RFU_1$. **Note:** Choose time points, which occur after the initial lag phase (roughly 5-10 min in our experience) and during the linear range of probe development (usually within 60 min, samples with extremely low ALT activity may require longer). **If Sample Background Control reading is significant then subtract the Sample Background Control reading from its paired sample readings to get the corrected sample reading.** Apply the sample $\Delta OD_{570\text{ nm}}$ or ΔF values to the Pyruvate Standard Curve to obtain B nmol of pyruvate generated in the sample well during the reaction time ($\Delta t = t_2 - t_1$). Calculate the ALT activity of the test samples using the following equation:

$$\text{Sample ALT Activity} = \frac{B}{(t_2 - t_1) \times V} = \text{nmol/min/ml} = \text{mU/ml}$$


Where: B is the amount of pyruvate in the sample well, calculated from the Standard Curve (in nmol)

t_1 is the time of the first reading (in min)

t_2 is the time of the second reading (in min)

V is the original sample volume added into the reaction well (in ml)

One unit is defined as the amount of ALT which generates 1.0 μmol of pyruvate per minute at 37°C.



Figures: (a) Pyruvate Standard Curve (colorimetric). (b) Pyruvate Standard Curve (fluorometric). (c) ALT Activity (mU/ml) in 3T3 Cell Lysate (2.5 μl, 12 mg/ml protein). Pooled Human Serum (2.5 μl), and Jurkat Cell Lysate (1.25 μl, 4 mg/ml protein). Assays were performed according to the kit protocol.

VII. Related Products:

Alanine Aminotransferase (ALT/SGPT) Activity Assay Kit (K752)	EZScreen™ β-Lactamase Activity Assay Kit - 384 Well (K954)
EZScreen™ NAD ⁺ /NADH Assay Kit - 384 Well (K958)	EZScreen™ Lactate Assay Kit - 384 Well (K951)
EZScreen™ Triglyceride Assay Kit - 384 Well (K952)	EZScreen™ Glycogen Assay Kit - 384 Well (K960)
EZScreen™ Free Fatty Acid Assay Kit - 384 Well (K956)	EZScreen™ ATP Assay Kit - 384 Well (K959)

FOR RESEARCH USE ONLY! Not to be used on humans.