




Cholinesterase Activity Assay Kit (Colorimetric)

rev 06/21

(Catalog # K975-100; 100 assays; Store at -20 °C)

I. Introduction:

Cholinesterase (ChE) consists of a group of enzymes that hydrolyze choline esters. There are two ChE isoenzymes in blood: acetylcholinesterase (AChE; EC 3.1.1.7), also known as erythrocytes or true ChE, which is found mainly in red blood cells; and butyrylcholinesterase (BChE; EC 3.1.1.8), also known as plasma ChE or pseudo-ChE, which is present in plasma. Blood AChE or BChE activity would be selectively reduced by exposing them to poisonous chemical agents, insecticides such as organophosphates or carbamates, anesthetics, and a variety of therapeutic drugs including donepezil or rivastigmine which are used for treating Alzheimer's diseases. Therefore, Blood Cholinesterases (ChE=AChE+BChE) are potential biomarkers of suppressed and/or increased central and peripheral nervous system activity and tools for confirming possible therapeutics. Since plasma BChE and erythrocyte AChE can be selectively inhibited by certain insecticides or drugs, quantification of both isoenzymes' activities is important. **BioVision's cholinesterase activity kit** combines the specific AChE and BChE substrates and a selective BChE inhibitor to measure and distinguish AChE and BChE activities in Whole Blood samples without separating plasma from erythrocytes. The principle is based on the ability of AChE and BChE to hydrolyze their respective substrates and produce thiocholine. Thiocholine reacts with 5,5'-dithiobis(2-nitrobenzoic acid) (DTNB) generating a yellow chromophore (TNB) that can be quantified at 412 nm. It is simple, easy to implement, and useful in clinical research to monitor exposure to anti-ChE compounds in Blood Samples.



II. Applications:

- Measurement of ChE activity in Biological Fluids
- Screening of ChE inhibitors in Biological Samples

III. Sample Type:

- Biological Fluids: Blood

IV. Kit Contents:

Components	K975-100	Cap Code	Part Number
ChE Assay Buffer	100 ml	NM	K975-100-1
AChE Substrate	1 vial	Purple	K975-100-2
BChE Substrate (in DMSO)	100 μ l	Blue	K975-100-3
Acetylcholinesterase	1 vial	Orange	K975-100-4
Butyrylcholinesterase	1 vial	Green	K975-100-5
BChE Inhibitor	1 vial	White	K975-100-6
DTNB	2 vials	Red	K975-100-7
TNB Standard	1 vial	Amber	K975-100-8

V. User Supplied Reagents and Equipment:

- 96-well clear plate with flat bottom
- Multi-well spectrophotometer

VI. Storage Conditions and Reagent Preparation:

Store the kit at -20 °C, protected from light. Briefly centrifuge all small vials prior to opening. Read the entire protocol before performing the assay.

- **ChE Assay Buffer:** Store at 4 °C or -20 °C. Bring to room temperature (RT) before use.
- **AChE Substrate:** Reconstitute the vial in 100 μ l ChE Assay Buffer. Store at -20 °C. Use within two months.
- **BChE Substrate:** Store at -20 °C, protected from light. Bring to RT before use.
- **Acetylcholinesterase:** Reconstitute the vial in 100 μ l of ChE Assay Buffer. Divide into aliquots and store at -20 °C. Use within two months.
- **Butyrylcholinesterase:** Reconstitute the vial in 20 μ l ChE Assay Buffer. Store at -20 °C. Use within two months.
- **BChE Inhibitor:** Reconstitute the vial in 150 μ l **dH₂O**. Vortex intensively at RT to facilitate solubilization. Divide into aliquots and store at -20°C. Bring to RT before use. Use it within two months.
- **DTNB Solution:** Dissolve 1 vial of DTNB in 625 μ l ChE Assay Buffer. Each vial can be used to carry out up to 50 reactions. Dissolve the vial contents when needed. Store at -20 °C. Use within two months.
- **TNB Standard:** Dissolve the vial in 1 ml of ChE Assay Buffer to generate 2.5 mM TNB Standard. Use within two months.

VII. Cholinesterase Assay Protocol:

1. BChE inhibitor, Sample Preparations:

a. BChE inhibitor: Dilute BChE Inhibitor 15-fold (i.e. Dilute 10 μ l BChE Inhibitor with 140 μ l ChE Assay Buffer).

b. **Blood sample:** Prepare a 40-200 fold dilution of Blood in **dH₂O**. *Record Dilution Factor*. Add 10-20 μ l of Diluted Blood into 3

