


Aromatase (CYP19A) Inhibitor Screening Kit (Fluorometric)

rev 06/19

(Catalog # K984-100; 100 Reactions; Store at -20°C)

I. Introduction:

Aromatase (CYP19A, EC 1.14.14.14) is a member of the cytochrome P450 monooxidase (CYP) family of microsomal xenobiotic metabolism enzymes. Aromatase plays a critical role in steroidogenesis, catalyzing the conversion of androgenic hormones into estrogens. The enzyme is expressed in high levels in reproductive tissues, placenta, brain and adipose tissue and is responsible for mammalian sexual dimorphism and development of secondary sexual characteristics. Inhibitors of aromatase are used to treat estrogen-dependent breast cancer, as estrogens promote the expression of peptide growth factors responsible for tumorigenesis. Aromatase activity and expression can be affected by many organic environmental pollutants such as pesticides and plasticizers. Such compounds, known as endocrine disruptors, are suspected of causing precocious puberty, obesity, infertility and various cancers. BioVision's Aromatase Inhibitor Screening Kit enables rapid screening of drugs and other small molecules for compound-Aromatase interaction in a reliable, high-throughput fluorescence-based assay. The assay utilizes a fluorogenic substrate that is converted into a highly fluorescent metabolite detected in the visible range (Ex/Em = 488/527 nm), ensuring a high signal-to-background ratio with little interference by autofluorescence. The kit contains a complete set of reagents sufficient for performing 100 reactions in a 96-well plate format.



II. Applications:

- Rapid, high-throughput screening and characterization of drugs and novel ligands for interaction with human aromatase.
- Development of structure-activity relationship (SAR) models to predict aromatase inhibition liability of novel compounds and analogues.
- Prediction of aromatase-dependent endocrine disruption potential for novel compounds.

III. Kit Contents:

Components	K984-100	Cap Code	Part Number
Aromatase Assay Buffer	100 ml	NM	K984-100-1
Fluorescence Standard (1 mM)	50 µl	Yellow	K984-100-2
Aromatase Inhibitor (Letrozole)	1 vial	Blue	K984-100-3
NADPH Generating System (100X)	1 vial	Green	K984-100-4
β-NADP ⁺ Stock (100X)	1 vial	Amber	K984-100-5
Aromatase Substrate	1 vial	Red	K984-100-6
Recombinant Human Aromatase	2 vials	Violet	K984-100-7

IV. User Supplied Reagents and Equipment:

- Multi-well fluorescence microplate reader
- Precision multi-channel pipette and reagent reservoir
- Anhydrous (reagent grade) acetonitrile and DMSO
- White 96-well plates with flat bottom

V. Storage Conditions and Reagent Preparation:

Store kit at -20°C and protect from light. Briefly centrifuge all small vials prior to opening. Allow the Aromatase Assay Buffer to warm to room temperature prior to use. Read entire protocol before performing the assay procedure.

- **Fluorescence Standard (1 mM):** Provided as a 1 mM stock solution in DMSO. Prior to use, warm to room temperature and vortex until fully dissolved. Store at -20°C, stable for at least 3 freeze/thaw cycles.
- **Aromatase Inhibitor (Letrozole):** Reconstitute in 55 µl of acetonitrile and vortex until fully dissolved to yield a 1 mM stock solution. To obtain a 5 µM working solution of letrozole (5X final concentration), add 5 µl of the 1 mM stock solution to 995 µl of Aromatase Assay Buffer. The 5 µM working solution should be stored at -20°C and is stable for 2 freeze/thaw cycles. The stock solution is stable for 2 months at -20°C.
- **NADPH Generating System (100X):** Reconstitute with 220 µl Aromatase Assay Buffer, aliquot and store at -20°C. Avoid repeated freeze/thaw cycles and keep on ice while in use.
- **β-NADP⁺ Stock (100X):** Dissolve in 110 µl Aromatase Assay Buffer and vortex thoroughly to yield a 100X stock solution of NADP⁺. Store at -20°C, stable for at least 3 freeze/thaw cycles.
- **Aromatase Substrate:** Reconstitute with 55 µl anhydrous reagent-grade acetonitrile and vortex until fully dissolved to obtain a 1 mM stock solution. Store at -20°C. Allow the vial to warm to room temperature before opening and promptly retighten cap after use to avoid absorption of airborne moisture.

$= (RFU_2 - RFU_1)$ and $\Delta T = (T_2 - T_1)$. Calculate the rate of change in fluorescence over time according to the equation below. Subtract the rate of the no substrate/background control (BC) well from the rates of each of the no inhibitor/solvent control (R_{SC}) and test compound (R_{TC}) wells to determine background-corrected reaction rates (denoted by R) for each well:

$$R = \frac{\Delta F - \Delta F_{BC}}{\Delta T}$$

Calculate the percent inhibition due to the test ligand or positive inhibition control using the following equation:

$$\% \text{ Relative Inhibition} = \frac{R_{SC} - R_{TC}}{R_{SC}} \times 100\%$$

Notes:

- The Aromatase Substrate undergoes rapid photobleaching in aqueous solutions. This photophysical property may give rise to a sharp non-linear phase in first few minutes of the reaction progress curves (a lag phase). When calculating ΔF values, it is important to choose time points that occur *after* this initial lag phase, during the linear range of the reaction. *In our experience, the linear phase begins roughly 5-10 mins after the initiation of the reaction.*
- If the background control (BC) well rate calculation yields a negative value, subtraction of the BC value may be ignored.
- If desired, reaction rate calculations can also be expressed in terms of pmoles of fluorescent metabolite formed per unit time per unit amount of protein by interpolation from the standard curve. Each well will contain a total of 50 μg of protein when the Recombinant Human Aromatase is used at the proportions suggested in the kit protocol.




Figure: (a) Standard curve of Aromatase Substrate metabolite fluorescence. One mole of fluorescence standard corresponds to the metabolism of one mole of Aromatase Substrate. (b) Reaction kinetics of recombinant human aromatase enzyme at 37°C in the presence and absence of the indicated aromatase inhibitors (the solvent control reaction contained assay buffer with 0.2% acetonitrile). (c) Dose-response curves for various aromatase ligands of differing structural and mechanistic classes: the potent competitive inhibitor letrozole, the mechanism-based irreversible inhibitor exemestane, the phytoestrogen-like synthetic flavonoid α -naphthoflavone and the antifungal ketoconazole (a competitive inhibitor of several CYP isoforms). For dose-response curves, percent activity was calculated for each concentration of inhibitor by comparison to activity of reactions containing no inhibitor. For each inhibitor, IC_{50} values were derived by 4-parameter logistic curve fitting with each point representing the mean \pm SEM of at least four replicates. Assays were performed according to the kit protocol.

VII. RELATED PRODUCTS:

Microsome Isolation Kit (K249)
Cytochrome P450 Reductase Activity Kit (K700)
Cytochrome P450 3A4 Activity Assay Kit (K701)
Cytochrome P450 3A4 Inhibitor Screening Kit (K702)
Cytochrome P450 2D6 Activity Assay Kit (K703)
Cytochrome P450 2D6 Inhibitor Screening Kit (K704)
Cytochrome P450 2C19 Activity Assay Kit (K848)
Cytochrome P450 2C19 Inhibitor Screening Kit (K849)
Cytochrome P450 1A2 Activity Assay Kit (K893)
Cytochrome P450 1A2 Inhibitor Screening Kit (K894)

Cytochrome P450 2C9 Activity Assay Kit (K895)
Cytochrome P450 2C9 Inhibitor Screening Kit (K896)
Aromatase (Human) ELISA Kit (K3599)
Exemestane (9473)
EZCyp™ Active human Cytochrome P450 3A4 (7872)
EZCyp™ Active human Cytochrome P450 2D6 (7873)
EZCyp™ Active human Cytochrome P450 2C19 (7874)
EZCyp™ Active human Cytochrome P450 2C9 (7875)
EZCyp™ Active human Cytochrome P450 1A2 (7876)
UGT Activity Assay / Ligand Screening Kit (K692)

FOR RESEARCH USE ONLY! Not to be used on humans.