



**Notes:**

- a. Taurine concentration varies over a wide range depending on the sample. For unknown samples, we recommend doing a pilot experiment & testing several doses to ensure the readings are within the Standard Curve range.
  - b. For samples having high protein content, we recommend deproteinizing the samples (tissue or cell lysate or biological fluids) using 10 kDa Spin Column (Cat. # 1997 or equivalent). Add sample to the spin column, centrifuge at 10,000 X g, 4°C for 10 min. collect and use the filtrate.
  - c. To ensure accurate determination of Taurine in the test samples or for samples having low concentrations of Taurine, we recommend spiking samples with a known amount of Taurine Standard (e.g. 15 nmol).
  - d. Thiol present in biological samples might show high background signal. To quantify the signal contribution from taurine-generated sulfite only, add formaldehyde to 5 mM final concentration *before adding the probe* (step 4). This will complex with the sulfite and prevent signal generation. The difference of signal (with and without formaldehyde) will correspond to the signal from sulfite only.
- 2. Taurine Standard Curve Preparation:** Add 0, 2, 4, 6, 8 and 10 µl of 5 mM Taurine into a series of wells in a 96-well plate and adjust the final volume to 180 µl/well with WS to generate 0, 10, 20, 30, 40 and 50 nmol/well of Taurine Standards respectively.
- Note:** For unknown samples, we suggest testing several dilutions to ensure the readings are within the Standard Curve range.
- 3. Enzyme mix:** To each well of Sample(s) and Taurine Standards, add 20 µl of Enzyme Mix, mix properly, and then incubate the 96 well plate at 30°C for 30 mins. For BC, replace Enzyme mix with 20 µl of Taurine Assay Buffer.
- 4. Sulfite Probe:** add 30 µl of the Sulfite Probe to each well containing Sample, Background Control, and Taurine Standards Mix well, incubate for 5 minutes at 30°C.
- 5. Measurement:** Measure absorbance in an endpoint mode at 415 nm using a microplate reader.
- 6. Calculation:** Subtract 0 Standard reading from all readings. Plot the Taurine-Standard Curve and obtain the slope of the curve ( $\Delta$ Absorbance/nmol). If the background control reading is significant then subtract the background control reading from sample reading.

$$[\text{Taurine}] = B \times D / V = \text{nmol/ml}$$

Where:

**B** = Taurine in sample based on Std. curve slope (nmol)


**V** = the sample volume added into the reaction well (ml)

**D** = Sample dilution factor (D=1 when samples are undiluted)

**Note:** For spiked samples, correct for any sample interference by subtracting the sample reading from spiked sample reading.

$$\text{For spiked samples Taurine amount in sample well} = \frac{(\text{OD}_{\text{Sample corrected}})}{(\text{OD}_{\text{Sample+Tau Standard}}) - (\text{OD}_{\text{Sample (corrected)}})} * \text{Taurine spike (nmole)}$$

Taurine Molecular weight: 125.15.



**Figure:** (a) **Taurine Standard Curve** (5 nmoles-25 nmoles), error bars indicate SD (n=3). (b) **Taurine in urine and energy drink:** Taurine present in Donor 1 urine sample was 134 nmoles/ mg Creatinine (Cr) and in Donor 2 was 592 nmoles/ mg Creatinine. Energy drink was diluted 10x with water before performing the assay.

**VIII. RELATED PRODUCTS:**

- Taurine Dioxygenase (P1071)
- Cysteine Assay Kit (Fluorometric) (K558)
- Glycine Assay Kit (Fluorometric) (K589)
- Glutamate Colorimetric Assay Kit (K629)
- Sarcosine Colorimetric/Fluorometric Assay Kit (K636)
- Albumin (Albuminuria) Fluorometric Assay Kit (K550)
- Creatinine Colorimetric/Fluorometric Assay Kit (K627)

**FOR RESEARCH USE ONLY! Not to be used on humans.**